



# Treatment Tables




# TREATMENT

# TABLES

Intensa offers a collection of treatment tables *to meet any need* in the exam room.

The tables are offered with an adjustable or fixed back and contoured or flat top. Each style of table can have up to two pass-through drawers. All of our tables come in the standard Foghorn Powder Coat Finish and have a weight capacity of 500 lbs. We also offer an extensive collection of vinyls to add color to the treatment tables.



**i** Special Feature

*With seamless upholstery for a smooth, easy-to-clean surface, our comfortable foam tops are contoured and curved to the body's natural shape.*

CHOOSE FROM A VARIETY  
OF STYLES AND OPTIONS

## 415 SERIES

### Standard Features & Specifications

Shelf | Adjustable back | Contour top  
Overall Length 76" | Height 32" | Width 28"



415



415-D



415-DD

### Special Feature

*Our pass-through drawers can be opened from either side of the table, so you'll always have easy access to supplies no matter where your table sits. Drawer stops add further flexibility, allowing the pass-through drawers to become standard drawers that will not accidentally knock into anything on the other side.*

## 414 SERIES

### Standard Features & Specifications

Shelf | Adjustable back | Flat top  
Overall Length 76" | Height 32" | Width 28"



414



414-DD



414-D

## 413 SERIES

### Standard Features & Specifications

Shelf | Fixed back | Flat top  
Overall Length 76" | Height 32" | Width 28"



413-DD



413-D



413

## 412 SERIES

### Standard Features & Specifications

Adjustable back | Contour top  
Overall Length 76" | Height 32" | Width 28"



### Special Feature

*The adjustable back option offers maximum flexibility and ease of access for healthcare professionals while providing a stable, comfortable surface to support patients.*

412



412-DD



412-D

## 411 SERIES

### Standard Features & Specifications

Adjustable back | Flat top  
Overall Length 76" | Height 32" | Width 28"



411



411-DD




411-D

## 410 SERIES

### Standard Features & Specifications

Fixed back | Flat top  
Overall Length 76" | Height 32" | Width 28"



410



410-DD



410-D



Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262  
888.636.2660 | sales@intensa.net | [intensa.net](http://intensa.net)

*Manufactured in the USA*

