



TREATMENT ACCESSORIES

140 SERIES MULTI PURPOSE CART

144 Folding Utility Cart - Gray Epoxy Coated Baskets *each sold separately*
Ideal for delivering medical records and supplies.

144 Features:

- Multi purpose, rolling utility cart - *Sold Separately*
- Durable 1 inch (2.54 cm) diameter chrome plated steel frame
- 3 inch (7.62 cm) ball bearing rubber wheels
- Folds for easy storage
- Gray epoxy coated frame support (for basket not Included)
- Assembly Required
- Folding cart ships separately and must be ordered separately.

144-Basket Features:

- Gray epoxy coated wire baskets - *Sold Separately*
- Wire baskets are ship separately
- Baskets must be ordered separately.

144 Product Info

144 Options:

	item no.
Rolling Utility Cart	144
Wire Basket [†] - 6 inch <i>Order Separately</i>	144-6
Wire Basket [†] - 12 inch <i>Order Separately</i>	144-12

sold and ship separately

144 Specifications:

	inch	cm		inch	cm			
Overall Width	19.50	49.53	Rubber Wheels ball bearing	3.00	7.26	Linen Bag Included		no
Overall Length	34.00	86.36	Steel Frame Tubing	1.00	2.54	Ships Assembled		no
Overall Height	40.00	101.6	Steel Frame Finish	chrome plated			lb	kg
			Wire Basket Finish	gray epoxy			Shipping Weight	23

Gray Epoxy Coated Baskets[†]

	inch	cm		inch	cm		inch	cm
Single Height [†]	144-6		Width	14.50	36.83	Length	24.50	62.23
			Height	6.00	15.24			
Double Height [†]	144-12		Width	14.50	36.83	Length	24.50	62.23
			Height	12.00	30.48			

[◊] Ships: Ground at no additional charge [†] Both Cart and Basket must be ordered and ship separately



Multi Purpose, Chrome Plated Steel Cart shown with Gray Epoxy Coated Baskets - each basket and cart are sold and ship separately



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

Color variations may occur due to monitor or printer settings.

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net



Manufactured in Taiwan