

## ACCESSORIES

### 170 SERIES

# DUAL TOP OVERBED TABLE

#### Overbed Table with Heavy Duty Steel, Low Profile Base and Thermofoil Dual Top

This Overbed Table combines an attractive, ergonomic design with features that reduce infection control issues by minimizing areas in which dust and dirt can accumulate. Manufactured with the lowest base in the industry, our Overbed tables move easily under any hospital bed or recliner.

#### 179 Features:

- Supports up to 100 lb of evenly distributed weight
- Exposed steel base with shock-absorbent gas spring cylinder
- Seamless spill guard top with integrated cup holders
- Laminate or Thermofoil tops are customizable in shape, finish, and edge
- Low profile design allows overbed table to be used with limited clearance beds

#### 179 Product Info

#### 179 Options: \*\*

	item no.	kidney
Thermofoil	179TF-y-xx	179TFK-y-xx


options available at additional cost

#### 179 Specifications:

	inch	cm		inch	cm
Overall Width	32.00	81.28	Base Width	32.00	81.28
Overall Depth	22.00	55.88	Base Depth	18.00	45.72
Overall Height <i>minimum</i>	29.00	73.66	Base Profile Height	1.94	4.93
Overall Height <i>maximum</i>	44.00	111.76	Top Width	32.00	81.28
			Top Depth	17.00	43.18

\*-xx in the product number is the placeholder for top finish choice    †-y in the product number is the placeholder for base paint finish choice  
Solid Surface top available - contact customer service    \*Ships: Freight request additional shipping quote from Intensa

#### 170 OBT Top Finishes (xx)



Candlelight | CL    Chocolate Pear | CP    Hard Rock Maple | HM    Medium Cherry | MC    Medium Oak | MO


#### 170 Paint Finishes (y)



Almond | A    Graphite Gray | G    Nevada | N



	lb	kg
Weight Capacity EVENLY DISTRIBUTED	100	45.36
Shipping Weight	43	19.50



#### 179 Accessories:

	item number
Vanity w/Mirror Basket 10" x 6" x 9" inch (25.40 x 15.24 x 22.86 cm)	179-VM
Basket	179-BK

accessories available at additional cost



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in Canada

