

## EXAM TABLES

# 420 SERIES FIVE DRAWER EXAM TABLE

### Exam Table with Three Pass Through Drawers, Two Storage Drawers and Stirrups

An exam table with pass through drawers, storage drawers and 5 position stirrups. All metal construction of the base table provides for longevity. of the base provides for longevity. Intensa offers an extensive collection of vinyls to add color to the exam room.

#### 420 Features:

- Three Pass Through Drawers & Two Storage Drawers
- Pull Out Leg Rest with Integrated Treatment Tray
- Five Position Continuous-Lock Stirrups
- Powder Coat Finish Foghorn
- Rubber Grip Pull Out Step
- Pelvic Tilt Standard
- Paper Roll Holder

 CF Stinson Pre-Approved Fabrics

 420 Product Info




#### 420 Options: †

	item no.
Five Drawer Exam Table	460-xx
Five Drawer Exam Table <small>with integrated electrical outlet</small>	460-E0-xx
Five Drawer Exam Table <small>with drawer warmer with integrated outlet</small>	460-DW-xx

options available at additional cost

#### 420 Specifications:


	inch	cm		inch	cm		yd	m
Overall Width	27.00	68.58	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement <sup>1</sup>	3	2.74
Overall Length <small>leg rest extended</small>	72.00	182.88		19.00	48.26	Pull Out Leg Rest	standard	
Overall Height	32.00	81.28		20.00	50.80	Treatment Tray	integrated	
Drawer Width <sup>7</sup>	18.00	45.72	Foot Step Width	18.00	45.72	 FREIGHT*	lb	kg
Drawer Length <sup>7</sup>	16.50	41.91	Foot Step Length	17.00	43.18	Foot Step Capacity EVENLY DISTRIBUTED	800	362.87
Drawer Height <sup>7</sup>	4.00	10.16	Stirrups	five position continuous lock		Weight Capacity EVENLY DISTRIBUTED	500	226.80
Pelvic Tilt	standard		Powder Coat Finish	foghorn		Shipping Weight	335	151.95

\*Ships: Freight request additional shipping quote from Intensa †-xx in the product number is the placeholder for top finish choice.  
<sup>1</sup> Submit COM sample for approval <sup>7</sup> Interior dimensions



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

