

EXAM & TREATMENT TABLES

RT404 SERIES

REPLACEMENT TOP

Replacement Top with Foot Pad for Ritter® by Midmark® 404 Exam Tables

Intensa offers a replacement top program for your existing exam tables. We provide tops that meet certain manufacture model number specifications. The new top is made of completely new foam and wood substrate and can be upholstered in a variety of vinyl colors and patterns. Revitalize/renew your exam rooms without the expense of new tables!

RT404 Features:

- Replaces Ritter® by Midmark® 404 Exam Table
- Foot pad included
- Includes hook and loop attachment

 [CF Stinson Pre-Approved Fabrics](#)

 [Replacement Tops Overview](#)


 [RT404 Product Info](#)

RT404 Options:

	Item no. †
Replacement Top	RT404-xx



RT404 Specifications:


	inch	cm		inch	cm		yd	m
Overall Width	27.50	69.85	Foot Rest Width	20.50	52.07	COM Requirement ¹	3	2.74
Overall Length	60.00	152.40	Foot Rest Length	17.50	44.45		lb	kg
			Hook & Loop Attachment	included		Shipping Weight	35	15.88

[‡] Ships: Ground at no additional charge † -xx in the product number is the placeholder for vinyl choice ¹ Submit COM sample for approval



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

