

## TREATMENT TABLES

# 411 SERIES EXAM ROOM TABLE

**Exam Room Treatment Table with Adjustable Back**  
Intensa offers a collection of treatment tables to meet any need in the exam room. The tables are offered with an upholstered back and an option for zero, one or two drawers. Intensa offers an extensive collection of vinyls to add color to the treatment table.

### 411 Features:

- Powder coat frame finish - Foghorn
- Optional one or two pass through drawers
- Adjustable back

 CF Stinson Pre-Approved Fabrics

 411 Product Info



411

### 411 Options:

	item number <sup>†</sup>
No Drawers	411-xx
Single Drawer	411D-xx
Double Drawer	411DD-xx

*options available at additional cost*



411D

### 411 Specifications:

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder	18.00	45.72	COM Requirement <sup>1</sup>	4	3.66
Overall Length	76.00	193.04	Widths Available	19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width <sup>7</sup>	17.75	45.09	Foam Thickness	1.75	4.45	Weight Capacity	500	226.80
Drawer Length <sup>7</sup>	19.25	48.90				EVENLY DISTRIBUTED		
Drawer Height <sup>7</sup>	5.00	12.70	Powder Coat Finish		foghorn	Shipping Weight	190	86.18

<sup>\*</sup> Ships: Freight request additional shipping quote from Intensa    <sup>†</sup> -xx in the product number is the placeholder for vinyl choice    <sup>1</sup> Submit COM sample for approval  
<sup>2</sup> Adjust Shipping Weight +10 lb (+4.54 kg) per drawer    <sup>3</sup> Interior Dimensions



411DD



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

Manufactured in the USA

