

TREATMENT TABLES

412 SERIES EXAM ROOM TABLE

Exam Room Treatment Table with Adjustable Back and Contour Top
Intensa offers a collection of treatment tables to meet any need in the exam room. The tables are offered with an upholstered back and an option for zero, one or two drawers. Intensa offers an extensive collection of vinyls to add color to the treatment table.

412 Features:

- Powder coat frame finish - Foghorn
- Optional one or two pass through drawers
- Adjustable back
- Contour top

 CF Stinson Pre-Approved Fabrics

 412 Product Info



412

412 Options:

	item no.*
No Drawers	412-xx
Single Drawer	412D-xx
Double Drawer	412DD-xx

options available at additional cost




412D



412DD

412 Specifications:


	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹	4	3.66
Overall Length	76.00	193.04		19.00	48.26	 FREIGHT*		
Overall Height	32.00	81.28		20.00	50.80		Weight Capacity EVENLY DISTRIBUTED	500
Drawer Width ³	17.75	45.09	Foam Thickness	1.75	4.45	Shipping Weight ²	170	81.65
Drawer Length ³	19.25	48.90						
Drawer Height ³	5.00	12.70	Powder Coat Finish	foghorn				

*Ships: Freight request additional shipping quote from Intensa †-xx in the product number is the placeholder for vinyl choice ¹ Submit COM sample for approval
² Adjust Shipping Weight +10 lb (+4.54 kg) per drawer ³ Interior Dimensions



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

