

TREATMENT TABLES

415 SERIES EXAM ROOM TABLE

Exam Room Treatment Table with Shelf, Adjustable Back and Contour Top
Intensa offers a collection of treatment tables to meet any need in the exam room. The tables are offered with an upholstered back and an option for zero, one or two drawers. Intensa offers an extensive collection of vinyls to add color to the treatment table.

415 Features:

- Powder coat frame finish - Foghorn
- Optional one or two pass through drawers
- Adjustable back
- Contour top

 CF Stinson Pre-Approved Fabrics

 415 Product Info



415

415 Options:

	item number [†]
No Drawers	415-xx
Single Drawer	415D-xx
Double Drawer	415DD-xx

options available at additional cost




415D



415DD

415 Specifications:


	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80		Weight Capacity EVENLY DISTRIBUTED	500
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	Shipping Weight	210	95.25
Drawer Length ⁷	19.25	48.90						
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn				

^{*}Ships: Freight request additional shipping quote from Intensa [†]-xx in the product number is the placeholder for vinyl choice ¹ Submit COM sample for approval
² Adjust Shipping Weight +10 lb (+4.54 kg) per drawer ³ Interior Dimensions



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

