

## SEATING

# 950 SERIES PHYSICIAN EXAM STOOL

### 950 Features:

- Single lever adjustment release
- Pneumatic height adjustment
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 **CF Stinson** Pre-Approved Fabrics

 **950 Options** Detail and Photos

Product Info:

 **951**  **958**




### 950 Options:

base with standard seat height	item no. <sup>†</sup>	procedure arm <sup>‡</sup>	with footring <sup>‡</sup>	with backrest <sup>‡</sup>	crescent back <sup>‡</sup>
<b>Black Nylon Composite</b>	951-xx	951PA-xx	951-FR-xx	951B-xx	951CB-xx
<b>Five Leg Tubular Steel<sup>4</sup></b>	958-xx	958PA-xx	standard	958B-xx	958CB-xx

*Footring and Glides*


*options available at additional cost*

### 950 Specifications:

	inch	cm		inch	cm		yd	m
<b>Seat Diameter</b>	16.00	40.64	<b>Base Diameter</b> <i>with casters</i>	24.00	60.96	<b>COM Requirement<sup>1</sup></b>	1	0.92
<b>Overall Height</b> <i>with backrest -low</i>	33.00	83.82	<b>Height Adjustment</b> <i>method</i>	pneumatic			lb	kg
<b>Overall Height</b> <i>with backrest -high</i>	39.00	99.06	<b>Height Release</b> <i>hardware</i>	lever		<b>Weight Capacity</b> EVENLY DISTRIBUTED	350	158.76
						<b>Shipping Weight</b>	20	9.07

### Adjustable Height Seat<sup>\*\*</sup>


	inch	cm		inch	cm
<b>Standard</b>	95x-xx		<b>Seat Height -low</b>	18.50	46.99
<b>7 [seven]</b>	95x-7-xx		<b>Seat Height -low</b>	20.50	52.07
<b>10 [ten]</b>	95x-10-xx		<b>Seat Height -low</b>	24.50	62.23
			<b>Seat Height -high</b>	24.00	60.98
			<b>Seat Height -high</b>	28.25	71.76
			<b>Seat Height -high</b>	35.00	88.90

<sup>♦</sup> Ships: Ground at no additional charge <sup>†</sup> -xx in the product number is the placeholder for vinyl choice  
<sup>‡</sup> Replacement Cylinders may be ordered separately at additional cost   
<sup>4</sup> 20 inch (50.80cm) footring; casters add 2" (5.08 cm) in height


<sup>2</sup> Specification Variances

item	seat height		ship weight	
	inch	cm	lb	kg
958	-1	-2.54	+3	+1.36



 GF Health Products, Inc. is an  
 intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
 Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA  
 +1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

Manufactured in the USA

