



Product Specification & Pricing Catalog

published *July 2026*





Put a little **COLOR**
in your exam room with Intensa!



[California Prop 65 Warnings](#)

[Intensa Limited Warranties](#)

[Intensa Terms & Conditions](#)

[Green Initiative](#)





Table of Contents



<i>Physician Exam</i> Seating	4
<i>Laboratory</i> Seating	10
<i>Hydraulic Surgical</i> Seating	19
<i>Clinical Care</i> Recliners	20
<i>Blood Draw / Phlebotomy</i> Seating	24
<i>Exam Room</i> Treatment Tables	25
<i>Replacement</i> Exam Table Tops	28
<i>Patient / Guest</i> Seating	30
<i>Overbed</i> Tables	31
<i>Step</i> Stools	32
<i>Exam</i> Lights	33
<i>IV</i> Poles	33
<i>Ordering</i> Info & Notices	34
<i>Stool</i> Options	34
<i>Recliner</i> Options & Accessories	35
<i>OBT Finish</i> Swatches, Care & Cleaning	36
<i>Approved Vinyl Upholstery</i> Swatches	37
<i>Vinyl Upholstery</i> Tier, Care, & Link Chart	42
<i>Vinyl Upholstery</i> Care, Cleaning & Disinfection	43
<i>Upholstery</i> Composition, Finish, & Performance Comparison	44
<i>Product</i> Index	45

🔗 [click above listing jumps to page](#)

Purchase Orders

please remit purchase orders to intensa@grahamfield.com

COM Requirements

submit required sample size indicated to: Intensa COM Department, 1810 S. Elm Street, High Point, NC 27260

Pricing & Availability

current at time of publication & subject to change without notice. Intensa is not responsible for typographical errors

Contact Customer Service

for most current information, pricing, freight quotes, and support +1.336.884.4096 | intensa@grahamfield.com

Intensa products are labeled in accordance with California Proposition 65 Warnings for products potentially sold to California consumers



WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov



WARNING: This product can expose you to chemicals including Formaldehyde (gas) which is known to the State of California to cause cancer. For more information go to - www.P65Warnings.ca.gov/furniture



WARNING: This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP) which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to - www.P65Warnings.ca.gov/furniture

© 2015, GF Health Products, Inc. All Rights Reserved | Intensa, Inc. is a wholly owned subsidiary of GF Health Products, Inc. | Stinson and Zeron are registered trademarks of C. F. Stinson, LLC. | Intensa, Dura-Gel, and SOCK-IT-GUARD are registered trademarks of GF Health Products, Inc. | Spradling, Neo, EcoSense, PBG+, Permaguard, Permablok, Permablok Plus, and Sureguard are trademarks of Spradling Holdings, Inc. | Crypton is a registered trademark of Crypton LLC | Tente is a registered trademark of TENTE International | Formula 409 is a registered trademark of The Clorox Company. | Fantastic is a registered trademark of S.C. Johnson & Son, Inc.



993 Physician Stool with Chrome 360° Hand Release and Polished Aluminum Base



Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$539	\$539	\$619	\$639	\$652	\$698	\$723	\$739	\$746	\$781	\$834
	inch	cm		inch	cm		yd	m			
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92			
Overall Height -low	18.00	45.72	Height Adjustment method	pneumatic			lb	kg			
Overall Height -high	23.00	55.88	Height Release hardware	360° hand			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76		
Optional Footring	18.00	45.72					Shipping Weight ⁸	23	11.34		

Options:

[details & photos on page 34*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		993-FR-xx	\$110	not available	

992 Physician Stool with Chrome 360° Hand Release and Brushed Aluminum Base with Toe Caps



Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$529	\$529	\$609	\$629	\$642	\$688	\$713	\$729	\$736	\$771	\$824
	inch	cm		inch	cm		yd	m			
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92			
Overall Height -low	17.00	43.18	Height Adjustment method	pneumatic			lb	kg			
Overall Height -high	22.00	55.88	Height Release hardware	360° hand			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76		
Optional Footring	18.00	45.72					Shipping Weight ⁸	23	11.34		

Options:

[details & photos on page 34*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		992-FR-xx	\$110	not available	

991 Physician Stool with Chrome 360° Hand Release and Black Nylon Composite Base



Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$442	\$442	\$522	\$542	\$555	\$601	\$626	\$642	\$649	\$684	\$737
	inch	cm		inch	cm		yd	m			
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92			
Overall Height -low	18.00	45.72	Height Adjustment method	pneumatic			lb	kg			
Overall Height -high	23.00	55.88	Height Release hardware	360° hand			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76		
Optional Footring	18.00	45.72					Shipping Weight ⁸	23	11.34		

Options:

[details & photos on page 34*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		991-FR-xx	\$110	not available	

983 Physician Stool with 360° Foot Release and Polished Aluminum Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$744	\$744	\$824	\$844	\$857	\$903	\$928	\$944	\$951	\$986	\$1,039

	inch	cm		inch	cm		yd	m
--	------	----	--	------	----	--	----	---

Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	33.75	85.73	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	38.75	98.43	Height Release hardware	360° foot		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
						Shipping Weight⁸	25	11.34

Adjustable Height Seat ^{††12}	add		inch	cm		inch	cm
5 [five] \$	983-xx	\$0			Seat Height -low	19.75	50.17
7 [seven]	983-7-xx	\$38			Seat Height -low	23.50	59.69
10 [ten]	983-10-xx	\$48			Seat Height -low	25.75	65.41

Options:
[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		not available		Full 983B-xx	\$235
								Crescent 983CB-xx	\$315

[Product Info](#)

982 Physician Stool with 360° Foot Release and Brushed Aluminum Base with Toe Caps



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$747	\$747	\$827	\$847	\$860	\$906	\$931	\$947	\$954	\$989	\$1,042

	inch	cm		inch	cm		yd	m
--	------	----	--	------	----	--	----	---

Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	32.50	82.55	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	37.50	95.25	Height Release hardware	360° foot		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
						Shipping Weight⁸	25	11.34

Adjustable Height Seat ^{††12}	add		inch	cm		inch	cm
5 [five] \$	982-xx	\$0			Seat Height -low	19.00	48.26
7 [seven]	982-7-xx	\$38			Seat Height -low	22.50	59.69
10 [ten]	982-10-xx	\$48			Seat Height -low	24.75	62.87

Options:
[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		not available		Full 982B-xx	\$235
								Crescent 982CB-xx	\$315

[Product Info](#)

981 Physician Stool with 360° Foot Release and Black Nylon Composite Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$675	\$675	\$755	\$775	\$788	\$834	\$859	\$875	\$882	\$917	\$970

	inch	cm		inch	cm		yd	m
--	------	----	--	------	----	--	----	---

Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	33.75	85.73	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	38.75	98.43	Height Release hardware	360° foot		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
						Shipping Weight⁸	25	11.34

Adjustable Height Seat ^{††12}	add		inch	cm		inch	cm
5 [five] \$	981-xx	\$0			Seat Height -low	19.75	50.17
7 [seven]	981-7-xx	\$38			Seat Height -low	23.50	59.69
10 [ten]	981-10-xx	\$48			Seat Height -low	25.75	65.41

Options:
[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		not available		Full 981B-xx	\$235
								Crescent 981CB-xx	\$315

[Product Info](#)

973 Physician Stool with D-Handle Release and Polished Aluminum Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$514	\$514	\$594	\$614	\$627	\$673	\$698	\$714	\$721	\$756	\$809

	inch	cm
Seat Diameter	16.00	40.64
Overall Height with backrest -low	33.50	85.09
Overall Height with backrest -high	38.50	97.79
Optional Footring	18.00	45.72

	inch	cm
Base Diameter with casters	24.00	60.96
Height Adjustment method	pneumatic	
Height Release hardware	D-ring	

	yd	m
COM Requirement¹⁰	1	0.92

	lb	kg
Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Shipping Weight⁸	22	9.07

Adjustable Height Seat ^{††12}	add
5 [five] \$	973-xx \$0
7 [seven]	973-7-xx \$38
10 [ten]	973-10-xx \$48

	inch	cm
Seat Height -low	18.50	46.99
Seat Height -low	20.50	52.07
Seat Height -low	24.50	62.23

	inch	cm
Seat Height -high	24.00	60.96
Seat Height -high	28.25	71.76
Seat Height -high	35.00	88.90

Options:
[details & photos on page 34*](#)

Casters	add
page 34	varies

T-Style Armrest[†]	add
not available	

Procedure Arm^{†‡}	add
973PA-xx	\$450

Footring[†]	add
973-FR-xx	\$110

Backrest^{†‡}	add
Full 973B-xx	\$235
Crescent 973CB-xx	\$315

[Product Info](#)

972 Physician Stool with D-Handle Release and Brushed Aluminum Base with Toe Caps



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$505	\$505	\$585	\$605	\$618	\$664	\$689	\$705	\$712	\$747	\$800

	inch	cm
Seat Diameter	16.00	40.64
Overall Height with backrest -low	32.50	82.55
Overall Height with backrest -high	37.50	95.25
Optional Footring	18.00	45.72

	inch	cm
Base Diameter with casters	24.00	60.96
Height Adjustment method	pneumatic	
Height Release hardware	D-ring	

	yd	m
COM Requirement¹⁰	1	0.92

	lb	kg
Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Shipping Weight⁸	22	9.07

Adjustable Height Seat ^{††12}	add
5 [five] \$	972-xx \$0
7 [seven]	972-7-xx \$38
10 [ten]	972-10-xx \$48

	inch	cm
Seat Height -low	17.50	44.45
Seat Height -low	19.50	49.53
Seat Height -low	23.50	59.69

	inch	cm
Seat Height -high	23.00	58.42
Seat Height -high	27.25	69.22
Seat Height -high	34.00	86.36

Options:
[details & photos on page 34*](#)

Casters	add
page 34	varies

T-Style Armrest[†]	add
not available	

Procedure Arm^{†‡}	add
972PA-xx	\$450

Footring[†]	add
972-FR-xx	\$110

Backrest^{†‡}	add
Full 972B-xx	\$235
Crescent 972CB-xx	\$315

[Product Info](#)

971 Physician Stool with D-Handle Release and Black Nylon Composite Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$432	\$432	\$512	\$532	\$545	\$591	\$616	\$632	\$639	\$674	\$727

	inch	cm
Seat Diameter	16.00	40.64
Overall Height with backrest -low	33.50	85.09
Overall Height with backrest -high	38.50	97.79
Optional Footring	18.00	45.72

	inch	cm
Base Diameter with casters	24.00	60.96
Height Adjustment method	pneumatic	
Height Release hardware	D-ring	

	yd	m
COM Requirement¹⁰	1	0.92

	lb	kg
Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Shipping Weight⁸	22	9.07

Adjustable Height Seat ^{††12}	add
5 [five] \$	971-xx \$0
7 [seven]	971-7-xx \$38
10 [ten]	971-10-xx \$48

	inch	cm
Seat Height -low	18.50	46.99
Seat Height -low	20.50	52.07
Seat Height -low	24.50	62.23

	inch	cm
Seat Height -high	24.00	60.96
Seat Height -high	28.25	71.76
Seat Height -high	35.00	88.90

Options:
[details & photos on page 34*](#)

Casters	add
page 34	varies

T-Style Armrest[†]	add
not available	

Procedure Arm^{†‡}	add
971PA-xx	\$450

Footring[†]	add
971-FR-xx	\$110

Backrest^{†‡}	add
Full 971B-xx	\$235
Crescent 971CB-xx	\$315

[Product Info](#)

963 Physician Stool with Single Lever Release and Polished Aluminum Base



Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$510	\$510	\$590	\$610	\$623	\$669	\$694	\$710	\$717	\$752	\$805
	inch	cm									
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92			
Overall Height with backrest -low	33.50	85.09	Height Adjustment method	pneumatic							
Overall Height with backrest -high	38.50	97.79	Height Release hardware	single lever							
Optional Footring	18.00	45.72									

Options:

[details & photos on page 34*](#)

[Product Info](#)

	add		add		add		add		add
Casters	varies	T-Style Armrest [†]	not available	Procedure Arm ^{†‡}	963PA-xx \$450	Footring [†]	963-FR-xx \$110	Backrest ^{†‡}	Full 963B-xx \$235
page 34									Crescent 963CB-xx \$315

962 Physician Stool with Single Lever Release and Brushed Aluminum Base with Toe Caps



Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$498	\$498	\$578	\$598	\$611	\$657	\$682	\$698	\$705	\$740	\$793
	inch	cm									
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92			
Overall Height with backrest -low	32.50	82.55	Height Adjustment method	pneumatic							
Overall Height with backrest -high	37.50	95.25	Height Release hardware	single lever							
Optional Footring	18.00	45.72									

Options:

[details & photos on page 34*](#)

[Product Info](#)

	add		add		add		add		add
Casters	varies	T-Style Armrest [†]	not available	Procedure Arm ^{†‡}	962PA-xx \$450	Footring [†]	962-FR-xx \$110	Backrest ^{†‡}	Full 962B-xx \$235
page 34									Crescent 962CB-xx \$315

961 Physician Stool with Single Lever Release and Black Nylon Composite Base



Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$412	\$412	\$492	\$512	\$525	\$571	\$596	\$612	\$619	\$654	\$707
	inch	cm									
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92			
Overall Height with backrest -low	33.50	85.09	Height Adjustment method	pneumatic							
Overall Height with backrest -high	39.00	99.06	Height Release hardware	single lever							
Optional Footring	18.00	45.72									

Options:

[details & photos on page 34*](#)

[Product Info](#)

	add		add		add		add		add
Casters	varies	T-Style Armrest [†]	not available	Procedure Arm ^{†‡}	961PA-xx \$450	Footring [†]	961-FR-xx \$110	Backrest ^{†‡}	Full 961B-xx \$235
page 34									Crescent 961CB-xx \$315

958 Industrial Stool with Single Lever Release, Footring, 5 Leg Black Tubular Steel Base, and Glides
**Specs:**

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$634	\$634	\$714	\$734	\$747	\$793	\$818	\$834	\$841	\$876	\$929

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with glides	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	33.00	83.82	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	39.00	99.06	Height Release hardware	single lever		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Footring	20.00	50.80				Shipping Weight⁸	23	10.43

Adjustable Height Seat ^{††12}		add		inch	cm		inch	cm
5 [five] \$	958-xx	\$0	Seat Height-low	17.50	44.45	Seat Height-high	23.00	58.42
7 [seven]	958-7-xx	\$38	Seat Height-low	19.50	49.53	Seat Height-high	26.50	67.31
10 [ten]	958-10-xx	\$48	Seat Height-low	23.50	59.69	Seat Height-high	34.00	86.36

958 Options:
[details & photos on page 34*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm ^{†‡}	add	Footring [†]	add	Backrest ^{†‡}	add
not available		not available		958PA-xx	\$450	standard	\$0	Full 958B-xx	\$235
								Crescent 958CB-xx	\$315

951 Physician Stool with Single Lever Release and Black Nylon Composite Base
**Specs:**

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$411	\$411	\$491	\$511	\$524	\$570	\$595	\$611	\$618	\$653	\$706

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	33.00	83.82	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	39.00	99.06	Height Release hardware	single lever		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Optional Footring	18.00	45.72				Shipping Weight⁸	20	9.07

Adjustable Height Seat ^{††12}			add		inch	cm		inch	cm
5 [five] ^S	951-xx	\$0	Seat Height -low		18.50	46.99	Seat Height -high	24.00	60.98
7 [seven]	951-7-xx	\$38	Seat Height -low		20.50	52.07	Seat Height -high	28.25	71.76
10 [ten]	951-10-xx	\$48	Seat Height -low		24.50	62.23	Seat Height -high	35.00	88.90

951 Options:
[details & photos on page 34*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm ^{†‡}	add	Footring [†]	add	Backrest ^{†‡}	add
page 34	varies	not available		951PA-xx	\$450	951-FR-xx	\$110	Full 951B-xx	\$235
								Crescent 951CB-xx	\$315

928 Urethane Seat Industrial Stool with Footring and Black Tubular 5-Leg Steel Base with Glides
**Specs:**

list	
\$591	 available in black urethane only

	inch	cm		inch	cm		lb	kg
Seat Diameter	14.00	35.56	Base Diameter with glides	24.00	60.96			
Overall Height -low	15.00	38.10	Height Adjustment method	pneumatic		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Overall Height -high	19.75	50.17	Height Release hardware	single lever		Shipping Weight⁸	23	10.43
Footring	20.00	50.80						

Adjustable Height Seat ^{††12}		add		inch	cm		inch	cm
5 [five] [§]	928	\$0	Seat Height ^{-low}	15.00	38.10	Seat Height ^{-high}	19.75	50.17
7 [seven]	928-7	\$38	Seat Height ^{-low}	16.50	41.91	Seat Height ^{-high}	23.25	59.06
10 [ten]	928-10	\$48	Seat Height ^{-low}	20.50	52.07	Seat Height ^{-high}	31.00	78.74

928 Options:
[details & photos on page 34*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm ^{†‡}	add	Footring [†]	add	Backrest ^{†‡}	add
not available		not available		928PA-xx	\$450	standard	\$0	not available	

878 Lab Stool with D-Handle Release, Footring, and Black Tubular 5-Leg Steel Base with Glides



 Product Info

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$881	\$881	\$961	\$981	\$994	\$1,040	\$1,065	\$1,081	\$1,088	\$1,123	\$1,176
†	inch		cm	inch		cm	yd		m	
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰		1	0.92	
Overall Height with backrest -low	38.00	96.52	Backrest Height	15.75	45.72	Height Release hardware		D-ring		
Overall Height with backrest -high	48.00	121.92	Base Diameter with glides	24.00	60.96			lb	kg	
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76	
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸		32	14.51	

Options:



Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	878-7-xx	\$0		19.25	48.90		26.50	67.31	
10 [ten] [§]	878-xx standard	\$0		23.50	59.69		34.00	86.36	

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
not available		not available		not available		standard	\$0	Crescent 878CB-xx	\$315

873 Lab Stool with D-Handle Release, Footring, and Polished Aluminum Base



 Product Info


Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10	
\$856	\$856	\$936	\$956	\$969	\$1,015	\$1,040	\$1,056	\$1,063	\$1,098	\$1,151	
†											
		inch	cm			inch	cm			yd	m
Seat Diameter		16.00	40.64	Backrest Width		9.75	24.77	COM Requirement ¹⁰		1	0.92
Overall Height with backrest -low		37.00	93.98	Backrest Height		15.75	45.72	Height Release hardware		D-ring	
Overall Height with backrest -high		47.00	119.38	Base Diameter with glides		24.00	60.96				
Seat Foam thickness		3.00	7.62	Height Adjustment method		pneumatic				lb	kg
Footring		20.00	50.80	Backrest Adjustment direction		lateral		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
								Shipping Weight ⁸		28	12.70

Options:



Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	873-7-xx	\$0		20.50	52.07		27.50	69.85	
10 [ten] [§]	873-xx standard	\$0		24.50	62.23		35.00	88.90	

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
 varies		not available		not available		standard	\$0	Crescent 873CB-xx	\$315

872 Lab Stool with D-Handle Release, Footring, and Brushed Aluminum Base with Toe Caps



 Product Info


Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10	
\$890	\$890	\$970	\$990	\$1,003	\$1,049	\$1,074	\$1,090	\$1,097	\$1,132	\$1,185	
		inch	cm			inch	cm			yd	m
Seat Diameter		16.00	40.64	Backrest Width		9.75	24.77	COM Requirement ¹⁰		1	0.92
Overall Height <i>with backrest -low</i>		36.00	91.44	Backrest Height		15.75	45.72	Height Release <i>hardware</i>		D-ring	
Overall Height <i>with backrest -high</i>		46.00	116.84	Base Diameter <i>with glides</i>		24.00	60.96			lb	kg
Seat Foam <i>thickness</i>		3.00	7.62	Height Adjustment <i>method</i>		pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
Footring		20.00	50.80	Backrest Adjustment <i>direction</i>		lateral		Shipping Weight ⁸		28	12.70

Options:



Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	872-7-xx	\$0		20.50	52.07		27.50	69.85	
10 [ten] [§]	872-xx standard	\$0		23.50	59.69		34.00	86.36	

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
 varies		not available		not available		standard	\$0	Crescent 872CB-xx	\$315

871 Lab Stool with D-Handle Release, Footring, and Black Nylon Composite Base

 Product Info

 details & photos on page 34*
Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$810	\$810	\$890	\$910	\$923	\$969	\$994	\$1,010	\$1,017	\$1,052	\$1,105
	inch		cm		inch		cm		yd		m
Seat Diameter	16.00		40.64		9.75		24.77		1		0.92
Overall Height with backrest -low	37.00		93.98		15.75		45.72		D-ring		
Overall Height with backrest -high	47.00		119.38		24.00		60.96		lb		kg
Seat Foam thickness	3.00		7.62		pneumatic				350		158.76
Footring	20.00		50.80		lateral				27		12.25

Options:

Adjustable Height Seat ^{††12}			add		inch	cm		inch	cm
7 [seven]	871-7-xx		\$0	Seat Height -low	19.50	49.53	Seat Height -high	27.50	69.85
10 [ten]	871-xx	standard	\$0	Seat Height -low	24.50	62.23	Seat Height -high	35.00	88.90
Casters	add			T-Style Armrest [†]	add		Procedure Arm [†]	add	
 page 34	varies			not available			not available		
							Footring [†]	add	
							standard	\$0	
							Backrest [†]	add	
							Crescent 871CB-xx		\$315

860 SERIES

Laboratory Seating

868 Lab Stool with Single Lever Release, Footring, 5 Leg Black Tubular Steel Base, and Glides

 Product Info


 details & photos on page 34*
Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$825	\$825	\$905	\$925	\$938	\$984	\$1,009	\$1,025	\$1,032	\$1,067	\$1,120
	inch		cm		inch		cm		yd		m
Seat Diameter	18.00		45.72		9.75		24.77		1		0.92
Overall Height with backrest -low	38.00		96.52		15.75		45.72		single lever		
Overall Height with backrest -high	48.00		121.92		24.00		60.96		lb		kg
Seat Foam thickness	3.00		7.62		pneumatic				350		158.76
Footring	20.00		50.80		lateral				32		14.51

Options:

Adjustable Height Seat ^{††12}			add		inch	cm		inch	cm
7 [seven]	868-7-xx		\$0	Seat Height -low	19.50	49.53	Seat Height -high	27.50	69.85
10 [ten]	868-xx	standard	\$0	Seat Height -low	23.50	59.69	Seat Height -high	34.00	86.36
Casters	add			T-Style Armrest [†]	add		Procedure Arm [†]	add	
not available				not available			not available		
							Footring [†]	add	
							standard	\$0	
							Backrest [†]	add	
							Crescent 868CB-xx		\$315

CB | optional crescent backrest★ | standard full backrest



863 Lab Stool with Single Lever Release, Footring, and Polished Aluminum Base


[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$852	\$852	\$932	\$952	\$965	\$1,011	\$1,036	\$1,052	\$1,059	\$1,094	\$1,147

	inch	cm		inch	cm		yd	m
Specs: Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	37.00	93.98	Backrest Height	15.75	45.72	Height Release hardware	single lever	
Overall Height with backrest -high	47.00	119.38	Base Diameter with glides	24.00	60.96		lb	kg
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	18.00	45.72	Backrest Adjustment direction	lateral		Shipping Weight ⁸	32	14.51

Options:

[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		standard	\$0	Crescent 863CB-xx	\$315

862 Lab Stool with Single Lever Release, Footring, and Brushed Aluminum Base with Toe Caps


[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$882	\$882	\$962	\$982	\$995	\$1,041	\$1,066	\$1,082	\$1,089	\$1,124	\$1,177

	inch	cm		inch	cm		yd	m
Specs: Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	36.00	91.44	Backrest Height	15.75	45.72	Height Release hardware	single lever	
Overall Height with backrest -high	46.00	116.84	Base Diameter with glides	24.00	60.96		lb	kg
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	18.00	45.72	Backrest Adjustment direction	lateral		Shipping Weight ⁸	28	12.70

Options:

[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
 page 34	varies	not available		not available		standard	\$0	Crescent 862CB-xx	\$315

861 Lab Stool with Single Lever Release, Footring, and Black Nylon Composite Base


[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$802	\$802	\$882	\$902	\$915	\$961	\$986	\$1,002	\$1,009	\$1,044	\$1,097

	inch	cm		inch	cm		yd	m
Specs: Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	37.00	93.98	Backrest Height	15.75	45.72	Height Release hardware	single lever	
Overall Height with backrest -high	47.00	119.38	Base Diameter with glides	24.00	60.96		lb	kg
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	18.00	45.72	Backrest Adjustment direction	lateral		Shipping Weight ⁸	27	12.25

Options:

[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	not available		not available		standard	\$0	Crescent 861CB-xx	\$315

848DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Black Tubular 5-Leg Steel Base with Glides



[Product Info](#)

[details & photos on page 34*](#)

848DC Options:

Specs:	list		\$758		 available in black urethane only	
	†	inch	cm	inch	cm	
	Seat Width	18.00	45.72	Backrest Width	16.00	40.64
	Seat Depth	17.00	43.18	Backrest Height	12.00	30.48
	Overall Height with backrest -low	37.00	93.98	Base Diameter with casters	24.00	60.96
Overall Height with backrest -high		47.50	120.65	Height Adjustment method	pneumatic	
Footring		20.00	50.80	Backrest Adjustment direction	height	
Adjustable Height Seat ^{†12}		add				
7 [seven]		848DC-7	\$0	Seat Height-low	17.50	44.45
10 [ten] [§]		848DC standard	\$0	Seat Height-low	21.00	53.34
Casters		add				
not available						
T-Style Armrest		add				
848DC-T		\$187				
Procedure Arm		add				
not available						
Footring		add				
standard		\$0				
Backrest		add				
Full standard		\$0				

843DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Polished Aluminum Base



[Product Info](#)

[details & photos on page 34*](#)

843DC Options:

Specs:	list		\$755		 available in black urethane only	
	†	inch	cm	inch	cm	
	Seat Width	18.00	45.72	Backrest Width	16.00	40.64
	Seat Depth	17.00	43.18	Backrest Height	12.00	30.48
	Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96
Overall Height with backrest -high		49.00	124.46	Height Adjustment method	pneumatic	
Footring		20.00	50.80	Backrest Adjustment direction	height	
Adjustable Height Seat ^{†12}		add				
7 [seven]		843DC-7	\$0	Seat Height-low	18.50	46.99
10 [ten] [§]		843DC standard	\$0	Seat Height-low	22.50	57.15
Casters		add				
page 34		varies				
T-Style Armrest		add				
843DC-T		\$187				
Procedure Arm		add				
not available						
Footring		add				
standard		\$0				
Backrest		add				
Full standard		\$0				

842DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Brushed Aluminum Base with Toe Caps



[Product Info](#)

[details & photos on page 34*](#)

842DC Options:

Specs:	list		\$789		 available in black urethane only	
	†	inch	cm	inch	cm	
	Seat Width	18.00	45.72	Backrest Width	16.00	40.64
	Seat Depth	17.00	43.18	Backrest Height	12.00	30.48
	Overall Height with backrest -low	37.50	97.79	Base Diameter with casters	24.00	60.96
Overall Height with backrest -high		48.00	124.46	Height Adjustment method	pneumatic	
Footring		20.00	50.80	Backrest Adjustment direction	height	
Adjustable Height Seat ^{†12}		add				
7 [seven]		842DC-7	\$0	Seat Height-low	17.50	44.45
10 [ten] [§]		842DC standard	\$0	Seat Height-low	21.50	54.61
Casters		add				
page 34		varies				
T-Style Armrest		add				
842DC-T		\$187				
Procedure Arm		add				
not available						
Footring		add				
standard		\$0				
Backrest		add				
Full standard		\$0				

841DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Black Nylon Composite Base



[Product Info](#)

[details & photos on page 34*](#)

list
\$748

[i](#) available in black urethane only

Specs:

	inch	cm		inch	cm		
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96	Height Adjustment seat & back tilt	dual lever
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic			
Footring	20.00	50.80	Backrest Adjustment direction	height			
						Weight Capacity ¹ EVENLY DISTRIBUTED	350 158.76
						Shipping Weight ⁸	43 19.50

841DC Options:

Adjustable Height Seat ^{†12}	add	inch	cm	inch	cm
7 [seven]	841DC-7	\$0		Seat Height-low	18.50 46.99
10 [ten] [§]	841DC standard	\$0		Seat Height-low	22.50 57.15
				Seat Height-high	25.50 57.15
				Seat Height-high	33.00 83.82

Casters	add	T-Style Armrest	add	Procedure Arm	add	Footring	add	Backrest	add
page 34	varies	841DC-T	\$187	not available		standard	\$0	Full standard	\$0

848SC Self Skin Ergonomic Chair with Footring and Black Tubular 5-Leg Steel Base with Glides



[Product Info](#)

[details & photos on page 34*](#)

list
\$693

[i](#) available in black urethane only

Specs:

	inch	cm		inch	cm		
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane
Overall Height with backrest -low	37.50	95.25	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever
Overall Height with backrest -high	47.50	120.65	Height Adjustment method	pneumatic			
Footring	20.00	50.80	Backrest Adjustment direction	lateral			
						Weight Capacity ¹ EVENLY DISTRIBUTED	350 158.76
						Shipping Weight ⁸	39 17.69

848SC Options:

Adjustable Height Seat ^{†12}	add	inch	cm	inch	cm
7 [seven]	848DC-7	\$0		Seat Height-low	17.50 44.45
10 [ten] [§]	848DC standard	\$0		Seat Height-low	21.00 53.34
				Seat Height-high	24.50 62.23
				Seat Height-high	31.50 80.01

Casters	add	T-Style Armrest	add	Procedure Arm	add	Footring	add	Backrest	add
not available		848SC-T	\$187	not available		standard	\$0	Full standard	\$0

843SC Self Skin Ergonomic Chair with Footring and Polished Aluminum Base



[Product Info](#)

[details & photos on page 34*](#)

list
\$690

[i](#) available in black urethane only

Specs:

	inch	cm		inch	cm		
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic			
Footring	20.00	50.80	Backrest Adjustment direction	lateral			
						Weight Capacity ¹ EVENLY DISTRIBUTED	350 158.76
						Shipping Weight ⁸	39 17.69

843SC Options:

Adjustable Height Seat ^{†12}	add	inch	cm	inch	cm
7 [seven]	843SC-7	\$0		Seat Height-low	18.50 46.99
10 [ten] [§]	843SC standard	\$0		Seat Height-low	22.50 57.15
				Seat Height-high	25.50 57.15
				Seat Height-high	33.00 83.82

Casters	add	T-Style Armrest	add	Procedure Arm	add	Footring	add	Backrest	add
page 34	varies	843SC-T	\$187	not available		standard	\$0	Full standard	\$0

842SC Self Skin Ergonomic Chair with Footring and Brushed Aluminum Base with Toe Caps


[Product Info](#)

 list
\$725

available in black urethane only

Specs:

	inch	cm		inch	cm		
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane
Overall Height with backrest -low	37.50	95.25	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever
Overall Height with backrest -high	48.00	121.92	Height Adjustment method	pneumatic			
Footring	20.00	50.80	Backrest Adjustment direction	lateral			
						Weight Capacity ¹ EVENLY DISTRIBUTED	350 158.76
						Shipping Weight ⁸	39 17.69

842SC Options:

[details & photos on page 34*](#)

Adjustable Height Seat ^{†12}		add		inch	cm		inch	cm	
7 [seven]	842SC-7	\$0	Seat Height-low	17.50	44.45	Seat Height-high	24.50	62.23	
10 [ten] [§]	842SC <i>standard</i>	\$0	Seat Height-low	21.50	54.61	Seat Height-high	32.00	81.28	
T-Style Armrest	add	Procedure Arm	add	Footring		add	Backrest		add
842SC-T	\$187	not available		standard		\$0	Full standard		\$0

841SC Self Skin Ergonomic Chair with Footring, and Black Nylon Composite Base


[Product Info](#)

 list
\$616

available in black urethane only

Specs:


	inch	cm		inch	cm		
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic			
Footring	20.00	50.80	Backrest Adjustment direction	lateral			
						Weight Capacity ¹ EVENLY DISTRIBUTED	350 158.76
						Shipping Weight ⁸	38 17.24

841SC Options:

[details & photos on page 34*](#)

Adjustable Height Seat ^{†12}		add		inch	cm		inch	cm
7 [seven]	841SC-7	\$0	Seat Height-low	18.50	46.99	Seat Height-high	25.50	57.15
10 [ten] [§]	841SC standard	\$0	Seat Height-low	22.50	57.15	Seat Height-high	33.00	83.82
T-Style Armrest	add	Procedure Arm	add	Footring	add	Backrest	add	
841SC-T	\$187	not available		standard	\$0	Full standard		\$0

T | adjustable urethane t-style armrests option★



838 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Black Tubular 5-Leg Steel Base with Glides


 Product Info

 details & photos on page 34*

Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$1,075	\$1,075	\$1,155	\$1,175	\$1,188	\$1,234	\$1,259	\$1,275	\$1,282	\$1,317	\$1,370
	inch	cm		inch	cm			yd	m		
Seat Width	20.00	50.80	Backrest Width	19.00	48.26	COM Requirement ¹⁰		1	0.92		
Seat Depth	18.00	45.72	Backrest Height	16.00	40.64	Height Release hardware		dual lever			
Overall Height with backrest -low	39.00	99.06	Base Diameter with casters	24.00	60.96			lb	kg		
Overall Height with backrest -high	49.50	125.73	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76		
Footring	20.00	50.80	Backrest Adjustment direction	height		Shipping Weight ⁸		45	20.41		

Options:

Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	838-7-xx	\$0		19.50	49.53		26.50	67.31	
10 [ten]	838-xx standard	\$0		22.50	57.15		33.00	83.82	
Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
not available		838-T-xx	\$187	not available		standard	\$0	Full standard	\$0

833 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Polished Aluminum Base


 Product Info

 details & photos on page 34*

Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$1,065	\$1,065	\$1,145	\$1,165	\$1,178	\$1,224	\$1,249	\$1,265	\$1,272	\$1,307	\$1,360
	inch	cm		inch	cm			yd	m		
Seat Width	20.00	50.80	Backrest Width	19.00	48.26	COM Requirement ¹⁰		1	0.92		
Seat Depth	18.00	45.72	Backrest Height	16.00	40.64	Height Release hardware		dual lever			
Overall Height with backrest -low	39.00	99.06	Base Diameter with casters	24.00	60.96			lb	kg		
Overall Height with backrest -high	49.50	125.73	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76		
Footring	20.00	50.80	Backrest Adjustment direction	height		Shipping Weight ⁸		45	20.41		

Options:

Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	833-7-xx	\$0		20.00	50.80		27.00	68.58	
10 [ten]	833-xx standard	\$0		24.00	60.96		34.50	87.63	
Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	833-T-xx	\$187	not available		standard	\$0	Full standard	\$0

832 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Brushed Aluminum Base with Toe Caps


 Product Info

 details & photos on page 34*

Specs:

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$1,054	\$1,054	\$1,134	\$1,154	\$1,167	\$1,213	\$1,238	\$1,254	\$1,261	\$1,296	\$1,349
	inch	cm		inch	cm			yd	m		
Seat Width	20.00	50.80	Backrest Width	19.00	48.26	COM Requirement ¹⁰		1	0.92		
Seat Depth	18.00	45.72	Backrest Height	16.00	40.64	Height Release hardware		dual lever			
Overall Height with backrest -low	38.00	96.52	Base Diameter with casters	24.00	60.96			lb	kg		
Overall Height with backrest -high	48.50	123.19	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76		
Footring	20.00	50.80	Backrest Adjustment direction	height		Shipping Weight ⁸		45	20.41		

Options:

Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	832-7-xx	\$0		19.00	48.26		26.00	66.04	
10 [ten]	832-xx standard	\$0		23.00	58.42		33.50	85.09	
Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	832-T-xx	\$187	not available		standard	\$0	Full standard	\$0

831 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Black Nylon Composite Base
**Specs:**

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$950	\$950	\$1,030	\$1,050	\$1,063	\$1,109	\$1,134	\$1,150	\$1,157	\$1,192	\$1,245
		inch	cm			inch	cm			yd	m
Seat Width		20.00	50.80	Backrest Width		19.00	48.26	COM Requirement ¹⁰		1	0.92
Seat Depth		18.00	45.72	Backrest Height		16.00	40.64	Height Release hardware		dual lever	
Overall Height with backrest -low		39.00	99.06	Base Diameter with casters		24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		lb	kg
Overall Height with backrest -high		49.50	125.73	Height Adjustment method		pneumatic		Shipping Weight ⁸		45	20.41
Footring		20.00	50.80	Backrest Adjustment direction		height					

Options:

Adjustable Height Seat ^{††12}			add				inch	cm			inch	cm
7 [seven]	831-7-xx	\$0	Seat Height-low		20.00	50.80	Seat Height-high		27.00	68.58		
10 [ten] [§]	831-xx standard	\$0	Seat Height-low		24.00	60.96	Seat Height-high		34.50	87.63		
T-Style Armrest [†]	add	Procedure Arm [†]		add	Footring [†]		add	Backrest [†]		add		
831-T-xx	\$187	not available			standard		\$0	Full standard		\$0		

 Product Info

[details & photos on page 34★](#)

810 SERIES

Laboratory Seating

818 Ergonomic Lab Chair with Footring, and Black Tubular 5-Leg Steel Base with Glides
**Specs:**

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$944	\$944	\$1,024	\$1,044	\$1,057	\$1,103	\$1,128	\$1,144	\$1,151	\$1,186	\$1,239
		inch	cm			inch	cm			yd	m
Seat Width		20.00	50.80	Backrest Width		16.00	40.64	COM Requirement ¹⁰		1	0.92
Seat Depth		18.00	45.72	Backrest Height		18.00	45.72	Height Release hardware		single lever	
Overall Height with backrest -low		41.50	105.41	Base Diameter with casters		24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		lb	kg
Overall Height with backrest -high		52.00	132.08	Height Adjustment method		pneumatic		Shipping Weight ⁸		45	20.41
Footring		20.00	50.80	Backrest Adjustment direction		lateral					

Options:

Adjustable Height Seat ^{††12}		add			inch	cm			inch	cm
7 [seven]	818-7-xx	\$0	Seat Height-low		19.50	49.53	Seat Height-high		26.50	67.31
10 [ten] [§]	818-xx standard	\$0	Seat Height-low		22.50	57.15	Seat Height-high		33.00	83.82
T-Style Armrest [†]	add	Procedure Arm [†]		add	Footring [†]		add	Backrest [†]		add
818-T-xx	\$187	not available			standard		\$0	Full standard		\$0

 Product Info

[details & photos on page 34★](#)
813 Ergonomic Lab Chair with Footring, and Polished Aluminum Base
**Specs:**

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$924	\$924	\$1,004	\$1,024	\$1,037	\$1,083	\$1,108	\$1,124	\$1,131	\$1,166	\$1,219
		inch	cm			inch	cm			yd	m
Seat Width		20.00	50.80	Backrest Width		16.00	40.64	COM Requirement ¹⁰		1	0.92
Seat Depth		18.00	45.72	Backrest Height		18.00	45.72	Height Release hardware		single lever	
Overall Height with backrest -low		42.00	106.68	Base Diameter with casters		24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		lb	kg
Overall Height with backrest -high		52.50	133.35	Height Adjustment method		pneumatic		Shipping Weight ⁸		43	19.50
Footring		20.00	50.80	Backrest Adjustment direction		lateral					

Options:

Adjustable Height Seat ^{††12}		add		inch	cm		inch	cm
7 [seven] [§]	813-7-xx	\$0	Seat Height - low	20.00	50.80	Seat Height -high	27.00	68.58
10 [ten]	813-xx standard	\$0	Seat Height -low	24.00	60.96	Seat Height -high	34.50	87.63
T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add	
813-T-xx	\$187	not available		standard	\$0	Full standard	\$0	

 Product Info

[details & photos on page 34★](#)

812 Ergonomic Lab Chair with Footring and Brushed Aluminum Base with Toe Caps



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$938	\$938	\$1,018	\$1,038	\$1,051	\$1,097	\$1,122	\$1,138	\$1,145	\$1,180	\$1,233

Specs: †

	inch	cm		inch	cm		yd	m
Seat Width	20.00	50.80	Backrest Width	16.00	40.64	COM Requirement ¹⁰	1	0.92
Seat Depth	18.00	45.72	Backrest Height	18.00	45.72	Height Release hardware	single lever	
Overall Height with backrest -low	41.00	104.14	Base Diameter with casters	24.00	60.96		lb	kg
Overall Height with backrest -high	51.50	130.81	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸	43	19.50

Options: 812

[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	812-T-xx	\$187	not available		standard	\$0	Full standard	\$0

[Product Info](#)

811 Ergonomic Lab Chair with Footring, and Black Nylon Composite Base



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$832	\$832	\$912	\$932	\$945	\$991	\$1,016	\$1,032	\$1,039	\$1,074	\$1,127

Specs: †

	inch	cm		inch	cm		yd	m
Seat Width	20.00	50.80	Backrest Width	16.00	40.64	COM Requirement ¹⁰	1	0.92
Seat Depth	18.00	45.72	Backrest Height	18.00	45.72	Height Release hardware	single lever	
Overall Height with backrest -low	42.00	106.68	Base Diameter with casters	24.00	60.96		lb	kg
Overall Height with backrest -high	52.50	133.35	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸	42	19.05

Options: 811

[details & photos on page 34*](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 34	varies	811-T-xx	\$187	not available		standard	\$0	Full standard	\$0

[Product Info](#)



xx1 | black nylon composite base



xx2 | brushed aluminum base with toe caps

xx3 | polished aluminum base

xx8 | black tubular five-leg steel base and glides



748 Hydraulic Surgical Stool with Silver Column, Foot Release and Polished Aluminum Base


[Product Info](#)
[details & photos on page 34*](#)

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$1,178	\$1,178	\$1,258	\$1,278	\$1,291	\$1,337	\$1,362	\$1,378	\$1,385	\$1,420	\$1,473
	inch	cm			inch	cm				yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters		24.00	60.96	COM Requirement ¹⁰		1	0.92	
Overall Height - low	35.50	90.17	Height Adjustment method		hydraulic				lb	kg	
Overall Height - high	42.50	107.95	Rotation Lock		on seat		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76	
Height Release hardware	foot control		Powder Coat Finish		silver		Shipping Weight ⁸		45	20.41	

Options: 748

Casters	add	T-Style Armrest [‡]	add	Procedure Arm [‡]	add	Footring [‡]	add	Backrest [‡]	add
 page 34	varies	not available		not available		not available		Full 748B-xx	\$235
								Crescent 748CB-xx	\$315

744 Hydraulic Surgical Stool with Black Column, Foot Release and Black Aluminum Base


[Product Info](#)
[details & photos on page 34*](#)

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$1,178	\$1,178	\$1,258	\$1,278	\$1,291	\$1,337	\$1,362	\$1,378	\$1,385	\$1,420	\$1,473
	inch	cm			inch	cm				yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters		24.00	60.96	COM Requirement ¹⁰		1	0.92	
Overall Height - low	35.50	90.17	Height Adjustment method		hydraulic				lb	kg	
Overall Height - high	42.50	107.95	Rotation Lock		on seat		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76	
Height Release hardware	foot control		Powder Coat Finish		black		Shipping Weight ⁸		45	20.41	

Options: 744

Casters	add	T-Style Armrest [‡]	add	Procedure Arm [‡]	add	Footring [‡]	add	Backrest [‡]	add
 page 34	varies	not available		not available		not available		Full 744B-xx	\$235
								Crescent 744CB-xx	\$315

HSS500 Surgeons Stool with Foot-Controlled Hydraulics, 3-Spindle Base, and 4 Dual Tente® Casters


[Product Info](#)


Dual-density Foam Seat with Dura-Gel® Insert

Specs: †

list											
\$6,420											
	inch	cm		inch	cm					degrees	
Overall Width ³	25.50	64.77	Base Width ³	24.00	60.96	Armrest Rotation		360			
Overall Depth ³	32.25	81.92	Base Depth ³	28.50	72.39	Armrest Support		180			
Overall Height - low	21.00	53.34	Seat Height - low	21.50	54.61	Backrest Tilt		20			
Overall Height - high	28.00	71.12	Seat Height - high	28.50	72.39						
Armrest Width ³	4.75	12.07	Backrest Width	14.00	35.56			lb	kg		
Armrest Length ³	12.00	30.48	Backrest Length	8.25	20.96	Weight Capacity ¹ EVENLY DISTRIBUTED		300	136.08		
Armrest Height - low	7.00	17.78	Backrest Height - low	10.50	26.67	Shipping Weight ⁸		74	33.57		
Armrest Height - high	15.00	38.10	Backrest Height - high	12.50	31.75						
			Tente® Casters four sets - dual	3.00	7.62						

HSS500 Options:

Saddle Seat [§]			add	Wedge Seat			add	Water Drop Seat			add	Reverse Water Drop			add
HSS500-S			\$0	HSS500-W			\$0	HSS500-WD			\$0	HSS500-RWD			\$0
		inch	cm			inch	cm			inch	cm			inch	cm
	Width ³	17.00	43.18		Width ³	13.50	34.29		Width ³	15.75	40.00		Width ³	15.75	40.00
	Depth ³	17.65	44.83		Depth ³	18.00	45.72		Depth ³	18.00	45.72		Depth ³	18.00	45.72

³ variance: ± 0.5 inch

IN601P Pivot-Arm Clinical Care Recliner with 3 Positions, Trendelenburg, and Upholstered Arms

Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10	
\$4,308	\$4,308	\$4,787	\$4,907	\$4,987	\$5,263	\$5,414	\$5,507	\$5,548	\$5,758	\$6,078	
†		inch	cm			inch	cm			yd	m
Overall Width	34.50	87.63	Seat Width		22.50	57.15	COM Requirement ¹⁰		7	6.40	
Overall Width <i>side table up</i>	54.00	137.16	Seat Depth		19.00	48.26			inch	cm†	
Overall Depth	38.00	96.52	Seat Height		22.00	55.88	Fold Down Side Tables		10 x 19	25 x 48	
Overall Depth <i>trendelenburg</i>	72.00	182.88	Arm Width		3.00	7.62					
Overall Height	51.50	130.81	Width <i>between arms</i>		25.00	63.50			lb	kg	
Back Width	30.50	77.47	Arm Height <i>from floor</i>		28.50	72.39	Weight Capacity ¹ EVENLY DISTRIBUTED		400	181.44	
Back Height <i>above seat</i>	30.00	76.20	Tente® Casters		5.00	12.70	Shipping Weight ⁸		197	89.36	

IN601P Options & Accessories:

details & photos on page 35*

Product Info

Heat & Massage†			add	IV Mounting			add	Headrest†			add	Extended†			add	Folding Tables			add
None	IN601Pxx	\$0		Mount	FR601IVMNT ⁵	\$64		Skirt	HRSxx	\$78		Legrest	IN601PLxx ^{6S}	\$0		Fold Down Side	\$	\$0	
Both	IN601PHxx	\$1,740		Pole	2101A	\$64		Cover	HRC601xx	\$149		Legs	not available			Swing-Away	SAFT20-601	\$447	
Heat Only	IN601PHOxx	\$1,494						Pillow	FR5659xxUS	\$98									

IN601U Pivot-Arm Clinical Care Recliner with 3 Positions, Trendelenburg, and Urethane Capped Arms

Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$4,480	\$4,480	\$4,959	\$5,079	\$5,159	\$5,435	\$5,586	\$5,679	\$5,720	\$5,930	\$6,250
†	inch	cm			inch	cm			yd	m
Overall Width	34.50	87.63	Seat Width		22.50	57.15	COM Requirement ¹		7	6.40
Overall Width <i>side table up</i>	54.00	137.16	Seat Depth		19.00	48.26			inch	cm†
Overall Depth	38.00	96.52	Seat Height		22.00	55.88	Fold Down Side Tables		10 x 19	25 x 48
Overall Depth <i>trendelenburg</i>	72.00	182.88	Arm Width		3.00	7.62				
Overall Height	51.50	130.81	Width <i>between arms</i>		25.00	63.50			lb	kg
Back Width	30.50	77.47	Arm Height <i>from floor</i>		28.50	72.39	Weight Capacity ¹⁰ EVENLY DISTRIBUTED		400	181.44
Back Height <i>above seat</i>	30.00	76.20	Tente® Casters		5.00	12.70	Shipping Weight ⁸		197	89.36

IN601U Options & Accessories:

details & photos on page 35*

Product Info

Heat & Massage†			add	IV Mounting			add	Headrest†			add	Extended†			add	Folding Table			add
None	IN601Uxx	\$0		Mount	FR601IVMNT ⁵	\$64		Skirt	HRSxx	\$78		Legrest	IN601ULxx ⁶	\$0		Fold Down Side	\$	\$0	
Both	IN601UHxx	\$1,740		Pole	2101A	\$64		Cover	HRC601xx	\$149		Legs	not available			Swing-Away	SAFT20-601	\$447	
Heat Only	IN601UHOxx	\$1,494						Pillow	FR5659xxUS	\$98									



IN597G Ortho-Biotic II Recliner with Infinite Back Positioning, Trendelenburg, and Wood Capped Arms



Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$4,788	\$4,788	\$5,267	\$5,387	\$5,467	\$5,743	\$5,894	\$5,987	\$6,028	\$6,238	\$6,558
†	inch	cm			inch	cm			yd	m
Overall Width	33.00	83.82	Seat Width		23.00	58.42	COM Requirement ¹⁰		6	5.49
Overall Width <i>side table up</i>	45.00	114.30	Seat Depth		20.00	50.80			inch	cm†
Overall Depth	37.00	93.98	Seat Height		21.00	53.34	Fold Down Side Table ⁴⁵		10 x 19	25 x 48
Overall Depth <i>trendelenburg</i>	73.00	185.42	Arm Width		3.50	8.89				
Overall Height	46.00	116.84	Width <i>between arms</i>		24.00	60.96			lb	kg
Back Width	23.00	58.42	Arm Height <i>from floor</i>		26.75	67.95	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
Back Height <i>above seat</i>	30.00	76.20	Tente® Casters		4.00	10.16	Shipping Weight ⁸		210	95.25

IN597G Options & Accessories:

details & photos on page 35*

Product Info

Heat & Massage [†]	add	IV Mounting	add	Headrest [†]	add	Extended	Folding Table	add
None IN597Gxx	\$0	Mount FR597GIVMNT ⁵	\$48	Skirt HRSxx	\$78	Legrest not available	Left Side Single FR597GTBLMNT ⁴⁵	\$102
Both IN597GHxx	\$1,852	Pole 2101A	\$64	Cover not available		Legs not available	Swing-Away SAFT20-597	\$447
Heat Only IN597GHOxx	\$1,604			Pillow standard [§]	\$0			

IN597P Ortho-Biotic II Recliner with Infinite Back Positioning, Trendelenburg, and Upholstered Arms



Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$4,646	\$4,646	\$5,125	\$5,245	\$5,325	\$5,601	\$5,752	\$5,845	\$5,886	\$6,096	\$6,416
†	inch	cm			inch	cm			yd	m
Overall Width	33.00	83.82	Seat Width		23.00	58.42	COM Requirement ¹⁰		6	5.49
Overall Width <i>side table up</i>	45.00	114.30	Seat Depth		20.00	50.80				
Overall Depth	37.00	93.98	Seat Height		21.00	53.34	Fold Down Side Table		10 x 19	25 x 48
Overall Depth <i>trendelenburg</i>	73.00	185.42	Arm Width		3.50	8.89				
Overall Height	46.00	116.84	Width <i>between arms</i>		24.00	60.96			lb	kg
Back Width	23.00	58.42	Arm Height <i>from floor</i>		26.75	67.95	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
Back Height <i>above seat</i>	30.00	76.20	Tente® Casters		4.00	10.16	Shipping Weight ⁸		210	95.25

IN597P Options & Accessories:

details & photos on page 35*

Product Info

Heat & Massage [†]	add	IV Mounting	add	Headrest [†]	add	Extended	Folding Table	add
None IN597Pxx	\$0	Mount FR597GIVMNT ⁵	\$48	Skirt HRSxx	\$78	Legrest not available	Left Side Single FR597GTBLMNT ⁴⁵	\$102
Both IN597PHxx	\$1,852	Pole 2101A	\$64	Cover not available		Legs not available	Swing-Away SAFT20-597	\$447
Heat Only IN597PHOxx	\$1,604			Pillow standard [§]	\$0			



IN587W Extra Wide Clinical Care Recliner with 3 Positions and Trendelenburg

Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$3,592	\$3,592	\$4,071	\$4,191	\$4,271	\$4,547	\$4,698	\$4,791	\$4,832	\$5,042	\$5,362
†	inch	cm			inch	cm			yd	m
Overall Width	33.00	83.82	Seat Width		24.00	60.96	COM Requirement ¹⁰		6	5.49
Overall Width <i>side table up</i>	53.50	135.89	Seat Depth		20.50	52.07			inch	cm†
Overall Depth	45.50	115.57	Seat Height		22.00	55.88	Fold Down Side Table		10 x 19	25 x 48
Overall Depth <i>trendelenburg</i>	78.00	198.12	Arm Width		4.50	11.43				
Overall Height	50.00	127.00	Width <i>between arms</i>		23.25	59.06			lb	kg
Back Width	30.50	77.47	Arm Height <i>from floor</i>		28.00	71.39	Weight Capacity ¹ EVENLY DISTRIBUTED		450	204.12
Back Height <i>above seat</i>	33.50	85.09	Tente® Casters		5.00	12.70	Shipping Weight ⁸		143	64.86

IN587W
Options
& Accessories:

details & photos on page 35★

Product Info

Heat & Massage [‡]			add		IV Mounting		add		Headrest [‡]		add		Extended		add		Folding Table			add	
None	IN587Wxx	\$0	Mount	5782RG	\$54	Skirt	HRSxx	\$78	Legrest	not available		Fold Down Side ^{\$}		\$0							
Both	IN587WHxx	\$1,740	Pole	2101A	\$64	Cover	HRC587xx	\$128	Legs 4" ¹¹	FR5699G	\$40	Swing-Away		SAFT10	\$447						
Heat Only	IN587WHOxx	\$1,494				Pillow	FR5659xxUS	\$98													

IN587WD Extra Wide Drop-Arm Clinical Care Recliner with 3 Positions and Trendelenburg

Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$3,900	\$3,900	\$4,379	\$4,499	\$4,579	\$4,855	\$5,006	\$5,099	\$5,140	\$5,350	\$5,670
†	inch	cm			inch	cm			yd	m
Overall Width	35.75	90.81	Seat Width		24.00	60.96	COM Requirement ¹⁰		6	5.49
Overall Width <i>side table up</i>	58.25	147.96	Seat Depth		20.50	52.07			inch	cm [†]
Overall Depth	44.00	111.76	Seat Height		22.00	55.88	Fold Down Side Table		10 x 19	25 x 48
Overall Depth <i>trendelenburg</i>	78.00	198.12	Arm Width		4.50	11.43	Tente® Casters		5.00	12.70
Overall Height	50.00	127.00	Width <i>between arms</i>		26.50	67.31	 [★]		lb	kg
Back Width	30.50	77.47	Drop Arm Height <i>from floor, arm up</i>		27.50	69.85	Weight Capacity ¹ EVENLY DISTRIBUTED		450	204.12
Back Height <i>above seat</i>	33.50	85.09	Drop Arm Height <i>from floor, arm down</i>		21.50	54.61	Shipping Weight ⁸		145	65.77

IN587WD
Options

details & photos on page 35★

Product Info

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended	add	Folding Table	add
None	IN587WDxx	\$0	Mount 5782RG	\$54	Skirt HRSxx	\$78	Legrest not available	Fold Down Side \$	\$0
Both	IN587WDHxx	\$1,742	Pole 2101A	\$64	Cover HRC587xx	\$128	Legs not available	Swing-Away SAFT10	\$447
Heat Only	IN587WDHOxx	\$1,496			Pillow FR5659xxUS	\$98			



IN566G Deluxe Clinical Care Recliner with 3 Positions and Trendelenburg

 Product Info

IN566G Options & Accessories:
 details & photos on page 35*
Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$3,086	\$3,086	\$3,565	\$3,685	\$3,765	\$4,041	\$4,192	\$4,285	\$4,326	\$4,536	\$4,856
†	inch	cm		inch	cm				yd	m
Overall Width	31.25	79.38	Seat Width	20.50	52.07	COM Requirement ¹⁰			6	5.49
Overall Width <i>side table up</i>	51.50	130.81	Seat Depth	19.00	48.26				inch	cm†
Overall Depth	37.50	95.25	Seat Height	22.25	56.52	Fold Down Side Table			10 x 19	25 x 48
Overall Depth <i>trendelenburg</i>	75.00	190.50	Arm Width	4.50	11.43					
Overall Height	49.50	125.73	Width <i>between arms</i>	21.50	54.61	 *			lb	kg
Back Width	27.50	69.85	Arm Height <i>from floor</i>	28.50	72.39	Weight Capacity ¹ EVENLY DISTRIBUTED			400	181.44
Back Height <i>above seat</i>	32.00	81.28	Tente® Casters	5.00	12.70	Shipping Weight ⁸			140	64.41



Heat & Massage [‡]	add	IV Mounting	add	Headrest [‡]	add	Extended	add	Folding Table	add
None IN566Gxx	\$0	Mount 5782RG	\$54	Skirt HRSxx	\$78	Legrest	not available	Fold Down Side [§]	\$0
Both IN566GHxx	\$1,666	Pole 2101A	\$64	Cover HRC566xx	\$134	Legs 4" ¹¹ FR5699G	\$40	Swing-Away SAFT10	\$447
Heat Only IN566GHOxx	\$1,420			Pillow FR5659xxUS	\$98				

IN566DG Deluxe Drop-Arm Clinical Care Recliner with 3 Positions and Trendelenburg

 Product Info

IN566DG Options & Accessories:
 details & photos on page 35*
Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$3,518	\$3,518	\$3,997	\$4,117	\$4,197	\$4,473	\$4,624	\$4,717	\$4,758	\$4,968	\$5,288
†	inch	cm		inch	cm				yd	m
Overall Width	32.50	82.55	Seat Width	20.50	52.07	COM Requirement ¹⁰			6	5.49
Overall Width <i>side table up</i>	52.50	133.35	Seat Depth	19.00	48.26				inch	cm†
Overall Depth	38.00	96.52	Seat Height	22.50	57.15	Fold Down Side Table			10 x 19	25 x 48
Overall Depth <i>trendelenburg</i>	75.00	190.50	Arm Width	4.50	11.43	Tente® Casters			5.00	12.70
Overall Height	50.50	128.27	Width <i>between arms</i>	22.00	55.88				lb	kg
Back Width	27.50	69.95	Drop Arm Height <i>from floor, arm up</i>	28.00	71.12	Weight Capacity ¹ EVENLY DISTRIBUTED			400	181.44
Back Height <i>above seat</i>	30.00	76.20	Drop Arm Height <i>from floor, arm down</i>	21.50	54.61	Shipping Weight ⁸			142	64.41



Heat & Massage [‡]	add	IV Mounting	add	Headrest [‡]	add	Extended	add	Folding Table	add
None IN566DGxx	\$0	Mount 5782RG	\$54	Skirt HRSxx	\$78	Legrest	not available	Fold Down Side [§]	\$0
Both IN566DGHxx	\$1,666	Pole 2101A	\$64	Cover HRC566xx	\$134	Legs	not available	Swing-Away SAFT10	\$447
Heat Only IN566DGHOxx	\$1,420			Pillow FR5659xxUS	\$98				



vinyl: Oatmeal 8802

IN577RG Multi-Care Recliner with 3 Positions and Trendelenburg**Specs:** †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$2,962	\$2,962	\$3,441	\$3,561	\$3,641	\$3,917	\$4,068	\$4,161	\$4,202	\$4,412	\$4,732
	inch	cm		inch	cm					yd	m
Overall Width	31.50	80.01	Seat Width	21.50	54.61	COM Requirement ¹⁰				6	5.49
Overall Width side table up	52.00	132.08	Seat Depth	19.50	49.53						
Overall Depth	40.50	102.87	Seat Height	20.50	52.07	Fold Down Side Table				10 x 19	25 x 48
Overall Depth trendelenburg	73.50	186.69	Arm Width	4.50	11.43						
Overall Height	48.00	121.92	Width between arms	21.50	54.61	Weight Capacity ¹ EVENLY DISTRIBUTED				350	158.76
Back Width	21.50	54.61	Arm Height from floor	27.50	69.85						
Back Height above seat	34.00	86.36	Tente® Casters	5.00	12.70	Shipping Weight ⁸				129	58.51

IN577RG**Options
& Accessories:**

details & photos on page 35★

Product Info

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended	add	Folding Table	add
None IN577RGxx	\$0	Mount FR5782RG	\$54	Skirt HRSxx	\$78	Legrest	not available	Fold Down Side \$	\$0
Both IN577RGHxx	\$1,666	Pole 2101A	\$64	Cover HRC577xx	\$115	Legs 4" 11 FR5699G	\$40	Swing-Away SAFT10	\$447
Heat Only IN577RGHOxx	\$1,420			Pillow FR5659xxUS	\$98				

**600 SERIES****Blood Draw / Phlebotomy Seating****620** Upholstered Bariatric Blood Draw / Phlebotomy Chair with Adjustable Height, Padded, Arm Rest and Wallsaver Wireform**Specs:** †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$1,840	\$1,840	\$2,040	\$2,090	\$2,123	\$2,238	\$2,301	\$2,340	\$2,357	\$2,644	\$2,578
	inch	cm		inch	cm					yd	m
Overall Width	35.00	88.90	Seat Width	25.50	64.77	COM Requirement ¹⁰				2.5	2.29
Overall Depth	28.00	71.12	Seat Depth	18.75	47.63						
Overall Height	37.50	95.25	Seat Height	22.00	55.88	Fold Down Side Table				10 x 19	25 x 48
Width Between Arms	25.00	63.50	Back Width	26.00	66.04						
Arm Height ⁸ - low	26.50	67.31	Back Height	16.00	40.64	Weight Capacity ¹ EVENLY DISTRIBUTED				800	362.87
Arm Height ⁸ - high	32.00	81.28	Powder Coat Finish	foghorn							
						Shipping Weight ⁸				55	24.95

620 Options: ★

Single Flip Arm†\$	add	Dual Flip Arms†	add
620-xx	\$0	620-DA-xx	\$172



460 High / Low Power Exam Table with Power Assisted Back, Hand Control, and Stirrups

 Product Info
Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$8,777	\$8,777	\$9,017	\$9,077	\$9,117	\$9,255	\$9,330	\$9,377	\$9,397	\$9,502	\$9,662

	inch	cm								
Overall Width	27.00	68.58	Pelvic Tilt	standard	Power Assisted	back				
Overall Length leg rest extended	72.00	182.88	Pull Out Leg Rest	standard	UL 60601-1	yes				
Overall Height ³ - low	18.50	46.99	Treatment Tray	integrated					yd	m
Overall Height ³ - high	35.50	90.17	Stirrups	five position continuous lock	COM Requirement ¹⁰	3	2.74			
Paper Roll Holder Widths Available	18.00	45.72	Base	dual pedestal				lb	kg	
	19.00	48.26	Drive	DC electromechanical	Weight Capacity ¹ EVENLY DISTRIBUTED	600	272.16			
	20.00	50.80	Powder Coat Finish	foghorn	Shipping Weight ⁸	395	179.7			

Options: †

Model [†]	add	Foot Switch [†]	add	Assist Arms [†]	add	Physician Storage Drawer [†]	add
Base 460-xx	\$0	460-FS-xx	\$512	460-AA-xx	\$722	Drawer Only 460-PS-xx	\$397
Loaded 460-FL-xx	\$1,839					With Warmer 460-PS-DW-xx	\$608


420 Exam Table with 3 Pass-Through Drawers, 2 Storage Drawers, and Stirrups

 Product Info
Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$3,171	\$3,171	\$3,411	\$3,471	\$3,511	\$3,649	\$3,724	\$3,771	\$3,791	\$3,896	\$4,056

	inch	cm		inch	cm					
Overall Width	27.00	68.58	Drawer Width ⁷	18.00	45.72	COM Requirement ¹⁰	3	2.74		
Overall Length leg rest extended	72.00	182.88	Drawer Length ⁷	16.50	41.91					
Overall Height	32.00	81.28	Drawer Height ⁷	4.00	10.16					
Foot Step Width	18.00	45.72	Pelvic Tilt	standard				lb	kg	
Foot Step Length	17.00	43.18	Pull Out Leg Rest	standard	Foot Step Weight Capacity ¹ EVENLY DISTRIBUTED	800	362.87			
Paper Roll Holder Widths Available	18.00	45.72	Treatment Tray	integrated	Table Weight Capacity ¹ EVENLY DISTRIBUTED	500	226.80			
	19.00	48.26	Stirrups	five position continuous lock	Shipping Weight ⁸	335	151.95			
	20.00	50.80	Powder Coat Finish	foghorn						

Options: †

Model	Integrated Electrical Outlet [†]	add	Integrated Electrical Outlet [†]	add
Base 420-xx \$	Outlet Only 420-E0-xx	\$281	with Drawer Warmer 420-DW-xx	\$669



415 Exam Room Treatment Table with Shelf, Adjustable Back, and Contour Top


[Product Info](#)

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,733	\$2,733	\$3,052	\$3,132	\$3,186	\$3,370	\$3,470	\$3,532	\$3,560	\$3,700	\$3,913

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	 Weight Capacity¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90				Shipping Weight⁸	190	90.72
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn				

Options:

No Drawers ⁺	add	Single Pass-Through Drawer ⁺	add	Dual Pass-Through Drawer ⁺	add
415-xx ⁵	\$0	415D-xx	\$310	415DD-xx	\$620
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

414 Exam Room Treatment Table with Shelf and Adjustable Back


[Product Info](#)

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,706	\$2,706	\$3,025	\$3,105	\$3,159	\$3,343	\$3,443	\$3,505	\$3,533	\$3,673	\$3,886

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	 Weight Capacity¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90				Shipping Weight⁸	200	90.72
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn				

Options:

No Drawers ⁺	add	Single Pass-Through Drawer ⁺	add	Dual Pass-Through Drawer ⁺	add
414-xx	\$0	414D-xx ⁵	\$310	414DD-xx	\$620
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

413 Exam Room Treatment Table with Shelf and Flat Top


[Product Info](#)

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,190	\$2,190	\$2,509	\$2,589	\$2,643	\$2,827	\$2,927	\$2,989	\$3,017	\$3,157	\$3,370

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	 Weight Capacity¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90				Shipping Weight⁸	200	90.72
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn				

Options:

No Drawers ⁺	add	Single Pass-Through Drawer ⁺	add	Dual Pass-Through Drawer ⁺	add
413-xx	\$0	413D-xx	\$310	413DD-xx ⁵	\$620
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

41xD Single Drawer Option:



41xDD Double Drawer Option:



412 Exam Room Treatment Table with Contoured Adjustable Back



 Product Info

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,337	\$2,337	\$2,656	\$2,736	\$2,790	\$2,974	\$3,074	\$3,136	\$3,164	\$3,304	\$3,517

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	 Weight Capacity¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90				Shipping Weight⁸	180	81.65
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn				

Options:

No Drawers ⁺	add	Single Pass-Through Drawer ⁺	add	Dual Pass-Through Drawer ⁺	add
412-xx ⁵	\$0	412D-xx	\$310	412DD-xx	\$620
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

411 Exam Room Treatment Table with Adjustable Back



 Product Info

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,310	\$2,310	\$2,629	\$2,709	\$2,763	\$2,947	\$3,047	\$3,109	\$3,137	\$3,277	\$3,490

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	 Weight Capacity¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90				Shipping Weight⁸	180	81.65
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn				

Options:

No Drawers ⁺	add	Single Pass-Through Drawer ⁺	add	Dual Pass-Through Drawer ⁺	add
411-xx	\$0	411D-xx ⁵	\$310	411DD-xx	\$620
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

410 Exam Room Treatment Table with Flat Top



 Product Info

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,632	\$1,632	\$1,951	\$2,031	\$2,085	\$2,269	\$2,369	\$2,431	\$2,459	\$2,599	\$2,812

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	 Weight Capacity¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90				Shipping Weight⁸	180	81.65
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn				

Options:


No Drawers ⁺	add	Single Pass-Through Drawer ⁺	add	Dual Pass-Through Drawer ⁺	add
410-xx	\$0	410D-xx ⁵	\$310	410DD-xx ⁵	\$620
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

41xD Single Drawer Option:



41xDD Double Drawer Option:




RT460 Replacement Top with Foot Pad for Intensa 460 Exam Tables

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$790	\$790	\$1,030	\$1,090	\$1,130	\$1,268	\$1,343	\$1,390	\$1,410	\$1,515	\$1,675
	inch		cm		inch		cm		yd		m
Overall Width	27.50		69.85		Foot Rest Width		19.00		3		2.74
Overall Length	60.00		152.40		Foot Rest Length		18.00		lb		kg
Paper Roll Holder	included		Mounting Hardware		included		Shipping Weight ⁸		55		24.95

 Order Number[‡] RT460-xx

 Product Info

RT420 Replacement Top with Foot Pad for Intensa 420 Exam Tables

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$790	\$790	\$1,030	\$1,090	\$1,130	\$1,268	\$1,343	\$1,390	\$1,410	\$1,515	\$1,675
	inch		cm		inch		cm		yd		m
Overall Width	27.50		69.85		Foot Rest Width		19.00		3		2.74
Overall Length	60.00		152.40		Foot Rest Length		18.00		lb		kg
Paper Roll Holder	included		Mounting Hardware		included		Shipping Weight ⁸		55		24.95

 Order Number[‡] RT420-xx

 Product Info


RT404 Replacement Top with Foot Pad for 404 Exam Tables

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$650	\$650	\$890	\$950	\$990	\$1,128	\$1,203	\$1,250	\$1,270	\$1,375	\$1,535
	inch		cm		inch		cm		yd		m
Overall Width	27.50		69.85		Foot Rest Width		20.50		3		2.74
Overall Length	60.00		152.40		Foot Rest Length		17.50		lb		kg
					Hook & Loop Attachment		included		Shipping Weight ⁸		35 15.88

 Order Number[‡] RT404-xx


 Product Info

RT304 Replacement Top with Foot Pad for 304 Exam Tables

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$881	\$881	\$1,121	\$1,181	\$1,221	\$1,359	\$1,434	\$1,481	\$1,501	\$1,606	\$1,766
	inch		cm		inch		cm		yd		m
Overall Width	27.50		69.85		Foot Rest Width		20.00		3		2.74
Overall Length	55.00		139.70		Foot Rest Length		18.00		lb		kg
									Shipping Weight ⁸		55 24.95

 Order Number[‡] RT304-xx

 Product Info


RT204 Replacement Top with Foot Pad for 204 Exam Tables**Specs:** †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$881	\$881	\$1,121	\$1,181	\$1,221	\$1,359	\$1,434	\$1,481	\$1,501	\$1,606	\$1,766
	inch		cm		inch		cm		yd		m
Overall Width	27.50		69.85		Foot Rest Width		20.00		COM Requirement ¹⁰		3
Overall Length	55.00		139.70		Foot Rest Length		18.00		45.72		2.74
Lift Handles	not included		Top Attachment Knobs		not included		Shipping Weight ⁸		55		24.95

Order Number* RT204-xx




Does not include lift handles and top attachment knobs
If the OEM model has a seamless top, request paper roll holder when ordering

 Product Info
RT104 Replacement Top with Foot Pad for 104 Exam Tables**Specs:** †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$938	\$938	\$1,178	\$1,238	\$1,278	\$1,416	\$1,491	\$1,538	\$1,558	\$1,663	\$1,823
	inch		cm		inch		cm		yd		m
Overall Width	27.50		69.85		Foot Rest Width		20.00		COM Requirement ¹⁰		3
Overall Length	55.00		139.70		Foot Rest Length		18.00		45.72		2.74
Paper Roll Holder	included		Mounting Hardware		included		Shipping Weight ⁸		55		24.95

Order Number* RT104-xx

 Product Info
RT100 Replacement Top with Foot Pad for 100 Exam Tables**Specs:** †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$938	\$938	\$1,178	\$1,238	\$1,278	\$1,416	\$1,491	\$1,538	\$1,558	\$1,663	\$1,823
	inch		cm		inch		cm		yd		m
Overall Width	27.50		69.85		Foot Rest Width		20.00		COM Requirement ¹⁰		3
Overall Length	55.00		139.70		Foot Rest Length		18.00		45.72		2.74
Paper Roll Holder	included		Mounting Hardware		included		Shipping Weight ⁸		55		24.95

Order Number* RT100-xx

 Product Info


202 Wallsaver, Stacking, Patient / Guest Chair with Slant Arms


[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$586	\$586	\$666	\$686	\$699	\$745	\$770	\$786	\$793	\$828	\$881

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	22.50	57.15	Seat Width	18.00	45.72	COM Requirement ¹⁰	1	0.91
Overall Depth	22.00	55.88	Seat Depth	18.00	45.72			
Overall Height	32.50	82.55	Seat Height	18.00	45.72			
Back Width	18.00	45.72	Arm Width Inside	19.50	49.53		lb	kg
Back Height	12.50	31.75	Arm Height	26.00	66.04	Weight Capacity ¹ EVENLY DISTRIBUTED	250	113.40
Steel Tubing	14 gauge		Powder Coat Finish	black		Shipping Weight ⁸	28	12.70

 Order Number[‡] 202-xx

201 Wallsaver, Stacking, Patient / Guest Chair with Straight Arms


[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$560	\$560	\$640	\$660	\$673	\$719	\$744	\$760	\$767	\$802	\$855

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	23.50	59.69	Seat Width	18.00	45.72	COM Requirement ¹⁰	1	0.91
Overall Depth	22.00	55.88	Seat Depth	18.00	45.72			
Overall Height	32.50	82.55	Seat Height	18.00	45.72			
Back Width	18.00	45.72	Arm Width Inside	20.00	50.80		lb	kg
Back Height	12.50	31.75	Arm Height	26.00	66.04	Weight Capacity ¹ EVENLY DISTRIBUTED	250	113.40
Steel Tubing	14 gauge		Powder Coat Finish	black		Shipping Weight ⁸	28	12.70

 Order Number[‡] 201-xx

200 Wallsaver, Stacking, Patient / Guest Chair without Arms



[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$524	\$524	\$604	\$624	\$637	\$683	\$708	\$724	\$731	\$766	\$819

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	19.50	49.53	Seat Width	18.00	45.72	COM Requirement ¹⁰	1	0.91
Overall Depth	22.00	55.88	Seat Depth	18.00	45.72			
Overall Height	32.50	82.55	Seat Height	18.00	45.72		lb	kg
Back Width	21.00	53.34	Steel Tubing	14 gauge		Weight Capacity ¹ EVENLY DISTRIBUTED	250	113.40
Back Height	12.50	31.75	Powder Coat Finish	black		Shipping Weight ⁸	26	11.79

 Order Number[‡] 200-xx


179 Overbed Table with Heavy-Duty Steel, Low Profile Base, and Raised Edge Thermofoil Dual Top with Cup Holder

Product Info

finish options page 36*

179TFx-y-xx
\$1,676**Specs:**

	inch	cm
Overall Width	46.00	116.84
Overall Depth	19.00	48.26
Overall Height - low	28.00	71.12
Overall Height - high	44.00	111.76
Top Depth	17.00	43.18

	inch	cm
Base Width	32.50	82.85
Base Depth	19.00	48.26
Base Profile Height	1.94	4.93
Main Top Width	32.00	81.28
Flip Top Width	13.50	34.29

 Solid Surface top available
- contact customer service


	lb	kg
Weight Capacity ¹ EVENLY DISTRIBUTED	100	45.36
Shipping Weight ⁸	43	19.50

Options:

Thermofoil Top Shape [‡]	add
Rectangle 179TFR-y-xx\$	\$0
Kidney 179TFK-y-xx	\$135

Vanity with Mirror [‡]	add
179TFx-y-xx-VM\$	\$210

Wire Basket [‡]	add
179TFx-y-xx-BK\$	\$220

vanity with mirror and wire basket 177TFx-y-xx-VM-BK

177 Overbed Table with Heavy-Duty Steel, Low Profile Base and Laminate Top *or* Raised Edge Thermofoil Top

Product Info

finish options page 36*

laminate: 177-y-xx
\$1,460**Specs:**

	inch	cm
Overall Width	32.50	82.55
Overall Depth	19.00	48.26
Overall Height - low	28.00	71.12
Overall Height - high	44.00	111.76

	inch	cm
Base Width	32.50	82.85
Base Depth	19.00	48.26
Base Profile Height	1.94	4.93
Top Width	32.00	81.28
Top Depth	17.00	43.18

 Solid Surface top available
- contact customer service


	lb	kg
Weight Capacity ¹ EVENLY DISTRIBUTED	100	45.36
Shipping Weight ⁸	36	16.33

Options:

Thermofoil Top Shape [‡]	add
Rectangle 177TFR-y-xx\$	\$135
Kidney 177TFK-y-xx	\$135

Optional [‡]	add
Cup Holder	\$102

Vanity with Mirror [‡]	add
177TFx-y-xx-VM\$	\$210

Wire Basket [‡]	add
177TFx-y-xx-BK\$	\$220

vanity with mirror and wire basket 177TFx-y-xx-VM-BK

175 Overbed Table with Heavy-Duty Steel, Infinite Stop, Low Profile 'H' Base and Thermofoil Top

Product Info

finish options page 36*

base: A919000i\$
\$495top: OS3DLOBTF-xx\$
\$260

order separately / ships separately*

Specs:

	inch	cm
Overall Width	35.00	88.90
Overall Depth	22.00	55.88
Overall Height - low	29.00	73.66
Overall Height - high	43.00	109.20
Range Stops	automatic/infinite	

	inch	cm
Base Width	29.30	74.40
Base Depth	18.20	46.20
Base Profile Height	1.94	4.93
Top Width	35.00	88.90
Top Depth	18.00	45.72



	lb	kg
Weight Capacity ¹ EVENLY DISTRIBUTED	100	45.36
Shipping Weight ⁸	43	19.50

Options:


Thermofoil Top ^{‡*}	add
Flat Edge OS3DLOBTF-xx\$	\$0
Raised Edge OS3DLOBTR-xx	\$8

Talas Cherry Finish	add
OS3DLOBTx-TC\$	\$32


Talas Cherry Finish Option at Additional Charge

perceived texture and coloring in images may vary on finished product. monitor settings effect appearance.

Laminate Top - Rectangular




Thermofoil Top - Rectangular



Thermofoil Top - Kidney




Raised Edge



Cup Holder



Wire Basket

10 x 6 x 9 inch
25.40 x 15.24 x 22.86 cm

Vanity Drawer w/ Mirror



102BAR Bariatric Step Stool, Slip Resistant102BAR
\$220**Specs:** †

	inch	cm
Overall Width	20.00	50.80
Overall Length	26.50	67.31
Overall Height	8.00	20.32
Platform Width	14.75	37.47
Platform Length	19.00	48.26

	inch	cm
Step Height	8.00	20.32
Steel Tubing	1.00	2.54
Steel Tubing Finish	epoxy coating	
Frame Finish	gray / black	

Packaged	1/ea	
Ships Assembled	yes	
 Weight Capacity¹ EVENLY DISTRIBUTED	lb	kg
	1000	453.59
Shipping Weight⁸	9	4.08



Product Info

100BAR Bariatric Chrome Step Stool with Slip Resistant Mat100BAR
\$175**Specs:** †

	inch	cm
Overall Width	14.00	35.56
Overall Length	17.00	43.18
Overall Height	9.50	24.13
Platform Width	11.00	27.94
Platform Length	14.00	35.56

	inch	cm
Step Height	9.50	24.13
Steel Tubing	1.00	2.54
Steel Tubing Finish	chrome plated	
Frame Finish	gray / black	

Packaged	1/ea	
Ships Assembled	yes	
 Weight Capacity¹ EVENLY DISTRIBUTED	lb	kg
	500	226.80
Shipping Weight⁸	11	4.99



Product Info

100 Chrome Step Stool with Slip Resistant Mat100
\$99**Specs:** †

	inch	cm
Overall Width	14.00	35.56
Overall Length	17.00	43.18
Overall Height	9.00	22.86
Platform Width	11.00	27.94
Platform Length	14.00	35.56

	inch	cm
Step Height	9.00	22.86
Steel Tubing	1.00	2.54
Steel Tubing Finish	chrome plated	
Frame Finish	gray / black	

Packaged	2/cs	
Ships Assembled	yes	
 Weight Capacity¹ EVENLY DISTRIBUTED	lb	kg
	250	113.40
Shipping Weight⁸	14	6.35



Product Info



500 Pedestal Exam Light with Chrome-Plated Steel Base and Adjustable Goose Neck

500

\$233

Product Info

Specs: †

	inch	cm
Overall Height	74.00	187.96
Height Range -low	34.50	87.63
Height Range -high	53.50	135.89
Base Diameter	11.00	27.94
Gooseneck Arm	12.00	30.48
	ft	m
Power Cord ^U 3-wire grounded	8.00	2.44
UL Approved	yes	

Steel Shaft Finish	chrome-plated
Pharmacy Shade exterior finish	chrome-plated
Pharmacy Shade interior finish	white baked enamel
Gooseneck Finish	chrome-plated
Steel Shaft	adjustable
Height Adjustment method	SOK-IT-GUARD®
Lamp Head adjustable	360°
Light Bulb Included	no

	W	hp(M)
Max LED Equivalent	10	0.014
Max CFL Equivalent	13	0.018
Max Incandescent	60	0.082
Ships Assembled	no	
	lb	kg
Shipping Weight ⁸	34	15.42

130-2 Chrome-Plated Steel IV Pole with 4-Leg Steel Base, Rubber Caps, and 2-Hooks

130-2


\$119

Product Info

Specs: †

	inch	cm
Base Diameter	22.00	55.88
Overall Height -low	51.00	129.54
Overall Height -high	90.00	228.6
Rubber Wheels	2.00	5.08
Steel Tubing	1.00	2.54

Ram Style Hooks	2	Steel Tubing Finish	chrome-plated	
		Height Adjustment	twist-lock	
Packages Shipped	2		lb	kg
Ships Assembled	no	Shipping Weight ⁸	13	5.90

Footnotes Legend

★ ships freight: request shipping quote from customer service

◆ ships ground: no additional charge

web link - additional product information

page link - jumps to referenced page

‡ xx in the product number is the placeholder for vinyl/finish choice

† standard rounding method used

§ representative image shows this option

★ options available at additional cost & may alter weight specs - order requirements noted

f. this section of cleaning guide does not apply to Dolce fabric

d. procedure arm is not compatible/available with either backrest option

1. maximum weight capacity - evenly distributed - refer to user manual

2. includes resident, support surface, head/footboards, and accessories. max resident weight (evenly distributed) is 50 lb less - refer to user manual

3. tolerance: ± 0.5 inch (± 1.27 cm)

4. supports a single, left side, table only

5. requires factory installation include item with purchase

6. pullout footplate not compatible/available with extended legrest option

7. interior dimensions

8. approximate measurement, options may alter spec

9. product is not CFFA-HealthCare 201 Certified

10. submit required sample size indicated to:
Intensa COM Department, 1810 S. Elm Street, High Point, NC 27260.

11. adjust height spec: +4 inch (+10.16 cm)

12. replacement cylinders may be ordered separately at additional cost

A. accessories not included - sold separately

B. tolerance: ± 1 inch (± 2.54 cm)

C. call for pricing

D. side-by-side

F. by CF Stinson

L. limited availability on certain colors

M. multiple factors affect fabric durability and appearance retention

N. not pictured

P. details, options, and weight info on page indicated

R. recommendations do not serve as a product guarantee. use manufacturer instructions

S. the term "safe", pertains to surface effects and not in reference to human health

T. term and corresponding data refer to typical performance in specific tests

U. UL listed

Order Numbers: Multiple Options, Vinyl, & Finish Codes



Combining Options

Most single options show the full order number. When combining multiple options, pull each option code from the option order number [or as shown below]. To create the full order number, begin with the main item number (972), then each option code (CB) (-FR), a dash (-), seat height, a dash (-), vinyl/finish codes (see right), a dash (-), then wheel choice.

Order Number for 972 stool with crescent back, footring, medium seat height, with soft wheels in Tier 1 Violet vinyl : 972CB-7-FR-VI-SW.

Order Number for 972 stool with procedure arm, footring, with standard medium seat height & hard wheels in Tier 1 Violet vinyl : 972PA-FR-VI.

Customer's Own Material

Due to the nature of our product, some COM material may not perform or upholster to our standards. Intensa will not be responsible for materials that do not meet our standards. All customer material must be approved prior to manufacture. Submit required sample size indicated to:

Intensa COM Department, 1810 S. Elm Street, High Point, NC 27260.

Vinyl/Finish Codes

Most vinyl or finish codes follow the product code (see left) and start with a dash (-) followed by the color code.

Approved Customer Vinyls: add '-COM' followed by '-manufacturer name-pattern name-color name'.

- 972CB-7-FR-COM-DesignTex-Chasam-Blue-SW

Tier 1 Vinyls: add the two-letter color code. Example using Alabaster.

- 972CB-7-FR-AL-SW

Tiers 2 - 10 Vinyls: add the tier number '-T#' followed by '-manufacturer name-pattern name-color name'.

- 972CB-7-FR-T3-CFStinson-Tessuto-Lime-SW

OBT Finishes: add the base code letter (-y) followed by the two-letter top finish code (-xx). 175 does not have an optional base color.

- 179-G-MC • 175-SO

Head/Footboard Finishes: add the two-letter finish code at the end (-xx).

- H0336ST3DL-AR • F18R36ST-SO

Info & Notices

Pricing & Availability

Current at time of publication & subject to change without notice. Intensa is not responsible for typographical errors

Contact Customer Service

For most current information, pricing, freight quotes, and support
+1.336.884.4096 | intensa@grahamfield.com

Color & Texture Variations

Images and swatches are digital reproductions. perceived texture and coloring may vary on finished product. Monitor and printer settings effect appearance of color and texture. Order physical swatches using chart links [or](#) contact customer service.

Purchase Orders

Please remit purchase orders to intensa@grahamfield.com

Intensa Limited Warranties

Intensa Terms & Conditions

California Prop 65 Warnings


Green Initiative

Physician & Laboratory Seating Options

T-Style Arm Rest [†]	add	Rotating Procedure Arm ^{†‡}	add	Upholstered Backrest ^{†‡}	add	Crescent Backrest ^{†‡}	add	Chrome Footring [†]	add
### -T-xx	\$187	### PA-xx	\$450	### B-xx	\$235	### CB-xx	\$315	### -FR-xx	\$110
adjustable, urethane		adjustable, upholstered, rotates 360°		adjustable, black plastic outer cover		adjustable, upholstered		with twist lock feature	



Casters <i>set of five</i>	add
Hard Wheel standard	\$0
Soft Wheel Urethane Tread ### -SW (hard floors)	\$38
Soft Wheel Self Braking ### -SBCSW (roll weighted)	\$40
Soft Wheel Braking ### -BCSW (roll unweighted)	\$50
Stationary Glides ### -GLD (±2 inch height vs caster)	\$0



Gas Cylinders			add		inch	cm
5 [five] ^N	###-5	\$0		Height	5.50	14.00
7 [seven]	###-7	\$38		Height	7.87	20.00
10 [ten]	###-10	\$48		Height	10.26	27.00

^N not pictured



Universal Bed Socket IV Pole

Easily slides into IV Pole Attachment. Internal rubber bumpers ensure silent operation and reduce shock. Constructed of durable, anodized aluminum tubing. Adjustable height from 35 to 65 inches (88.90 - 165.10 cm).

Item No. 2101A **add \$64**

IV Pole Attachment

For use with Bed Socket IV Pole. FR601 & FR597 series require factory installation - add item with order.

Item No. FR601IVMNT⁵ (shown left)

■ IN5601 series

Item No. FR597GIVMNT^{4,5}

■ IN597 series

Item No. 5782RG (shown right)

■ IN566 ■ IN577 ■ IN587 series



Wall saver Arm Caps

Plastic border prevents scarring is a standard feature. No changes to the item number

■ IN566 ■ IN577 ■ IN587 series **add \$0**



Headrest Cover

Protects and extends the life of recliners vinyl. Provides added protection from hair oil, hairspray, coloring agents, etc.

Item No. HRC566xxx ■ IN566 series

Item No. HRC587xxx ■ IN587 series

Item No. HRC577xxx ■ IN577 series

Item No. HRC601xxx ■ IN601 series



Headrest Skirt

Easily mounts to recliner back and provides added protection to head-section vinyl from hair oil, hairspray, coloring agents, etc.

Item No. HRSxxx **add \$78**

■ IN566 ■ IN577 ■ IN587

■ IN597 ■ IN601 series



Fold Down Side Tables

Smooth surface edges makes cleaning and disinfection easier. Side tables are level with armrest in the up position.

Item No. FR597GTBLMNT^{4,5}

■ IN597 — single left side table only



Swing-Away Folding Table

360 degree swivel table. Scratch resistant table top is durable and easy to clean. Fold down and out of the way when not needed.

Item No. SAFT20-597 **add \$447**

■ IN597 series

Item No. SAFT20-601 **add \$447**

■ IN601 series

Item No. SAFT10 **add \$447**

■ IN566 ■ IN577 ■ IN587 series



Universal Headrest Pillow

The unique shape features an angled design that offers the patient a separate level of support depending on whether the small contour or large contour is placed under the patient's head. Hook and loop attachment allows for variable positioning along back of recliner.

Item No. FR5659xxxUS **add \$98**



Leg Extensions

Fits all round tube frame designs. Raises recliner by four inches (10.16 cm).

Item No. FR5699G **add \$40**

■ IN566 ■ IN577 ■ IN587 series

* ships freight if purchased with product
◇ ships ground where possible on individual purchase


175 Top Finishes [xx]




AR ARIA



CO COCOBALA



SC SHAKER CHERRY




SO SOLAR OAK



TC TALAS CHERRY




WC WILD CHERRY




WM WINDSOR MAHOGANY


177-179 Top Finishes [xx]



CL CANDLELIGHT



CP CHOCOLATE PEAR



DV DOLCE VITA



HM HARD ROCK MAPLE



MO MEDIUM OAK

175 Recommended Care & Cleaning^{R,S}

Clean the surface with a damp cloth or sponge and a mild soap or detergent. Rinse thoroughly with warm water and wipe dry.



Difficult stains such as coffee or tea can be removed using a mild household cleaner/detergent and a soft bristled brush, repeating as necessary.



If a stain persists, use a paste of baking soda and water and apply with a soft bristled brush. Light scrubbing for 10 to 20 strokes should remove most stains. Although baking soda is a low abrasive, excessive scrubbing or exerting too much force could damage the surface.



Stubborn stains that resist any of the above cleaning methods may require the use of undiluted household bleach or nail polish remover. Apply to the stain and let stand no longer than two minutes. Rinse thoroughly with warm water and wipe dry. May repeat if stain appears to be going away and the laminate color has not been affected.



After cleaning rinse surfaces using a water-dampened cloth, and dry thoroughly as some cleaning agents are corrosive by nature, and if not rinsed their residue can cause permanent damage.

Recommended Household Cleaners^{R,S}

Windex® • Glass Plus® • Fantastik® • Mr. Clean® • Formula 409® • Isopropyl Alcohol

Warning

Prolonged exposure to bleach will cause discoloration. Always rinse surfaces after cleaning! Failure to rinse after cleaning can cause damage; even if a small amount of cleaning solution remains on the surface.

179 - 177 Recommended Care & Cleaning^{R,S}

ALWAYS rinse overbed table with a water-dampened cloth to remove any cleaning agent residue, and dry thoroughly.



Recommend using a mild cleaning agent or detergent when cleaning overbed table. Some cleaning agents are corrosive by nature, and if not rinsed their residue can cause permanent damage.



If other hospital-grade cleaning agents are used, follow the instructions provided by the manufacturer.




After cleaning rinse surfaces using a water-dampened cloth, and dry thoroughly as some cleaning agents are corrosive by nature, and if not rinsed their residue can cause permanent damage.

A Few Notes of Caution

- Acidic or abrasive cleaners can damage laminate surfaces; do not use them.
- Hair, textile and food dyes can cause permanent stains. If spilled, wipe it up immediately with dish-washing detergent or an all-purpose cleaner.
- Steel wool and other abrasive pads will damage the laminate face. Do not use them for cleaning.
- Drain and Toilet bowl cleaners contain harsh chemicals that can cause permanent damage. If spilled, wipe up immediately, wash surface with soapy water and rinse several times.

177-179 OBT Base Colors [y]



A ALMOND



G GRAPHITE GREY



N NEVADA

swatches are digital reproductions. perceived texture and coloring may vary on finished product. monitor settings effect appearance

^R recommendations do not serve as a product guarantee. use manufacturer instructions ^S the term "safe" herein, pertains to surface effects and not in reference to human health



AL ALABASTER



AQ AQUA GREEN




AR ARTICHOKE



BL BLACK



BU BLUE




BT BRITISH TAN



DO DOE



DB DUSTY BLUE



GR GRAY




MO MOSS




NA NAVY




NO NORTHWOODS




OG OCEAN GREY




PM PALM




PL PEARL




PU PLUM




ST SADDLE TAN




SH SHADOW



WH SUPER WHITE




TL TEA LEAF



TE TEAL



VI VIOLET



WG WARM GRAY



WW WEDGEWOOD

swatches are digital reproductions. perceived texture and coloring may vary on finished product
monitor settings effect appearance. order physical swatches using chart links on page 42

Olympus TIER 6



OLY204 DUNE



OLY202 CORIANDER



OLY01 MELON



OLY213 LAUREL



OLY09 FERN



OLY210 ADMIRAL



OLY208 SHOYA


Natural World TIER 6



NTL221 TIMBER



NTL223 MARIGOLD



NTL222 JUNIPER




NTL228 WALNUT

Raffia TIER 5



SCP2270 ACADEMY




SCP2150 SHOW




SCP2130 FLICK



SCP2230 CINEMA



SCP2190 PREMIER



SCP2210 SILVER SCREEN



SCP2310 PRODUCER


Tessuto TIER 3



SEN224 VENETTO



213 BELLINI



SEN220 MALTA



219 MARBLE



SEN222 LIME



229 MURANO



234 PALAZZO



217 SENECA



SEN225 MINK



233 ADRIATIC



228 MARCO




231 MAJESTY




230 TRUFFLE


Galveston TIER 6




GV-004 TWILIGHT



GV-003 CACTUS




GV-010 CEDAR




GV-007 CREAM




GV-026 STUCCO




GV-011 OXBLOOD



GV-017 DESERT



GV-016 WINCHESTER



GV-014 STREAM

Belmont TIER 4



BL-002 CITRUS




BL-009 CLEMINTINE



BL-016 STONE



BL-007 SISAL



BL-000 GRAVEL



BL-004 OPAL



BL-027 PATINA



BL-024 TIDAL




BL-017 DUNE


Icon TIER 3



IC-003 CELERY



IC-010 COFFEE



IC-017 SAND



IC-013 GRASS




IC-014 OCEAN




IC-034 MARINE



IC-019 TERRA



IC-001 SANGRIA



IC-027 ALMOND

Midway TIER 2



MW-007 CEMENT



MW-009 COPPER



MW-013 ALGAE



MW-027 WHEAT




MW-017 SAND



MW-004 INDIGO



MW-002 SUNSHINE



MW-023 SPA



MW-000 STONE



MW-014 SKYLIGHT



MW-003 MOSS




MW-016 HAZE




MW-006 STORM

swatches are digital reproductions. perceived texture and coloring may vary on finished product
monitor settings effect appearance. order physical swatches using chart links on page 42

Guardian TIER 2



1325 ARMOR



1305 CONVOY




1323 MOAT




1319 SIREN




1324 SENTRY



1311 FENCE



1318 STEEL




1304 LATCH


Savoy TIER 2




9517 EUCALYPTUS



9524 BLACK



409 WARM TAUPE



9516 INDIGO


Cassara TIER 2



9611 AMETHYST



9613 CORK



9607 PATINA




9604 INDIGO


Modena TIER 1



2013 TITAN



2022 ELEPHANT



2002 BLUEBIRD



2008 SEAGRASS




2014 INK



2019 NUTMEG



2007 FRESH



2001 FARRO


Masiala TIER 2



321 BLUESTAR




320 DEMIN




311 UNDERWATER



313 WOODLET



304 CALM




306 BLACK WALNUT




309 BERRY BUSH




317 DAMASCUS

Dolce
TIER 2


6727 METEOR



6717 MOSS



6724 JET




6725 INDIGO

Patriot Plus
TIER 2


NEO




432 IMPERIAL BLUE



8602 OATMEAL




8604 CONCRETE




8567 DOVE




454 ROYAL BLUE




427 BLUERIDGE




8584 TAUPE




8516 CHESTNUT




8611 TERRA COTTA



8504 SAND




8533 PORT



8626 EUCALYPTUS

Renaissance
TIER 2

NEO




9214 STEEL BLUE




9201 COBBLESTONE




9205 GYPSUM



9208 WINEBERRY




9203 PYRITE



9207 ICE BLUE



9210 CYPRESS



9212 BLUE JADE



Color Variations swatches are digital reproductions. perceived texture and coloring may vary on finished product
monitor and printer settings will effect appearance of color and texture
order physical swatches using chart links on page 42

Pricing & Availability current at time of publication & subject to change without notice

Contact Customer Service for most current information, pricing, availability, freight quotes, and support
+1.336.884.4096 | intensa@grahamfield.com



Price Tier and Care Quick Reference | Care Codes detailed on next page | Link Icon connects to complete pattern selection, spec sheet, and care guide.

Stinson®	tier	care	link
Alamo	2	W	link
Amplify	2	W-S	link
Arbor	3	W	link
Ardent	2	W	link
Avant	3	W-S	link
Bandwidth	6	W	link
Blend	2	W	link
Blink Ombre	5	W	link
Box Chain	2	W	link
Burkshire	4	W-S	link
Cassidy 2.0	2	W-S	link
Cesto	3	W-S	link
Channel 2.0	3	W-S	link
Cirque 2.0	5	W	link
Citadel	3	W-S	link
Clean Slate 2.0	3	W	link
Connect 2.0 ^L	2	W	link
Core 2.0	2	W	link

Stinson®	tier	care	link
Criss Cross 2.0	2	W-S	link
Discovery	3	W-S	link
English Pub ⁹	3	W-S	link
Essential	2	W-S	link
Fizz 2.0	4	W	link
Flanders	3	W-S	link
Flip ⁹	6	W	link
Frequency 2.0	5	W	link
Fuse 2.0	2	W-S	link
Gabby	2	W-S	link
GeoMet ^L	2	W	link
Grille 2.0	3	W	link
Hopsack	3	W	link
Kusari 2.0 ⁹	3	W	link
Lara	2	W	link
Laredo 2.0 ⁹	2	W-S	link
Lattice	2	W	link
Level 2.0	3	W	link

Stinson®	tier	care	link
Loophole 2.0 ^L	3	W	link
Lunaria ⁹	6	W	link
Matrix ^{9L}	2	W-S	link
Meteor	6	W	link
Mica 2.0 ⁹	4	W	link
Moccasin 2.0	3	W	link
Montana 2.0	3	W	link
Natural World 2.0	6	W-S	link
Naugasoft ^{9L}	3	W-S	link
Nishiki 2.0	3	W-S	link
Obi Stripe	6	W	link
Orbit 2.0	3	W	link
Oxford 2.0	3	W	link
Painted Moon	5	W-S	link
Pebble Trail	5	W	link
Ping Pong	6	W	link
Pod	6	W	link
Polygon	6	W	link

Stinson®	tier	care	link
Promenade 2.0 ^L	5	W	link
Raffia 2.0	5	W	link
Rally 2.0 ⁹	2	W-S	link
Rialto ^{9L}	3	W-S	link
Riverstone	3	W	link
Rogue II ⁹	2	W-S	link
Round Trip	6	W	link
Route	2	W	link
Runway	3	W-S	link
Sand Wave	9	W-S	link
Satchel	3	W-S	link
Scout 2.0	3	W-S	link
Script 2.0	4	W	link
Sea of Clouds	5	W-S	link
Sentiero	3	W-S	link
Sentry-HC 2.0	3	W	link
Shantung 2.0	3	W	link
Sirocco	6	W	link

Stinson®	tier	care	link
Sleek 2.0 ⁹	3	W-S	link
Smoothie 2.0	3	W-S	link
Stat Con ⁹	6	W	link
Stitch in Time ^{9L}	1	W-S	link
Strand 2.0 ⁹	3	W-S	link
Stride	3	W	link
Sumi	6	W-S	link
Tendril 2.0 ^L	2	W	link
Tessuto 2.0	3	W	link
Tilt	2	W	link
Tinker	6	W	link
Topography	5	W-S	link
Touchstone	3	W-S	link
Vector	4	W	link
Walk the Line 2.0	6	W	link
Watercolor Grid	6	W	link
Weft	2	W	link
Willow Wave	6	W-S	link

Non-PVC ^F	tier	care	link
Arioso ^{9L}	7	W	link
Attaché ⁹	3	W	link
Bond ^{9L}	6	W	link
Chaos ⁹	10	W	link
Composite ⁹	10	W-S	link
Corium ^{9L}	7	W-S	link
Didi ⁹	6	W-S	link
Effervescence 2.0 ^{9L}	8	W-S	link
Etch ⁹	8	W	link
Fiberdrill ^{9L}	5	W	link
Fiddlesticks ⁹	10	W-S	link
Flare ⁹	10	W-S	link
Galileo ⁹	7	W-S	link
Impress ^{9L} w/StaKleen™	5	W-S	link
Jax ⁹	8	W	link
Katsura ⁹	7	W-S	link
Kid ⁹	5	W-S	link

Non-PVC ^F	tier	care	link
Mesa ^{9L}	C	W-S	link
Modena 2.0 ^{9L}	7	W-S	link
Petal Dots ⁹	8	W-S	link
Rain Pool ⁹	8	W-S	link
Savile ⁹	7	W-S	link
Sedona ^{9L}	C	W-S	link
Silicone-Avail ⁹	10	W-S	link
Silicone-Downtown ^{9L}	10	W-S	link
Silicone-Uptown ^{9L}	10	W-S	link
Stage ⁹	5	W-S	link
Stallion ⁹	6	W-S	link
Staple ⁹	2	W	link
Sundance ^{9L} w/StaKleen™	5	W-S	link
Tender ^{9L}	C	W-S	link
Terra ⁹	8	W	link
Volterra-HC ^{9L}	C	W-S	link
Wyatt ^{9L}	3	W-S	link

Patty Madden ^F	tier	care	link
Affinity 2.0 ^L	6	W	link
Bari 2.0	6	W	link
Bubblerrake 2.0	9	W	link
Cassiopeia ⁹	9	W-S	link
Eli 2.0 ^L	3	W	link
Hashtag	5	W	link
Intersect	9	W	link
Juxy 2.0	6	W	link
La Scalina ⁹	8	W	link
Lava 2.0	9	W	link
Lens ⁹	9	W-S	link
Live Wire	6	W	link
Martini 2.0	6	W	link
Median	6	W	link
Odeon	9	W	link
Olympus 2.0 ^L	6	W	link
Origami	9	W	link

Patty Madden ^F	tier	care	link
Orvilla ⁹	9	W	link
Parasail	8	W-S	link
Perspective 2.0 ^L	6	W	link
Prisma ^L	3	W	link
Raya ^L	6	W	link
Revolve 2.0	9	W	link
Ripple 2.0	9	W	link
Rivera ^{9L}	3	W	link
Scoprena	9	W	link
Swizzle	8	W-S	link
Twist	9	W	link
Verge	9	W	link
Vortexia 2.0	9	W	link
Waterwand 2.0 ^L	6	W	link
Wicker 2.0	6	W	link
Winnily ^L	5	W	link
Ybor ^L	3	W	link

Michael Graves ^F	tier	care	link
Entrada ^{9L}	2	W	link
Excursion ^{9L}	2	W	link
Latitude ^{9L}	2	W	link
Valise ^{9L}	2	W-S	link

 CF Stinson Portfolio
Approved Vinyl Upholstery

Spradling®	tier	care
Guardian	2	S
Savoy	2	S
Cassara	2	S
Modena	1	S
Masialas	2	S
Dolce ⁹	2	S
Patriot Plus	2	S
Renaissance	2	S

Mayer Fabrics	tier	care	link
Belmont	4	W-S	link
Caressa	3	CFC	link
Catwalk	4	CFC	link
Chaps	4	W-S	link
Cirque	4	CFC	link
Etch	2	CFC	link
Explorer	2	W-S	link
Galveston	6	CFC	link

Mayer Fabrics	tier	care	link
Icon	3	W-S	link
Imprint	2	CFC	link
Interlace	4	CFC	link
Key Largo	2	W-S	link
Kinsey	2	W-S	link
Ludlow	2	W-S	link
Lyric	2	W	link
Matka	3	W-S	link

Mayer Fabrics	tier	care	link
Maverick	2	W-S	link
Midway	2	W	link
Mingle	3	W-S	link
Mode	4	CFC	link
Noble	2	CFC	link
Outback	2	W-S	link
Ryder	2	W-S	link
Silverweave	4	CFC	link

Mayer Fabrics	tier	care	link
Sketch	3	W-S	link
Sonnet	4	W-S	link
Spartan	2	W-S	link
Truss	4	W-S	link
Verona	3	W-S	link
Zip	2	CFC	link

Intensa®	tier	care	link
Standard Vinyl ⁹	1	D	—
Venue ⁹	2	D	link
Vista	2	D	link

 link to pattern swatch, spec, and care ^F by CF Stinson ^L limited availability on certain colors ⁹ product is not CFFA-HealthCare ^{tier C} Call for pricing



**RECOMMENDED****Care & Cleaning**

tier 1 standard vinyl - industry standard colors

**Regular Cleaning & Maintenance**

Clean the soiled area with mild soap and water, then rinse with fresh water and wipe dry with a clean cloth.

**Food Stains / Oils**

Wipe affected area with a soft cloth or soft bristle brush using a non-abrasive cleaner, for example Formula 409® or Fantastik®. Rinse with fresh water and wipe dry.

**Difficult Stains**

Clean the soiled area with 50% a solution of isopropyl Alcohol and water. It is important to rinse the cleaned area with fresh water after applying the alcohol solution.

**Healthcare / Disinfection**

Clean the soiled area with a 10% solution of Household Bleach and water. It is important to rinse the cleaned area with fresh water after applying the bleach solution.

Example

Dirt
Dust
Grime

Ketchup
Chocolate
Coffee
Tea
Juice

Eyeshadow/
Mascara
Lipstick
Crayon
Grease
Mustard

Blood
Urine
Betadine

Disinfectants

In addition to cleaning, healthcare facilities maintain disinfecting regimens to reduce the spread of infections. There is a wide variety of available biocides that include bleaches, peroxides and quaternary ammonium products, among others. They should always be used at the recommended dilution, never in concentrated form. To prolong the life of coated fabrics, the disinfectants should be rinsed off after each use. It should be noted that misuse of disinfectants is the major source of surface deterioration.

For Color Transfer

Certain clothing and accessory dyes (such as those used on denim jeans) may migrate to lighter colors. This phenomenon is increased by humidity and temperature and is irreversible. Neither the manufacturer nor GF Health Products, Inc. is responsible for any damages caused in this manner. Please check compatibility when using this product in combination with painted or varnished surfaces.

Note: The cleaning products recommended were used in a laboratory setting. This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and referenced cleaning agents. Neither the manufacturer nor GF Health Products, Inc. has a relationship with the manufacturer of these cleaning products.

**Stinson
MAYER****RECOMMENDED****Care & Cleaning Codes**

use chart links on previous page for full details

W

Clean with water-based foam or water-based cleaners only. Do not use solvents. We recommend using a professional cleaning service.

Bleach solutions of up to 10% may be used for more difficult stains.

W-S

For water-based stains use water-based cleaners. For oil-based stains use solvent-based cleaners. We recommend using a professional cleaning service. Bleach solutions of up to 20% may be used for more difficult stains.

Crypton® Cleaners are recommended to clean many of Stinson's fabrics. These specifically developed cleaners are water-based, pH neutral, biodegradable, easy to use, and do not contain bleach, solvents or harsh chemicals. [Crypton® Cleaners](#)

S

Clean the whole surface in a circular motion using pH neutral liquid hand soap and water (1 part soap to 9 parts water). Wipe with a clean damp cloth to get rid of any excess soap. Repeat as needed and then let dry.

D

After disinfection, hand dish-washing soap and water or Formula 409® may be used with a soft cloth or damp sponge to remove chemical residue. After cleaning, rinse area with clean water, then dry with a clean, lint-free cloth. Water rinse recommended after each disinfectant application.

CFC

Clean with mild soap and water or Fantastik® or Formula 409® may be used with a soft cloth. More stubborn stains may be eliminated using a cotton cloth soaked in 70-91% of isopropyl alcohol. Wipe without spreading stain to surrounding area. Repeat with vigorous rubbing in a circular motion. For textured material, a soft brush is recommended. Rinse with fresh water and wipe dry.

Note: The cleaning products recommended were used in a laboratory setting. This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and referenced cleaning agents. Neither CF Stinson LLC, Mayer Fabrics, nor GF Health Products, Inc. has a relationship with the manufacturer of these cleaning products.

**SPRADLING®****RECOMMENDED****Care & Cleaning**

Clean within 24 hours to avoid permanent stains in the surface. After this time has elapsed, effectiveness drops to 50%.



We recommend cleaning at least once a week to extend the lifespan of the upholstery, which will prevent dirt and contamination buildup and permanent staining. In the event of a spill, wipe off any excess right away with a dry white cloth or fabric, but careful not to spread to remaining surface.



Clean the whole surface in a circular motion using pH neutral liquid hand soap and water (1 part soap to 9 parts water). Check the Approved Cleaners guide for recommended cleaning products. Wipe with a clean damp cloth to get rid of any excess soap. Repeat as needed and then let dry.



If there is any buildup that cannot be removed with soapy water, a mixture of alcohol and water (1 part isopropyl alcohol to 9 parts water) can be used as a last resort. Remember, using alcohol on a regular basis can affect the properties of the product.



Once a month, you can use a soft-bristle brush for a deeper cleaning. Remember, light colors need more care, more often.

Some products have more detailed procedures and/or restrictions; see specifics in below guides:

[Cleaning Guides specific for each finish.](#) [www.spradling.group]

[Guide to Healthcare Surface Disinfection for the Real World.](#) [www.spradling.group]

Disinfectants

There is a wide variety of available biocides that include bleaches, peroxides and quaternary ammonium products, among others. They should always be used at the recommended dilution, never in concentrated form. To prolong the life of coated fabrics, the disinfectants should be rinsed off after each use.

If you need to disinfect the product because it has come into contact with organic material (blood, urine or fecal matter), check the [disinfection guide](#) for appropriate cleaning products. Clean with a clean white cloth or fabric. If you choose to use a bleach solution, only use 0.5% concentration (1 part bleach to 10 parts water). Rinse and wipe with a clean damp cloth to get rid of any excess and dry.

For Color Transfer

Our products don't transfer onto other materials, but some clothing items (like jeans, microfibers and dark fabrics) could transfer their color to the coated fabrics, depending on how they are dry-cleaned. This happens more with dampness and higher temperatures and is irreversible. Neither GF Health Products, Inc. nor Spradling International, Inc. is responsible for any damages caused in this manner. Please check compatibility when using this product in combination with painted or varnished surfaces.

Note: The cleaning products recommended were used in a laboratory setting. This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and referenced cleaning agents. Neither Spradling International, Inc. nor GF Health Products, Inc. has a relationship with the manufacturer of these cleaning products.

f this section of cleaning guide does not apply to Dolce fabric



Composition & Specifications

PVC	■
PU	—
TPE	—
Polyester Backing (100%)	■
Polyester Backing Hi-Loft ² ™ (100%)	—
Polyester Backing Ultra Hi-Loft ² ™ (100%)	—
Polyester Backing Hi-Loft ² ™ RPET (100%)	—
Polyester / Cotton Backing	—
Polyester / Rayon Backing	—
Weight (oz/liner yd)	30
Roll Size (yd)	30
Roll Width (inch)	54

Stinson

	Olympus	Natural World	Raffia	Tessuto
✓	✓	✓	✓	✓
—	—	—	—	—
—	—	—	—	—
✓	✓	✓	✓	✓
—	—	—	—	—
—	—	—	—	—
—	—	—	—	—
—	—	—	—	—
32	30	30	30	30
30	30	30	30	30
54	54	54	54	54

MAYER

	Galveston	Belmont	Icon	Midway
—	—	—	—	—
■	■	■	■	■
—	—	—	—	—
■	■	■	■	■
—	—	—	—	—
—	—	—	—	—
—	—	—	—	—
30	29	31	28	
30	30	30	30	30
54	54	54	54	54

SPRADLING®

	NEO		eco ^{wise} LIVING					
	Dolce	Patriot Plus	Renaissance	Guardian	Savoy	Cassara	Modena	Masiala
—	—	■	■	■	■	—	■	■
■	■	—	—	—	—	—	—	—
—	—	—	—	—	—	■	—	—
—	—	—	—	—	—	■	—	—
—	—	■	—	—	■	—	■	—
—	—	—	—	■	—	—	—	—
—	—	—	—	—	—	—	—	■
■	■	—	—	—	—	—	—	—
15	34	32	30	30	22	26	31	
30	30	30	30	30	30	30	30	
54	54	54	54	54	54	54	54	

Flammability Codes^T

NFPA 260 Cover Fabric	Class 1	■
UFAC Fabric	Class 1	■
CAL TB 117-2013		■
FMVSS 302		—
EU/US: IMO	MSC.307 (88) (FTP2010)	—
	BIFMA	—
UK: BS 5852	Source 0	—
Does not contain flame retardants		■

OLY	NTW	RAF	TES
			
			
			
			
			
			
			
			

GV	BL	IC	MW
■	■	■	■
■	■	■	■
■	■	■	■
■	■	■	—
—	—	—	—
—	—	—	■
—	—	—	—
■	■	■	■

	DOL	PAT	REN	GUA	SAV	CAS	MOD	MAI
—		■	■	■	■	■	■	■
—		■	■	■	■	■	■	■
■		■	■	■	■	■	■	■
—		■	—	—	—	■	—	—
—	—	—	—	—	—	—	—	■
—	—	—	—	—	—	—	—	—
DOL	—	—	—	—	—	—	—	—
■	■	■	■	■	■	■	■	■

Finishes & Top Coats

	—
	—
	—
	—
Advanced Soil & Stain Resistant Top Coat	—
Zeron® Protective Finish	—

OLY	NTW	RAF	TES
—	—	—	—
—	—	—	—
—	—	—	—
—	—	—	—
■	■	■	■
■	—	—	—

GV	BL	IC	MW
—	—	—	—
—	—	—	—
—	—	—	—
—	—	—	—
■	■	■	■
—	—	—	—

DOL	PAT	REN	GUA	SAV	CAS	MOD	MAI
—	—	—	■	■	—	—	—
—	■	■	—	—	—	—	■
—	—	—	—	—	—	■	—
—	—	—	—	—	■	—	—
—	■	■	■	■	■	■	■
—	—	—	—	—	—	—	—

Performance & Environmental

Abrasion Resistance ^M	ASTM D-4157	■
Antibacterial Protection		—
Anti-Static	ASTM D257	—
Colorfast to Light		—
Colorfast to UV	QUV / ASTM G154-16 (hours)	1000
Crocking Resistance	AATCC 8-2007 [wet / dry] (k)	5/5
Mildew Resistance	ASTM G21	■
Stain Resistance		■
Stain Resistance	in Healthcare Environment CFFA142	—
Sulfide Stain Resistance	Fisher Body TM31-12	■
Denim Stain Resistance		—
Bleach Resistance	10% Bleach Solution Cleanable	■
Ink Resistant		—
CFFA-Healthcare-201: Certified		—
Greenguard/Greenguard Gold: Certified		—
SCS Global Indoor Advantage™	Gold	—
GRS – Global Recycled Standard		—
Hydrolytic Stability ISO 1419	Tropical Test Method C (7+ weeks)	—
REACH		■
RoHS		■
Wyzenbeek Cotton Test	(k)	250
Cold Crack	CFFA-6A-2007 (F)	-20°

OLY	NTW	RAF	TES
	—	—	
—	—	—	—
—	—	—	—
1000	1000	1000	1000
5/5	5/5	5/5	5/5
—	—	—	—
—	—	—	—
	—		
	—	—	—
—	—	—	—
—	—	—	—
—	—	—	—
—	—	—	—

GV	BL	IC	MW
■	■	■	■
—	—	—	—
200	200	300	650
—	—	—	—
4/4	4/4	5/5	4/4
—	—	—	—
■	■	■	■
—	■	—	■
—	—	—	—
■	■	—	—
■	■	■	■
■	■	■	■
—	■	—	■
■	■	■	■
—	—	—	—
—	—	—	—
—	—	—	—
■	■	■	—
■	■	■	■
■	■	■	■
—	-20°	-20°	-20°

DOL	PAT	REN	GUA	SAV	CAS	MOD	MAI
■	■	■	■	■	■	■	■
—	—	—	■	—	—	—	—
—	■	■	—	■	■	■	■
—	—	—	—	—	—	—	—
650	—	—	650	—	—	—	—
—	—	—	4/4	—	—	—	—
—	—	■	■	—	—	■	—
■	■	■	■	■	■	■	■
—	■	■	■	■	■	■	■
—	■	■	■	■	—	■	■
—	■	■	■	■	■	■	—
■	■	■	■	■	—	■	—
—	—	—	—	—	—	■	—
—	■	■	■	■	■	■	■
—	—	—	—	—	■	—	■
—	■	■	■	■	■	■	■
—	—	■	—	—	—	—	■
—	—	—	—	—	—	—	—
■	■	■	■	■	■	■	■
#10	#10	#10	#10	#10	#10	#10	#10
-10°	-20°	-20°	-10°	-10°	-10°	-10°	-10°

^M multiple factors affect fabric durability and appearance retention, including end-user application and proper maintenance. Wyzenbeek results above 100,000 double rubs have not been shown to be an indicator of increased lifespan

Put a little
COLOR
 in your
 exam
 room
 with **intensa**
 by graham-field



Product Index

click listing to jump to page

990 SERIES Physician Exam Stools	4
980 SERIES Physician Exam Stools	5
970 SERIES Physician Exam Stools	6
960 SERIES Physician Exam Stools	7
950 SERIES Physician Exam Stools	8
930 SERIES Physician Exam Stools	9
920 SERIES Physician Exam Stools	8
870 SERIES Laboratory Seating	10
860 SERIES Laboratory Seating	11
840 SERIES Laboratory Seating	13
830 SERIES Laboratory Seating	16
810 SERIES Laboratory Seating	17
740 SERIES Hydraulic Surgical Seating	19
HSS500 SERIES Hydraulic Surgical Seating	19
IN601 SERIES Pivot-Arm Clinical Care Recliners	20
IN597 SERIES Ortho-Biotic™ II Clinical Care Recliners	21
IN587 SERIES Extra Wide Clinical Care Recliners	22
IN566 SERIES Deluxe Clinical Care Recliners	23
IN577 SERIES Multi-Care Recliners	24
600 SERIES Blood Draw / Phlebotomy Seating	24
460 SERIES Exam Room Treatment Tables	25
420 SERIES Exam Room Treatment Tables	25
410 SERIES Exam Room Treatment Tables	26
RT400 SERIES Replacement Exam Table Tops	28
RT300 SERIES Replacement Exam Table Tops	29
RT200 SERIES Replacement Exam Table Tops	29
RT100 SERIES Replacement Exam Table Tops	29
200 SERIES Wallsaver Patient / Guest Seating	30
170 SERIES Overbed Tables	31
100 SERIES Step Stools	32
500 SERIES Exam Light	33
130 SERIES IV Pole	33
Footnote Legend	33
Order Numbering Multiple Options & Codes	34
Customer's Material Requirements	34
Additional Information Notices & Links	34
Physician & Laboratory Seating Options	34
Clinical Care Recliner Options & Accessories	35
Overbed Table Finish Swatches & Care	36
Tier 1 Vinyl Upholstery Standard Collection Swatches	37
Stinson® Vinyl Upholstery Pattern Collection Swatches	38
Mayer Vinyl Upholstery Texture Collection Swatches	39
Spradling® Vinyl Upholstery High Performance, Green Collection Swatches	40
Spradling® Vinyl Upholstery Performance Collection Swatches	41
Spradling® Vinyl Upholstery Premiere Collection Swatches	41
Approved Vinyl Upholstery Tier, Care, Link, & Swatch Chart	42
Vinyl Upholstery Recommended Care, Cleaning & Disinfection	43
Vinyl Upholstery Composition, Finish, & Performance Comparison Chart	44





A Better Healing Environment Better Results



published **July 2026**

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

GF1800088 RevU26

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net



© 2015, GF Health Products, Inc. All Rights Reserved | Intensa, Inc. is a wholly owned subsidiary of GF Health Products, Inc. | Stinson and Zeron are registered trademarks of C. F. Stinson, LLC.
| Intensa, Dura-Gel, and Sok-It-Guard are registered trademarks of GF Health Products, Inc. | Spradling, Neo, EcoSense, Hi-Loft, PBG+, Permablok, Permablock Plus, and Sureguard are trademarks of Spradling Holdings, Inc.
| Crypton is a registered trademark of Crypton LLC | Tente is a registered trademark of TENTE International | Formula 409 is a registered trademark of The Clorox Company. | Fantastic is a registered trademark of S.C. Johnson & Son, Inc.