



ACCESSORIES

170 SERIES AUTOMATIC OVERBED TABLE

175 Features:

- Supports up to 100 lb of evenly distributed weight
- Automatic/infinite stops within height travel range
- Modified "H" base allows use with chair or bed
- Low profile design allows overbed table to be used with limited clearance beds
- Tops are customizable in both finish and edge

175 Product Info

175 Options:

| | |
|------------------------------|---------------|
| base and top sold separately | item number † |
| Flat Edge 3DL Top† | OS3DLOBTF-xx |
| Raised Edge 3DL Top† | OS3DLOBTR-xx |
| Modified "H" 100 lb Base† | A919000i |

‡ -xx in the product number is the placeholder for top finish

† base and top sold separately



175 Specifications:


| | inch | cm | | inch | cm |
|--------------------------------|-------|--------|---------------------|-------|-------|
| Overall Width | 32.00 | 81.28 | Base Width | 29.30 | 74.40 |
| Overall Depth | 22.00 | 55.88 | Base Depth | 18.20 | 46.20 |
| Height Range <i>minimum</i> | 29.90 | 75.90 | Base Profile Height | 1.94 | 4.93 |
| Height Range <i>maximum</i> | 43.00 | 109.20 | Top Width | 35.00 | 88.90 |
| | | | Top Depth | 18.00 | 45.72 |



| | lb | kg |
|---------------------------------------|-----|-------|
| Weight Capacity EVENLY DISTRIBUTED | 100 | 45.36 |
| Shipping Weight | 43 | 19.50 |

*Ships: Freight request additional shipping quote from Intensa ‡ -xx in the product number is the placeholder for finish choice † base and top sold separately

175 Top Finishes:



Aria | AR

Wild Cherry | WC

Windsor Mahogany
WM

Cocobala | CO

Shaker Cherry | SC

Solar Oak | SO

Talas Cherry | TC



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

