

ACCESSORIES

170 SERIES OVERBED TABLE

Overbed Table with Heavy Duty Steel, Low Profile Base and Laminated Top

This Overbed Table combines an attractive, ergonomic design with features that reduce infection control issues by minimizing areas in which dust and dirt can accumulate. Manufactured with the lowest base in the industry, our Overbed tables move easily under any hospital bed or recliner.


177 Features:

- Supports up to 100 lb of evenly distributed weight
- Exposed steel base with shock-absorbent gas spring cylinder
- Seamless spill guard top with optional integrated cup holders
- Laminate or Thermofoil tops are customizable in shape, finish, and edge
- Low profile design allows overbed table to be used with limited clearance beds

 177 Product Info



177 Options: ††

| steel base with top | item no. | whiteboard | jellybean | kidney | |
|--|---------------|-------------|----------------|----------------|---|
| Laminate | 177-y-xx | 177-WB-y-xx | — | — | |
| Thermofoil Shape | 177TFR-y-xx | — | 177TFJ-y-xx | 177TFK-y-xx | |
| Thermofoil Shape <i>with Cup Holder</i> | 177TF-CH-y-xx | — | 177TFJ-CH-y-xx | 177TFK-CH-y-xx |  |

options available at additional cost

177 Specifications:

| | inch | cm | | inch | cm |
|----------------------------------|-------|--------|---------------------|-------|-------|
| Overall Width | 32.00 | 81.28 | Base Width | 32.00 | 81.28 |
| Overall Depth | 18.00 | 45.72 | Base Depth | 18.00 | 45.72 |
| Overall Height <i>minimum</i> | 29.00 | 73.66 | Base Profile Height | 1.94 | 4.93 |
| Overall Height <i>maximum</i> | 44.00 | 111.76 | Top Width | 32.00 | 81.28 |
| | | | Top Depth | 17.00 | 43.18 |



| | lb | kg |
|---------------------------------------|-----|-------|
| Weight Capacity EVENLY DISTRIBUTED | 100 | 45.36 |
| Shipping Weight | 41 | 18.60 |

177 Accessories:

| | item number |
|--|-------------|
| Vanity w/Mirror | 177-VM |
| Basket 10 x 6 x 9 inch (25.40 x 15.24 x 22.86 cm) | 177-BK |

accessories available at additional cost

170 OBT Top Finishes (xx)



Candlelight | CL Chocolate Pear | CP Hard Rock Maple | HM Medium Cherry | MC Medium Oak | MO

170 Paint Finishes (y)




Almond | A Graphite Gray | G Nevada | N



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in Canada

