



TREATMENT ACCESSORIES

130 SERIES IV POLE STAND

130 Features:

- Chrome Plated Steel Tubing
- Ram Style Hooks
- Assembly Required

Product Info:

130-2 134-2 134-4

130 Options:


	item no.
Chrome Plated Steel IV Pole <i>4 Leg Steel Base, Rubber Caps and 2 Hooks</i>	130-2
IV Pole Stand <i>5 Leg Flat Aluminum Base, Toe Caps and 2 Hooks</i>	134-2
IV Pole Stand <i>5 Leg Flat Aluminum Base, Toe Caps and 4 Hooks</i>	134-4

130 Specifications:

	inch	cm		inch	cm		
Overall Height <i>-low</i>	51.00	129.54	Casters	2.00	5.08	Packages Shipped	2
Overall Height <i>-high</i>	90.00	228.6	Steel Tubing	1.00	2.54	Ships Assembled	no
			Steel Tubing Finish	chrome plated			
			Hooks	ram style		Shipping Weight	14 lb / 6.35 kg

Adjustable Height IV Pole[‡]		inch	cm
Twist Lock	130-x	Steel Base Diameter	22.00 / 55.88
Turn Knob Lock	134-x	Aluminum Base Diameter	24.00 / 60.96


[‡] Ships: Ground at no additional charge † -x in the product number is the placeholder for number of hooks



130-2



134-2



134-4



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

Color variations may occur due to monitor or printer settings.

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in China

