

## EXAM TABLES

### 460 SERIES

# HIGH/LOW POWER EXAM TABLE

#### High/Low Power Exam Table with Power Assisted Back, Hand Control and Stirrups

A power table with power assisted back, storage drawer, and 5 position stirrups. All metal construction of the base provides for longevity. Intensa offers an extensive collection of vinyls to add color to the exam room.

#### 460 Features:

- Dual Pedestal Base with DC Electromechanical Drive
- Pull out Leg Rest with Integrated Treatment Tray
- 5 Position Continuous-Lock Stirrups
- Powder Coat Finish Foghorn
- Power Assisted Back
- Pelvic Tilt Standard

 CF Stinson Pre-Approved Fabrics

 460 Product Info

#### 460 Options: †

	item no.
Power Exam Table	460-xx
Power Exam Table <i>with Foot Switch</i>	460-FS-xx
Power Exam Table <i>with Assist Arms</i>	460-AA-xx
Power Exam Table <i>with Physician Storage only</i>	460-PS-xx
Power Exam Table <i>with Physician Storage &amp; Drawer Warmer</i>	460-PS-DW-xx
Power Exam Table <i>Fully Loaded</i>	460-FL-xx

options available at additional cost



#### 460 Specifications:

	inch	cm				
Overall Width	27.00	68.58	Pelvic Tilt	standard	Power Assisted	back
Overall Length <i>leg rest extended</i>	72.00	182.88	Pull Out Leg Rest	standard	UL 60601-1	yes
Overall Height <i>minimum</i>	18.50 ± 0.5	46.99 ± 1.27	Treatment Tray	integrated		yd m
Overall Height <i>maximum</i>	35.50 ± 0.5	90.17 ± 1.27	Stirrups	five position continuous lock	COM Requirement <sup>1</sup>	3 2.74
Paper Roll Holder Widths Available	18.00 19.00 20.00	45.72 48.26 50.80	Base	dual pedestal		lb kg
			Drive	DC electromechanical	Weight Capacity EVENLY DISTRIBUTED	600 272.16
			Powder Coat Finish	foghorn	Shipping Weight	395 179.7

\*Ships: Freight request additional shipping quote from Intensa †-xx in the product number is the placeholder for top finish choice  
<sup>1</sup>Submit COM sample for approval



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

