

EXAM & TREATMENT TABLES

RT204 SERIES

REPLACEMENT TOP

Replacement Top with Foot Pad for Ritter® by Midmark® 204 Exam Tables

Intensa offers a replacement top program for your existing exam tables. We provide tops that meet certain manufacture model number specifications. The new top is made of completely new foam and wood substrate and can be upholstered in a variety of vinyl colors and patterns. Revitalize/renew your exam rooms without the expense of new tables!

RT204 Features:

- Replaces Ritter® by Midmark® 204 Exam Table
- Foot pad included

 [CF Stinson Pre-Approved Fabrics](#)

 [Replacement Tops Overview](#)


 [RT204 Product Info](#)

RT204 Options:

	Item no.†
Replacement Top ^{2,3}	RT204-xx



RT204 Specifications: ^{2,3}

	inch	cm		inch	cm		yd	m
Overall Width	27.50	69.85	Foot Rest Width	20.00	50.80	COM Requirement ¹	3	2.74
Overall Length	55.00	139.70	Foot Rest Length	18.00	45.72			
Lift Handles ²	not included		Top Attachment Knobs ²	not included		Shipping Weight	55	24.95


^{*} Ships: Freight request additional shipping quote from Intensa [†] -xx in the product number is the placeholder for vinyl choice ¹ Submit COM sample for approval

² Lift Handles and Top Attachment Knobs for exam table year models 2009 and after are not included ³ If OEM model is seamless top, request paper roll holder when ordering



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

