

EXAM & TREATMENT TABLES

RT460 SERIES REPLACEMENT TOP

Replacement Top with Foot Pad for Intensa 460 Exam Tables

Intensa offers a replacement top program for your existing exam tables. We provide tops that meet certain manufacture model number specifications. The new top is made of completely new foam and wood substrate and can be upholstered in a variety of vinyl colors and patterns. Revitalize/renew your exam rooms without the expense of new tables!

RT460 Features:

- Replaces Intensa 460 Exam Table
- Paper roll holder and foot pad included
- Includes mounting hardware

 [CF Stinson Pre-Approved Fabrics](#)

 [Replacement Tops Overview](#)


 [RT460 Product Info](#)

RT460 Options:

	Item no. †
Replacement Top	RT460-xx



RT460 Specifications:


	inch	cm		inch	cm		yd	m
Overall Width	27.50	69.85	Foot Rest Width	19.00	48.26	COM Requirement ¹ <i>no railroad patterns</i>	3	2.74
Overall Length	60.00	152.40	Foot Rest Length	18.00	45.72	 + FREIGHT *	lb	kg
Paper Roll Holder	included		Mounting Hardware	included		Shipping Weight	55	24.95

* Ships: Freight request additional shipping quote from Intensa † -xx in the product number is the placeholder for vinyl choice ¹ Submit COM sample for approval



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

