

## LABORATORY

# 830 SERIES LABORATORY SEATING

### 830 Features:

- Single lever adjustment release with tilt
- Pneumatic height adjustment
- Standard back rest with vertical adjustment
- Standard 20 inch (50.80 cm) footring

 **CF Stinson Pre-Approved Fabrics**

 **830 Options Detail and Photos**

Product Info:


 831  832  833  838

### 830 Options:

base with standard seat height	item number †	with arm rest ‡
<b>Black Nylon Composite</b>	831-xx	831-T-xx
<b>Brushed Aluminum<sup>2</sup></b> and Toe Caps	832-xx	832-T-xx
<b>Polished Aluminum</b>	833-xx	833-T-xx
<b>Five Leg Tubular Steel<sup>4</sup></b> Footring and Glides	838-xx	838-T-xx

options available at additional cost

### 830 Specifications:

	inch	cm		inch	cm		yd	m
<b>Seat Width</b>	20.00	50.80	<b>Back Width</b>	19.00	48.26	<b>COM Requirement<sup>1</sup></b>	1	0.92
<b>Seat Depth</b>	18.00	45.72	<b>Back Height</b>	16.00	40.64		<b>lb</b>	<b>kg</b>
<b>Overall Height</b> with backrest -low	39.00	99.06	<b>Base Diameter</b> with casters	24.00	60.96		<b>Weight Capacity</b> EVENLY DISTRIBUTED	350
<b>Overall Height</b> with backrest -high	49.50	125.73	<b>Height Adjustment</b> method	pneumatic		<b>Shipping Weight</b>	45	20.41
<b>Height Release</b> hardware	lever		<b>Back Adjustment</b> direction	height				


### Adjustable Height Seat<sup>††</sup>

	inch	cm		inch	cm
<b>Standard</b>	83x-xx		<b>Seat Height -low</b>	24.00	60.96
<b>7 [seven]</b>	83x-7-xx		<b>Seat Height -high</b>	34.50	87.63
			<b>Seat Height -low</b>	20.00	50.80
			<b>Seat Height -high</b>	27.00	68.58

‡ Ships: Ground at no additional charge † -xx in the product number is the placeholder for vinyl choice † Replacement Cylinders may be ordered separately at additional cost † Submit COM sample for approval ‡ Adjust Height Specs: -1" (-2.54 cm) † Casters add 2" (5.08 cm) in height

 GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.



Manufactured in the USA

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA  
+1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

