

LABORATORY

840SC SERIES LABORATORY SEATING

All 840 series chairs are designed for the harshest environments. They have an all-steel inner frame, are wrapped in self-skin urethane, and are impervious to harsh chemicals.

840SC Features:

- Single lever adjustment release
- Pneumatic height adjustment
- Standard back rest with vertical adjustment
- Impact-resistant
- Standard 20 inch footring
- High density fire resistant foam impervious to harsh chemicals
- Available in black urethane only

 [840SC Options Detail and Photos](#)

Product Info:


 [841SC](#)  [842SC](#)  [843SC](#)  [848SC](#)

840SC Options:

base with standard seat height	item no [†]	with arm rest [†]
Black Nylon Composite	841SC-xx	841SC-T-xx
Brushed Aluminum <small>and Toe Caps</small>	842SC-xx	842SC-T-xx
Polished Aluminum	843SC-xx	843SC-T-xx
Five Leg Tubular Steel⁴ <small>Footring and Glides</small>	878SC-xx	878SC-T-xx

options available at additional cost

840 Specifications:

	inch	cm		inch	cm		
Seat Width	18.00	45.72	Back Width	16.00	40.64	Seat & Back substrate	steel
Seat Depth	17.00	43.18	Back Height	12.00	30.48	Seat & Back cover material	urethane
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96		lb kg
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic		Weight Capacity EVENLY DISTRIBUTED	350 158.76
Height Release hardware	lever		Back Adjustment direction	lateral		Shipping Weight	38 17.24

Adjustable Height Seat^{††}

	inch	cm		inch	cm
Standard	84xSC-xx		Seat Height -low	22.50	57.15
7 [seven]	84xSC-7-xx		Seat Height -high	33.00	83.82
			Seat Height -low	18.50	46.99
			Seat Height -high	25.50	64.77


[†] -xx in the product number is the placeholder for vinyl choice ^{††} Replacement Cylinders may be ordered separately at additional cost 
[♦] Ships: Ground at no additional charge [†] Submit COM sample for approval
⁴ 20 inch (50.80cm) footring; casters add 2" (5.08 cm) in height

Item	height		ship weight	
	inch	cm	lb	kg
841	—	—	-1	-0.45
842	-1	-2.54	—	—
848	+2	+5.08	+4	+1.81



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Manufactured in the USA

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

