

## LABORATORY

# 860 SERIES LABORATORY SEATING

### 860 Features:

- Single lever adjustment release
- Pneumatic height adjustment
- Lateral back rest adjustment
- Standard 18 inch footring<sup>3</sup>
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 **CF Stinson Pre-Approved Fabrics**

 **860 Options Detail and Photos**

Product Info:


 **861**  **862**  **863**  **868**

### 860 Options:

base with standard backrest & seat height	item no <sup>1</sup>	crescent back <sup>2</sup>
<b>Black Nylon Composite</b>	861-xx	861CB-xx
<b>Brushed Aluminum</b> <small>and Toe Caps</small>	862-xx	862CB-xx
<b>Polished Aluminum</b>	863-xx	863CB-xx
<b>Five Leg Tubular Steel<sup>4</sup></b> <small>Footring and Glides</small>	868-xx	868CB-xx

options available at additional cost

### 860 Specifications:

	inch	cm		inch	cm		yd	m
<b>Seat Diameter</b>	16.00	40.64	<b>Back Width</b>	9.75	24.77	<b>COM Requirement<sup>1</sup></b>	1	0.92
<b>Overall Height</b> <small>with backrest -low</small>	37.00	93.98	<b>Back Height</b>	15.75	45.72			
<b>Overall Height</b> <small>with backrest -high</small>	47.00	119.38	<b>Base Diameter</b> <small>with glides</small>	24.00	60.96		<b>Weight Capacity</b> EVENLY DISTRIBUTED	350
<b>Seat Foam</b> <small>thickness</small>	3.00	7.62	<b>Height Adjustment</b> <small>method</small>		pneumatic	<b>Shipping Weight</b>	32	14.51
<b>Height Release</b> <small>hardware</small>		lever	<b>Back Adjustment</b> <small>direction</small>		lateral			

### Adjustable Height Seat<sup>++</sup>

	inch	cm		inch	cm
<b>Standard</b>	86x-xx		<b>Seat Height -low</b>	24.50	62.23
<b>7 [seven]</b>	86x-7-xx		<b>Seat Height -high</b>	35.00	88.90
			<b>Seat Height -low</b>	19.50	49.53
			<b>Seat Height -high</b>	27.50	69.85

<sup>+</sup> -xx in the product number is the placeholder for vinyl choice

<sup>†</sup> Replacement Cylinders may be ordered separately at additional cost

<sup>1</sup> Submit COM sample for approval

<sup>4</sup> 20 inch (50.80cm) footring; casters add 2" (5.08 cm) in height

<sup>3</sup> Ships: Ground at no additional charge


### <sup>2</sup> Specification Variances

item	height		ship weight	
	inch	cm	lb	kg
861	—	—	-1	-0.45
862	-1	-2.54	—	—
868	+2	+5.08	+4	+1.81



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Manufactured in the USA

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA  
+1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

