

LABORATORY

920 SERIES LABORATORY SEATING

920 Features:

- Single lever adjustment release
- Pneumatic height adjustment
- Available in black urethane only
- Standard 20 inch (50.80 cm) footring

 [920 Options Detail and Photos](#)

 [928 Product Info](#)




920 Options:

base with standard seat height	item no. [†]	procedure arm [‡]
Five Leg Tubular Steel⁴ <i>Footring and Glides</i>	928-xx	928PA-xx

options available at additional cost

920 Specifications:

	inch	cm		inch	cm		lb	kg
Seat Diameter	14.00	35.56	Base Diameter with glides	24.00	60.96			
Overall Height -low	15.00	38.10	Height Adjustment method		pneumatic	Weight Capacity EVENLY DISTRIBUTED	350	158.76
Overall Height -high	19.75	50.17	Height Release hardware		lever	Shipping Weight	23	10.43

Adjustable Height Seat^{**}

	inch	cm		inch	cm		
Standard	928-xx	Seat Height -low	15.00	38.10	Seat Height -high	19.75	50.17
7 [seven]	928-7-xx	Seat Height -low	16.50	41.91	Seat Height -high	23.25	59.06
10 [ten]	928-10-xx	Seat Height -low	20.50	52.07	Seat Height -high	31.00	78.74

^{**}-xx in the product number is the placeholder for vinyl choice

[‡] Ships: Ground at no additional charge

[†] Replacement Cylinders may be ordered separately at additional cost


[‡] Submit COM sample for approval

⁴ Casters add 2" (5.08 cm) in height



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Manufactured in the USA

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

