

SEATING

930 SERIES PHYSICIAN EXAM STOOL

930 Features:

- Manual spin lift adjustment
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 CF Stinson Pre-Approved Fabrics

 930 Options Detail and Photos

Product Info:

 931  932  933

930 Options:

base with standard seat height	item no.†	with backrest‡	crescent back‡
Black Nylon Composite	931-xx	931B-xx	931CB-xx
Brushed Aluminum² <small>and Toe Caps</small>	932-xx	932B-xx	932CB-xx
Polished Aluminum	933-xx	933B-xx	933CB-xx

options available at additional cost



931




932



933

930 Specifications:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter <small>with casters</small>	24.00	60.96	COM Requirement[†]	1	0.92
Overall Height <small>with backrest -low</small>	33.00	83.82	Height Adjustment <small>method</small>	manual spin lift			lb	kg
Overall Height <small>with backrest -high</small>	39.00	99.06	Height Release <small>hardware</small>	none		Weight Capacity <small>EVENLY DISTRIBUTED</small>	350	158.76
						Shipping Weight	20	9.07

Adjustable Height Seat^{††}


	inch	cm		inch	cm
Standard	96x-xx		Seat Height -low	19.50	49.53
			Seat Height -high	23.00	58.42

†-xx in the product number is the placeholder for vinyl choice † Replacement Cylinders may be ordered separately at additional cost
 ‡ Ships: Ground at no additional charge † Submit COM sample for approval ‡ Adjust Height Specs: -1" (-2.54 cm)



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

