

SEATING

960 SERIES PHYSICIAN EXAM STOOL

960 Features:

- Single lever adjustment release
- Pneumatic height adjustment
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 **CF Stinson Pre-Approved Fabrics**

 **960 Options Detail and Photos**

Product Info:

 **961**  **962**  **963**




960 Options:

| base with standard seat height | item no. [†] | procedure arm [‡] | with footring [‡] | with backrest [‡] | crescent back [‡] |
|--|-----------------------|----------------------------|----------------------------|----------------------------|----------------------------|
| Black Nylon Composite | 961-xx | 961PA-xx | 961-FR-xx | 961B-xx | 961CB-xx |
| Brushed Aluminum² <i>and Toe Caps</i> | 962-xx | 962PA-xx | 962-FR-xx | 962B-xx | 962CB-xx |
| Polished Aluminum | 963-xx | 963PA-xx | 963-FR-xx | 963B-xx | 963CB-xx |

options available at additional cost

960 Specifications:

| | inch | cm | | inch | cm | | yd | m |
|---|-------|-------|-----------------------------------|-----------|-------|---|-----|--------|
| Seat Diameter | 16.00 | 40.64 | Base Diameter with casters | 24.00 | 60.96 | COM Requirement¹ | 1 | 0.92 |
| Overall Height with backrest -low | 33.50 | 85.09 | Height Adjustment method | pneumatic | |  | lb | kg |
| Overall Height with backrest -high | 38.50 | 97.79 | Height Release hardware | lever | | Weight Capacity EVENLY DISTRIBUTED | 350 | 158.76 |
| | | | | | | Shipping Weight | 21 | 53.34 |

Adjustable Height Seat^{††}


| | inch | cm | | inch | cm |
|------------------|-----------|----|--------------------------|-------|-------|
| Standard | 96x-xx | | Seat Height -low | 18.50 | 46.99 |
| 7 [seven] | 96x-7-xx | | Seat Height -low | 20.50 | 52.07 |
| 10 [ten] | 96x-10-xx | | Seat Height -low | 24.50 | 62.23 |
| | | | Seat Height -high | 24.00 | 60.98 |
| | | | Seat Height -high | 28.25 | 71.76 |
| | | | Seat Height -high | 35.00 | 88.90 |

[†]-xx in the product number is the placeholder for vinyl choice [‡] Replacement Cylinders may be ordered separately at additional cost [‡] Ships: Ground at no additional charge [‡] Submit COM sample for approval [‡] Adjust Height Specs: -1" (-2.54 cm)



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

