

SEATING

970 SERIES PHYSICIAN EXAM STOOL

970 Features:

- D-ring adjustment release
- Pneumatic height adjustment
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 **CF Stinson Pre-Approved Fabrics**

 **970 Options Detail and Photos**

Product Info:

 **971**  **972**  **973**




970 Options:

base with standard seat height	item no.†	procedure arm†	with footring†	with backrest†	crescent back†
Black Nylon Composite	971-xx	971PA-xx	971-FR-xx	971B-xx	971CB-xx
Brushed Aluminum² <i>and Toe Caps</i>	972-xx	972PA-xx	972-FR-xx	972B-xx	972CB-xx
Polished Aluminum	973-xx	973PA-xx	973-FR-xx	973B-xx	973CB-xx

options available at additional cost

970 Specifications:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter <i>with casters</i>	24.00	60.96	COM Requirement¹	1	0.92
Overall Height <i>with backrest -low</i>	33.50	85.09	Height Adjustment <i>method</i>	pneumatic			lb	kg
Overall Height <i>with backrest -high</i>	38.50	97.79	Height Release <i>hardware</i>	D-ring		Weight Capacity EVENLY DISTRIBUTED	350	158.76
						Shipping Weight	22	9.07

Adjustable Height Seat^{††}


		inch	cm		inch	cm	
Standard	97x-xx	Seat Height -low	18.50	46.99	Seat Height -high	24.00	60.96
7 [seven]	97x-7-xx	Seat Height -low	20.50	52.07	Seat Height -high	28.25	71.76
10 [ten]	97x-10-xx	Seat Height -low	24.50	62.23	Seat Height -high	35.00	88.90

†-xx in the product number is the placeholder for vinyl choice † Replacement Cylinders may be ordered separately at additional cost †
 † Ships: Ground at no additional charge † Submit COM sample for approval † Adjust Height Specs: -1" (-2.54 cm)



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

