

SEATING

980 SERIES PHYSICIAN EXAM STOOL

980 Features:

- 360° foot adjustment release
- Pneumatic height adjustment
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 **CF Stinson Pre-Approved Fabrics**

 **980 Options Detail and Photos**

Product Info:

 **981**  **982**  **983**


980 Options:

base with standard seat height	item number†	with backrest‡	crested back‡
Black Nylon Composite	981-xx	981B-xx	981CB-xx
Brushed Aluminum² <i>and Toe Caps</i>	982-xx	982B-xx	982CB-xx
Polished Aluminum	983-xx	983B-xx	983CB-xx

options available at additional cost



980 Specifications:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter <i>with casters</i>	24.00	60.96	COM Requirement¹	1	0.92
Overall Height <i>with backrest -low</i>	33.75	85.73	Height Adjustment <i>method</i>		pneumatic		lb	kg
Overall Height <i>with backrest -high</i>	38.75	98.43	Height Release <i>hardware</i>		360° foot	Weight Capacity EVENLY DISTRIBUTED	350	158.76
						Shipping Weight	25	11.34

Adjustable Height Seat^{††}


		inch	cm		inch	cm	
Standard	98x-xx	Seat Height -low	19.75	50.17	Seat Height -high	24.75	62.87
7 [seven]	98x-7-xx	Seat Height -low	23.50	59.69	Seat Height -high	30.00	76.20
10 [ten]	98x-10-xx	Seat Height -low	25.75	65.41	Seat Height -high	38.75	98.425

†-xx in the product number is the placeholder for vinyl choice †† Replacement Cylinders may be ordered separately at additional cost
¹ Submit COM sample for approval ² Adjust Height Specs: Overall: -1.25" (-3.175 cm); Standard Seat: -0.75" (-1.905 cm); 7 & 10 Seat -1" (-2.54 cm)
³ Ships: Ground at no additional charge



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Manufactured in the USA

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

