

SEATING

990 SERIES PHYSICIAN EXAM STOOL

990 Features:

- 360° Hand ring adjustment release
- Pneumatic height adjustment
- High density fire resistant foam
- Optional 18" (45.72 cm) footing
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 **CF Stinson Pre-Approved Fabrics**

 **990 Options Detail and Photos**

Product Info:


   **991 992 993**

990 Options:

base with standard seat height	item no. †	with footing ‡
Black Nylon Composite	991-xx	991-FR-xx
Brushed Aluminum² <i>and Toe Caps</i>	992-xx	992-FR-xx
Polished Aluminum	993-xx	993-FR-xx

options available at additional cost

990 Specifications:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter <i>with casters</i>	24.00	60.96	COM Requirement¹	1	0.92
Overall Height <i>-low</i>	18.00	45.72	Height Adjustment <i>method</i>		pneumatic		lb	kg
Overall Height <i>-high</i>	23.00	58.42	Height Release <i>hardware</i>		360° hand	Weight Capacity EVENLY DISTRIBUTED	350	158.76
Optional Footring	18.00	45.72				Shipping Weight	23	11.34

Adjustable Height Seat^{††}


		inch	cm		inch	cm	
Standard	99x-xx	Seat Height -low	18.00	45.72	Seat Height -high	23.00	58.42
7 [seven]	99x-7-xx	Seat Height -low	19.00	48.26	Seat Height -high	26.00	66.04
10 [ten]	99x-10-xx	Seat Height -low	22.00	55.88	Seat Height -high	32.00	81.28

† -xx in the product number is the placeholder for vinyl choice † Replacement Cylinders may be ordered separately at additional cost ‡ Ships: Ground at no additional charge † Submit COM sample for approval ‡ Adjust Height Specs -1" (-2.54 cm)



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.



Manufactured in the USA

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

