

SEATING

IN566DG SERIES

DELUXE DROP ARM RECLINER

The Deluxe Clinical Care recliner is ideal for clinical and acute care settings. It provides the best value while being engineered for maximum patient comfort, durability, and ease of use.

IN566DG Features:

- Three positions for patient comfort: Upright, First Recline, and Second Recline
- Fourth position: Trendelenburg (or Shock) position, Caregiver-activated using side pedal
- Tente® Swivel casters for ease of maneuverability. 3 Total Locking and 1 Directional Locking caster
- Wall saver armrest caps increase durability and protect walls from marring; also available with fully upholstered arms
- Removable upholstered side panels for easy cleaning and maintenance
- Two easy to use self-storing drop arms allow for a safer lateral transfer
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Rear push handle for improved patient transport
- Deluxe cushioning on seat and backrest for comfort

 CF Stinson Pre-Approved Fabrics

 IN566DG Product Info


IN566DG Options:

	item no.*
Standard Arm wall saver armrest tray	IN566DG-xx
Standard Arm wall saver armrest tray ² with heat and massage	IN566DGH-xx
Full Upholstered Arm	IN566DGxxPA
Full Upholstered Arm ² with heat and massage	IN566DGHxxPA

options available at additional cost



IN566DG Specifications:

	inch	cm		inch	cm		yd	m
Overall Width	32.50	82.55	Seat Width	20.50	52.07	COM Requirement ¹	6	5.49
Overall Depth	38.00	96.52	Seat Depth	19.00	48.26	 + FREIGHT *		
Overall Depth fully reclined	71.50	181.61	Seat Height	22.50	57.15		lb	kg
Overall Height	50.50	128.27	Back Height	30.00	76.20	Weight Capacity EVENLY DISTRIBUTED	350	158.76
Tente® Casters	5.00	12.70	Drop Arm Height from floor, arm up	28.00	71.12	Shipping Weight ²	142	64.41
			Drop Arm Height from floor, arm down	21.50	54.61			

*Ships: Freight request additional shipping quote from Intensa + -xx in the product number is the placeholder for vinyl choice

¹ Submit COM sample for approval ² Adjust Shipping Weight Specs: +2 lb (+ 0.91 kg)

IN566DG Accessories:


	item no.*
Headrest Pillow	FR5659xxUS
Headrest Skirt	HRSxx
Headrest Cover	HRC587xx
IV Pole	2101A
IV Mount	5782RG

accessories available at additional cost



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

