

## SEATING

# IN566G SERIES DELUXE CLINICAL RECLINER

The Deluxe Clinical Care recliner is ideal for clinical and acute care settings. It provides the best value while being engineered for maximum patient comfort, durability, and ease of use.

### IN566G Features:

- Three positions for patient comfort: Upright, First Recline, and Second Recline
- Fourth position: Trendelenburg (or Shock) position, Caregiver-activated using side pedal
- Tente® Swivel casters for ease of maneuverability. 3 Total Locking and 1 Directional Locking caster
- The wide, contoured seat and back, along with an integrated footrest, provide an enhanced level of patient comfort
- Wall saver armrest caps increase durability of the armrests and protect walls from marring
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Smooth and removable side panels for easy cleaning and maintenance
- Rear push handle for improved patient transport

 CF Stinson Pre-Approved Fabrics

 IN566G Product Info


### IN566G Options:

	item no.*
Deluxe Clinical Care	IN566G-xx
Deluxe Clinical Care <sup>2</sup> <i>with heat and massage</i>	IN566GH-xx

options available at additional cost



### IN566G Specifications:

	inch	cm		inch	cm	yd	m
Overall Width	31.25	79.38	Seat Width	20.50	52.07	COM Requirement <sup>1</sup>	6 5.49
Overall Depth	37.50	95.25	Seat Depth	19.00	48.26	 <b>+ FREIGHT</b> *	lb kg
Overall Depth <i>fully reclined</i>	71.50	181.61	Seat Height	22.50	57.15		
Overall Height	49.50	125.73	Back Height	30.00	76.20	Shipping Weight	139 63.05
Tente® Casters	5.00	12.70	Arm Height <i>from seat</i>	8.00	20.32		
			Arm Width <i>width between arms</i>	21.75	55.25		

\*Ships: Freight request additional shipping quote from Intensa \*-xx in the product number is the placeholder for vinyl choice  
<sup>1</sup> Submit COM sample for approval <sup>2</sup> Adjust Shipping Weight Specs: +2 lb (+ 0.91 kg) <sup>3</sup> adds 4 inch (10.16 cm) height

### IN566G Accessories:


	item no.*
Headrest Pillow	INFR5659xxUS
Headrest Skirt	INHRSxx
Headrest Cover	INHRC566xx
Extension Legs <sup>3</sup>	IN5699G
IV Pole	IN2101
IV Mount	IN5782RG

accessories available at additional cost



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

Manufactured in the USA

