

SEATING

IN577RG SERIES CLINICAL CARE RECLINER

The most versatile Intensa Clinical Care recliner is ideal for treatment within extended care, clinical and acute care settings. It is an exceptional chair that provides comfort, utility, and flexibility for multiple treatment applications.

IN577RG Features:

- Three positions for patient comfort: Upright, First Recline, and Second Recline can be activated by patient or Caregiver
- Fourth Trendelenburg (or Shock) position, Caregiver-activated using rear mounted hand trigger
- Tente® Swivel Casters for ease of maneuverability. 3 Total Locking and 1 Directional Locking caster
- Gas cylinder give smooth power assist during Trendelenburg activation
- Wall saver armrest caps increase durability and protect walls from marring
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Removable side panels for cleaning and maintenance
- Back-mounted arm pads for additional support
- Rear push handle for improved patient transport

 CF Stinson Pre-Approved Fabrics

 IN577RG Product Info


IN577RG Options:

	Item no.*
Arm Pad Recliner	IN577RG-xx
Arm Pad Recliner ² <i>with heat and massage</i>	IN577RGH-xx

options available at additional cost



IN577RG Specifications:

	inch	cm		inch	cm		yd	m
Overall Width	31.50	80.01	Seat Width	21.50	54.61	COM Requirement ¹	6	5.49
Overall Depth	39.00	99.06	Seat Depth	19.50	49.53	 + FREIGHT*		
Overall Depth <i>fully reclined</i>	65.00	165.10	Seat Height	20.50	52.07		lb	kg
Overall Height	48.00	121.92	Arm Height <i>from floor</i>	27.25	69.22	Weight Capacity EVENLY DISTRIBUTED	290	131.54
Tente® Casters	5.00	12.70	Back Height	29.00	73.66	Shipping Weight ²	129	58.51

*Ships: Freight request additional shipping quote from Intensa + -xx in the product number is the placeholder for vinyl choice
¹ Submit COM sample for approval ² Adjust Shipping Weight Specs: +2 lb (+ 0.91 kg) ³ adds 4 inch (10.16 cm) height

IN577RG Accessories:


	Item no.*
Headrest Pillow	FR5659xxUS
Headrest Skirt	HRSxx
Headrest Cover	HRC577xx
4" Extension Legs ³	5699G
IV Pole	2101A
IV Mount	5782RG

accessories available at additional cost



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

