

## SEATING

### IN587W SERIES

# WIDE CLINICAL CARE RECLINER

With 3 positions and trendelenburg, the Extra Wide Clinical Care Recliner is ideal for clinical and acute care settings. This recliner provides patients more maneuverability space while providing excellent value that maximizes patient comfort, product durability and ease of use.

#### IN587W Features:

- Three positions for patient comfort: Upright, First Recline, and Second Recline
- Fourth position: Trendelenburg (or Shock) position, Caregiver-activated using rear mounted hand trigger
- Tente® Swivel casters for ease of maneuverability. 3 Total Locking and 1 Directional Locking caster
- Dual gas cylinders give smooth power assist during Trendelenburg activation
- Extra-wide seat and back provides an enhanced level of patient comfort
- Two easy to use self-storing drop arms allow for a safer lateral transfer
- Wall saver armrest caps increase durability and protect walls from marring
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Removable side panels for easy cleaning
- Rear push handle for improved patient transport

 **CF Stinson** Pre-Approved Fabrics

 **IN587W** Product Info


#### IN587W Options:

	item no.†
Standard Arms	IN587W-xx
Standard Arms <sup>2</sup> with heat and massage	IN587WH-xx
Upholstered Arms	IN587WPA-xx
Upholstered Arms <sup>2</sup> with heat and massage	IN587WPAH-xx

options available at additional cost



#### IN587W Specifications:

	inch	cm		inch	cm		yd	m
Overall Width	33.00	83.82	Seat Width	23.50	59.69	COM Requirement <sup>1</sup>	6	5.49
Overall Depth	44.00	111.76	Seat Depth	20.50	52.07	 <b>+ FREIGHT</b> *		
Overall Depth fully reclined	71.00	180.34	Seat Height	22.00	55.88		lb	kg
Overall Height	49.50	125.73	Arm Height from floor	28.00	71.12	Weight Capacity EVENLY DISTRIBUTED	450	204.12
Tente® Casters	5.00	12.70	Back Height	31.50	76.20	Shipping Weight	145	65.77

\*Ships: Freight request additional shipping quote from Intensa <sup>†</sup>-xx in the product number is the placeholder for vinyl choice  
<sup>1</sup> Submit COM sample for approval <sup>2</sup> Adjust Shipping Weight Specs: +2 lb (+ 0.91 kg) <sup>3</sup> adds 4 inch (10.16 cm) height

#### IN587W Accessories:


	item no.†
Headrest Pillow	FR5659xxUS
Headrest Skirt	HRSxx
Headrest Cover	HRC587xx
4" Extension Legs <sup>3</sup>	FR5699G
IV Pole	2101A
IV Mount	5782RG

accessories available at additional cost



GF Health Products, Inc. is an  
ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

