

## SEATING

# IN587WD SERIES WIDE DROP ARM RECLINER

### Clinical Care Recliner Wide with 3 Positions and Trendelenburg

The Extra Wide Clinical Care Recliner is ideal for clinical and acute care settings. This recliner provides more space for patients to adjust their seated position while providing excellent value that maximizes patient comfort, product durability and ease of use.

### IN587WD Features:

- Three positions for patient comfort: Upright, First Recline, and Second Recline
- Fourth position: Trendelenburg (or Shock) position, Caregiver-activated using rear mounted hand trigger
- Tente® Swivel casters for ease of maneuverability. 3 Total Locking and 1 Directional Locking caster
- Dual gas cylinders give smooth power assist during Trendelenburg activation
- Extra-wide seat and back provides an enhanced level of patient comfort
- Wall saver armrest caps increase durability and protect walls from marring
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Removable upholstered side panels for easy cleaning and maintenance
- Two easy to use self-storing drop arms allow for a safer lateral transfer
- Rear push handle for improved patient transport


 **CF Stinson** Pre-Approved Fabrics

 **IN587WD** Product Info


### IN587WD Options:

|  | Item no.*    |
|--|--------------|
| Standard Arm <i>wall saver armrest tray</i>  | IN587WD-xx   |
| Standard Arm <i>wall saver armrest tray</i> <sup>2</sup><br><i>with heat and massage</i> | IN587WDH-xx  |
| Full Upholstered Arm   | IN587WDxxPA  |
| Full Upholstered Arm <sup>2</sup><br><i>with heat and massage</i>                        | IN587WDHxxPA |

options available at additional cost



### IN587WD Specifications:

|  | inch  | cm     |   | inch  | cm    |  | yd                                    | m     |
|--|-------|--------|---|-------|-------|--|---------------------------------------|-------|
| Overall Width                          | 33.00 | 83.82  | Seat Width                              | 23.50 | 59.69 | COM Requirement <sup>1</sup>   | 6                                     | 5.49  |
| Overall Depth                          | 44.00 | 111.76 | Seat Depth                              | 20.50 | 52.07 |  <b>+ FREIGHT</b> * | lb                                    | kg    |
| Overall Depth<br><i>fully reclined</i> | 71.00 | 180.34 | Seat Height                             | 22.00 | 55.88 |  | Weight Capacity<br>EVENLY DISTRIBUTED | 450   |
| Overall Height                         | 49.50 | 125.73 | Back Height                             | 31.50 | 76.20 | Shipping Weight  | 145                                   | 66.68 |
| Tente® Casters                         | 5.00  | 12.70  | Drop Arm Height<br>from floor, arm up   | 28.00 | 71.12 |  |                                       |       |
|  |       |        | Drop Arm Height<br>from floor, arm down | 21.50 | 54.61 |  |                                       |       |

\*Ships: Freight request additional shipping quote from Intensa <sup>1</sup>-xx in the product number is the placeholder for vinyl choice  
<sup>1</sup> Submit COM sample for approval <sup>2</sup> Adjust Shipping Weight Specs: +2 lb (+ 0.91 kg)


### IN587WD Accessories:

|                 | Item no.*  |
|-----------------|------------|
| Headrest Pillow | FR5659xxUS |
| Headrest Skirt  | HRSxx      |
| Headrest Cover  | HRC587xx   |
| IV Pole         | 2101A      |
| IV Mount        | 5782RG     |

accessories available at additional cost

 GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

Manufactured in the USA

