

SEATING

IN597 SERIES

ORTHO-BIOTIC II RECLINER

with Infinite Back Positioning and Trendelenburg The Premium recliner that is the most aesthetically pleasing and provides the ultimate comfort for a variety of applications in extended care, clinical and acute care settings.

IN597 Features:

- Infinite recline positions for patient comfort
- Trendelenburg (or Shock) position, Caregiver-activated using side pedal
- Tente® Swivel casters for ease of maneuverability. 3 Total Locking and 1 Directional Locking caster
- Visco (Memory) foam engineered into the seat for ultimate comfort
- Legrest-Recline lever activates independent movement from backrest
- Self-storing patient footplate for transport
- Rear push handle for improved patient transport
- Headrest pillow included for extra comfort

 **CF Stinson** Pre-Approved Fabrics

 **IN597** Product Info


IN597 Options:

	Item no.*
Wood Arms	IN597G-xx
Wood Arms ² with heat and massage	IN597GH-xx
Upholstered Arms	IN597P-xx
Upholstered Arms ² with heat and massage	IN597PH-xx

options available at additional cost



IN597 Specifications:

	inch	cm		inch	cm		yd	m
Overall Width	33.00	83.82	Seat Width	23.00	58.42	COM Requirement ¹	6	5.49
Overall Depth	36.50	92.71	Seat Depth	20.50	52.07	 + FREIGHT *		
Overall Depth fully reclined	64.00	162.56	Seat Height	21.00	53.34		lb	kg
Overall Height	45.50	115.57	Arm Height from floor	26.75	67.95	Weight Capacity EVENLY DISTRIBUTED	350	158.76
Tente® Casters	4.00	10.16	Back Height	31.00	78.74	Shipping Weight	210	95.25

*Ships: Freight request additional shipping quote from Intensa

¹ Submit COM sample for approval ² Adjust Shipping Weight Specs: +2 lb (+ 0.91 kg)

*-xx in the product number is the placeholder for vinyl choice

IN597 Accessories:


	Item no.*
Headrest Pillow	standard
Headrest Skirt	HRSxx
Side Table	FR597GTBLMNT
IV Pole	2101A
IV Mount	FR597GIVMNT

accessories available at additional cost



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

