

SEATING

IN601 SERIES PIVOT ARM RECLINER

Clinical Care Recliner with Pivot Arms

A premium recliner that focuses on being the most aesthetically pleasing, comfortable, and durable recliner for use across a large variety of clinical and acute care settings.

IN601 Features:

- Three positions: Upright, First Recline, and Second Recline
- Fourth Trendelenburg (or Shock) position, Caregiver-activated using side pedal
- Tente® Swivel casters for ease of maneuverability. 3 Total Locking and 1 Directional Locking caster
- Visco (Memory) foam engineered into the seat for ultimate comfort
- Wide contoured seat and back provides an enhanced level of patient comfort
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Two pivot arms for easy convenient cleaning
- Urethane capped or fully upholstered arms
- Self-storing patient footplate for safer patient transport
- Rear push handle for improved patient transport

 **CF Stinson Pre-Approved Fabrics**

 **IN601 Product Info**

IN601 Options:

Urethane Capped Arms

| Item no.† |
|-----------|
| IN601U-xx |

Urethane Capped Arms
with heat and massage

| |
|------------|
| IN601UH-xx |
|------------|

Fully Upholstered Arms

| |
|-----------|
| IN601P-xx |
|-----------|

Fully Upholstered Arms²
with heat and massage

| |
|------------|
| IN601PH-xx |
|------------|

**comprehensive
Service Parts Program**


options available at additional cost



IN601P-xx



IN601 Specifications:

| | inch | cm | | inch | cm | | yd | m |
|--|-------|--------|---------------------------------|-------|-------|---|-----|--------|
| Overall Width | 34.50 | 87.63 | Seat Width | 22.50 | 57.15 | COM Requirement ¹ | 6 | 5.49 |
| Overall Depth | 38.00 | 96.52 | Seat Depth | 19.00 | 48.26 |  + FREIGHT [*] | | |
| Overall Depth fully reclined | 68.00 | 172.72 | Seat Height | 22.00 | 55.88 | | lb | kg |
| Overall Height | 51.50 | 130.81 | Arm Height from floor | 27.00 | 68.58 | Weight Capacity EVENLY DISTRIBUTED | 400 | 181.44 |
| Tente® Casters | 5.00 | 12.70 | Back Height | 30.00 | 76.20 | Shipping Weight | 197 | 89.36 |

^{*}Ships: Freight request additional shipping quote from Intensa

[†]-xx in the product number is the placeholder for vinyl choice

¹Submit COM sample for approval

²Adjust Shipping Weight Specs: +2 lb (+ 0.91 kg)

IN601 Accessories:


| | Item no.† |
|------------------------|-------------|
| Headrest Pillow | FR5659xxUS |
| Headrest Skirt | HRSxx |
| IV Pole | 2101A |
| IV Mount | FR597GIVMNT |

accessories available at additional cost



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

