



HSS SERIES SURGICAL STOOL



HSS500-S SADDLE SEAT STYLE



HSS500-WD WATER DROP SEAT STYLE



HSS500-W WEDGE SEAT STYLE

BUILT FOR PERFORMANCE
AND COMFORT



HSS500 SURGICAL STOOL

IDEAL FOR EYE, ENT AND PLASTIC PROCEDURES




360° Rotating Urethane Arm Rest for Comfortable Arm Positioning

Adjustable Contoured Backrest

300 POUNDS
300 lb (136 kg)
Weight Capacity
EVENLY DISTRIBUTED

Dual-density Foam Seat with **DURA-GEL®** Insert for Enhanced Comfort



180° Arm Rest Support for more Stable Arm Positioning

Adjustable Height Range
21 - 28 inches (53 - 71 cm)

OPTIONAL ACCESSORY:
Riser Kit Increases Height Range an additional 2 inches (5 cm)
Item No: H081025

ABS Thermoformed Hood for Easy Cleaning

FOOT-CONTROLLED HYDRAULICS with Three-Spindle Base and Four 3 inch (8 cm) Dual Tente® Casters for Greater Stability



Exit
Loading
System

ENGINEERED TO PROVIDE
ENHANCED
SURGEON COMFORT

Color variations may occur due to monitor or printer settings.

© 2021, GF Health Products, Inc. All Rights Reserved.
Intensa and Hausted are registered trademarks of GF Health Products, Inc.
GF Health Products, Inc. is an ISO 13485:2016 Certified Company.



CE Partner4U BV
Esdoornlaan 13
3951 DB Maarn
The Netherlands
www.cepartner4u.eu

Country of Origin: United States

Intensa, Inc. | P.O. Box 5981 | High Point, NC 27262 | USA
+1 888.636.2660 | intensa@grahamfield.com | intensa.net

GF2100052 RevB23

