



Product Specification & Pricing Catalog

published *November 2025*





Put a little **COLOR** in your exam room with Intensa!



Footnotes:

- * ships freight: request shipping quote from customer service
 - ◇ ships ground: no additional charge
 - 🔗 hyperlink - additional product information and ordering
 - 🔗 page link - jumps to referenced page
 - ‡ xx in the product number is the placeholder for vinyl / finish choice
 - † standard rounding method used
 - § representative image shows this option
 - ★ options available at additional cost & may alter weight specs - order requirements noted
1. maximum weight capacity - evenly distributed - refer to user manual
 2. procedure arm is not compatible / available with either backrest option
 3. variance: ± 0.5 inch (± 1.27 cm)
 4. supports a single, left side, table only
 5. requires factory installation include item with purchase
 6. pullout footplate not compatible / available with extended legrest option
 7. interior dimensions
 8. approximate measurement, options may alter spec
 9. product is not CFFA-HealthCare 201 Certified
 10. submit required sample size indicated to:
Intensa COM Department, 1810 S. Elm Street, High Point, NC 27260.
 11. adjust height spec: +4 inch (+10.16 cm)
 12. replacement cylinders may be ordered separately at additional cost [🔗](#)





Table of Contents



<i>Physician Exam</i> Seating	4
<i>Laboratory</i> Seating	10
<i>Hydraulic Surgical</i> Seating	19
<i>Clinical Care</i> Recliners	20
<i>Blood Draw / Phlebotomy</i> Seating	24
<i>Exam Room</i> Treatment Tables	25
<i>Replacement</i> Exam Table Tops	28
<i>Patient / Guest</i> Seating	30
<i>Overbed</i> Tables	31
<i>Step</i> Stools	32
<i>Exam</i> Lights	33
<i>IV</i> Poles	33
<i>Recliner</i> Options	34
<i>Stool</i> Options	35
<i>Overbed Table Finish</i> Swatches & Care	36
<i>Tier 1 Vinyl Upholstery</i> Standard Collection Swatches	37
<i>Stinson® Vinyl Upholstery</i> Pattern Collection Swatches	38
<i>Mayer Vinyl Upholstery</i> Texture Collection Swatches	39
<i>Approved Vinyl Upholstery</i> Tier, Care, Link & Swatch Chart	40
<i>Spradling® Vinyl Upholstery</i> Collections Swatches	41
<i>Vinyl Upholstery</i> Recommended Care, Cleaning & Disinfection	43
<i>Upholstery</i> Composition, Finish, & Performance Comparison	44
<i>Product</i> Index	45

click listing to jump to page

Purchase Orders

please remit purchase orders to intensa@grahamfield.com

Pricing & Availability

current at time of publication & subject to change without notice

Contact Customer Service

for most current information, pricing, freight quotes, and support +1.336.884.4096 | intensa@grahamfield.com

Color Variations

perceived texture and coloring may vary on finished product. monitor and/or printer settings will affect appearance of color

Intensa products are labeled in accordance with California Proposition 65 Warnings for products potentially sold to California consumers

 **WARNING:** Cancer and Reproductive Harm - www.P65Warnings.ca.gov

 **WARNING:** This product can expose you to chemicals including Formaldehyde (gas) which is known to the State of California to cause cancer. For more information go to - www.P65Warnings.ca.gov/furniture

 **WARNING:** This product can expose you to chemicals including Di2-ethylhexylphthalate (DEHP) which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to - www.P65Warnings.ca.gov/furniture

993 Physician Stool with Chrome 360° Hand Release and Polished Aluminum Base



	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$518	\$518	\$598	\$618	\$631	\$677	\$702	\$718	\$725	\$760	\$813

† Specs:		inch	cm		inch	cm		yd	m
	Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height -low	18.00	45.72	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
Overall Height -high	23.00	55.88	Height Release hardware	360° hand		Shipping Weight ⁸	23	11.34	
Optional Footring	18.00	45.72							

Adjustable Height Seat ^{††12} add			inch	cm	inch	cm		
5 [five] [§]	993-xx	\$0	Seat Height -low	18.00	45.72	Seat Height -high	23.00	58.42
7 [seven]	993-7-xx	\$38	Seat Height -low	19.00	48.26	Seat Height -high	26.00	66.04
10 [ten]	993-10-xx	\$48	Seat Height -low	22.00	55.88	Seat Height -high	32.00	81.28

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		993-FR-xx	\$110	not available	

993 Options:
[details & photos on page 35*](#)

[Product Info](#)

992 Physician Stool with Chrome 360° Hand Release and Brushed Aluminum Base with Toe Caps



	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$509	\$509	\$589	\$609	\$622	\$668	\$693	\$709	\$716	\$751	\$804

† Specs:		inch	cm		inch	cm		yd	m
	Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height -low	17.00	43.18	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
Overall Height -high	22.00	55.88	Height Release hardware	360° hand		Shipping Weight ⁸	23	11.34	
Optional Footring	18.00	45.72							

Adjustable Height Seat ^{††12} add			inch	cm	inch	cm		
5 [five] [§]	992-xx	\$0	Seat Height -low	17.00	43.18	Seat Height -high	22.00	55.88
7 [seven]	992-7-xx	\$38	Seat Height -low	18.00	45.72	Seat Height -high	25.00	63.50
10 [ten]	992-10-xx	\$48	Seat Height -low	21.00	53.34	Seat Height -high	31.00	78.74

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		992-FR-xx	\$110	not available	

992 Options:
[details & photos on page 35*](#)

[Product Info](#)

991 Physician Stool with Chrome 360° Hand Release and Black Nylon Composite Base



	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$425	\$425	\$505	\$525	\$538	\$584	\$609	\$625	\$632	\$667	\$720

† Specs:		inch	cm		inch	cm		yd	m
	Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height -low	18.00	45.72	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
Overall Height -high	23.00	55.88	Height Release hardware	360° hand		Shipping Weight ⁸	23	11.34	
Optional Footring	18.00	45.72							

Adjustable Height Seat ^{††12} add			inch	cm	inch	cm		
5 [five] [§]	991-xx	\$0	Seat Height -low	18.00	45.72	Seat Height -high	23.00	58.42
7 [seven]	991-7-xx	\$38	Seat Height -low	19.00	48.26	Seat Height -high	26.00	66.04
10 [ten]	991-10-xx	\$48	Seat Height -low	22.00	55.88	Seat Height -high	32.00	81.28

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		991-FR-xx	\$110	not available	

991 Options:
[details & photos on page 35*](#)

[Product Info](#)

983 Physician Stool with 360° Foot Release and Polished Aluminum Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$744	\$744	\$824	\$844	\$857	\$903	\$928	\$944	\$951	\$986	\$1,039

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	33.75	85.73	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	lb	kg
Overall Height with backrest -high	38.75	98.43	Height Release hardware	360° foot			350	158.76
						Shipping Weight ⁸	25	11.34

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm
5 [five] [§]	983-xx	\$0		Seat Height -low	19.75 50.17
7 [seven]	983-7-xx	\$38		Seat Height -low	23.50 59.69
10 [ten]	983-10-xx	\$48		Seat Height -low	25.75 65.41

983 Options:
[details & photos on page 35*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		not available		Full 983B-xx	\$175
								Crescent 983CB-xx	\$315

982 Physician Stool with 360° Foot Release and Brushed Aluminum Base with Toe Caps



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$718	\$718	\$798	\$818	\$831	\$877	\$902	\$918	\$925	\$960	\$1,013

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	32.50	82.55	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	lb	kg
Overall Height with backrest -high	37.50	95.25	Height Release hardware	360° foot			350	158.76
						Shipping Weight ⁸	25	11.34

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm
5 [five] [§]	982-xx	\$0		Seat Height -low	19.00 48.26
7 [seven]	982-7-xx	\$38		Seat Height -low	22.50 59.69
10 [ten]	982-10-xx	\$48		Seat Height -low	24.75 62.87

982 Options:
[details & photos on page 35*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		not available		Full 982B-xx	\$175
								Crescent 982CB-xx	\$315

981 Physician Stool with 360° Foot Release and Black Nylon Composite Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$649	\$649	\$729	\$749	\$762	\$808	\$833	\$849	\$856	\$891	\$944

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	33.75	85.73	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	lb	kg
Overall Height with backrest -high	38.75	98.43	Height Release hardware	360° foot			350	158.76
						Shipping Weight ⁸	25	11.34

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm
5 [five] [§]	981-xx	\$0		Seat Height -low	19.75 50.17
7 [seven]	981-7-xx	\$38		Seat Height -low	23.50 59.69
10 [ten]	981-10-xx	\$48		Seat Height -low	25.75 65.41

981 Options:
[details & photos on page 35*](#)

[Product Info](#)

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		not available		Full 981B-xx	\$175
								Crescent 981CB-xx	\$315

973 Physician Stool with D-Handle Release and Polished Aluminum Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$494	\$494	\$574	\$594	\$607	\$653	\$678	\$694	\$701	\$736	\$789

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	33.50	85.09	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	38.50	97.79	Height Release hardware	D-ring		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Optional Footring	18.00	45.72				Shipping Weight⁸	22	9.07

Options:

[details & photos on page 35*](#)

Adjustable Height Seat^{††12}	add	inch	cm		inch	cm		inch	cm
5 [five] \$	973-xx	\$0		Seat Height -low	18.50	46.99	Seat Height -high	24.00	60.96
7 [seven]	973-7-xx	\$38		Seat Height -low	20.50	52.07	Seat Height -high	28.25	71.76
10 [ten]	973-10-xx	\$48		Seat Height -low	24.50	62.23	Seat Height -high	35.00	88.90

Casters	add	T-Style Armrest[†]	add	Procedure Arm^{†2}	add	Footring[†]	add	Backrest^{†2}	add
page 35	varies	not available		973PA-xx	\$450	973-FR-xx	\$110	Full	973B-xx \$175
								Crescent	973CB-xx \$315

[Product Info](#)

972 Physician Stool with D-Handle Release and Brushed Aluminum Base with Toe Caps



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$486	\$486	\$566	\$586	\$599	\$645	\$670	\$686	\$693	\$728	\$781

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	32.50	82.55	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	37.50	95.25	Height Release hardware	D-ring		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Optional Footring	18.00	45.72				Shipping Weight⁸	22	9.07

Options:

[details & photos on page 35*](#)

Adjustable Height Seat^{††12}	add	inch	cm		inch	cm		inch	cm
5 [five] \$	972-xx	\$0		Seat Height -low	17.50	44.45	Seat Height -high	23.00	58.42
7 [seven]	972-7-xx	\$38		Seat Height -low	19.50	49.53	Seat Height -high	27.25	69.22
10 [ten]	972-10-xx	\$48		Seat Height -low	23.50	59.69	Seat Height -high	34.00	86.36

Casters	add	T-Style Armrest[†]	add	Procedure Arm^{†2}	add	Footring[†]	add	Backrest^{†2}	add
page 35	varies	not available		972PA-xx	\$450	972-FR-xx	\$110	Full	972B-xx \$175
								Crescent	972CB-xx \$315

[Product Info](#)

971 Physician Stool with D-Handle Release and Black Nylon Composite Base



Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$415	\$415	\$495	\$515	\$528	\$574	\$599	\$615	\$622	\$657	\$710

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement¹⁰	1	0.92
Overall Height with backrest -low	33.50	85.09	Height Adjustment method	pneumatic			lb	kg
Overall Height with backrest -high	38.50	97.79	Height Release hardware	D-ring		Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76
Optional Footring	18.00	45.72				Shipping Weight⁸	22	9.07

Options:

[details & photos on page 35*](#)

Adjustable Height Seat^{††12}	add	inch	cm		inch	cm		inch	cm
5 [five] \$	971-xx	\$0		Seat Height -low	18.50	46.99	Seat Height -high	24.00	60.96
7 [seven]	971-7-xx	\$38		Seat Height -low	20.50	52.07	Seat Height -high	28.25	71.76
10 [ten]	971-10-xx	\$48		Seat Height -low	24.50	62.23	Seat Height -high	35.00	88.90

Casters	add	T-Style Armrest[†]	add	Procedure Arm^{†2}	add	Footring[†]	add	Backrest^{†2}	add
page 35	varies	not available		971PA-xx	\$450	971-FR-xx	\$110	Full	971B-xx \$175
								Crescent	971CB-xx \$315

[Product Info](#)

963 Physician Stool with Single Lever Release and Polished Aluminum Base



963 Options: details & photos on page 35*

Product Info

Casters	add
page 35	varies

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$491	\$491	\$571	\$591	\$604	\$650	\$675	\$691	\$698	\$733	\$786

Specs:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	33.50	85.09	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Overall Height with backrest -high	38.50	97.79	Height Release hardware	single lever		Shipping Weight ⁸	21	53.34
Optional Footring	18.00	45.72						

Adjustable Height Seat^{††12} add

	inch	cm		inch	cm
5 [five] \$	96x-xx	\$0	Seat Height -low	18.50	46.99
7 [seven]	96x-7-xx	\$38	Seat Height -low	20.50	52.07
10 [ten]	96x-10-xx	\$48	Seat Height -low	24.50	62.23

	inch	cm		inch	cm
Seat Height -high	24.00	60.98	Seat Height -high	28.25	71.76
Seat Height -high	28.25	71.76	Seat Height -high	35.00	88.90

T-Style Armrest [†]	add	Procedure Arm ^{†2}	add	Footring [†]	add	Backrest ^{†2}	add
not available		963PA-xx	\$450	963-FR-xx	\$110	Full 963B-xx	\$175
						Crescent 963CB-xx	\$315

962 Physician Stool with Single Lever Release and Brushed Aluminum Base with Toe Caps



962 Options: details & photos on page 35*

Product Info

Casters	add
page 35	varies

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$479	\$479	\$559	\$579	\$592	\$638	\$663	\$679	\$686	\$721	\$774

Specs:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	32.50	82.55	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Overall Height with backrest -high	37.50	95.25	Height Release hardware	single lever		Shipping Weight ⁸	21	53.34
Optional Footring	18.00	45.72						

Adjustable Height Seat^{††12} add

	inch	cm		inch	cm
5 [five] \$	96x-xx	\$0	Seat Height -low	17.50	44.45
7 [seven]	96x-7-xx	\$38	Seat Height -low	19.50	49.53
10 [ten]	96x-10-xx	\$48	Seat Height -low	23.50	59.69

	inch	cm		inch	cm
Seat Height -high	23.00	58.42	Seat Height -high	28.25	71.76
Seat Height -high	28.25	71.76	Seat Height -high	34.00	86.36

T-Style Armrest [†]	add	Procedure Arm ^{†2}	add	Footring [†]	add	Backrest ^{†2}	add
not available		962PA-xx	\$450	962-FR-xx	\$110	Full 962B-xx	\$175
						Crescent 962CB-xx	\$315

961 Physician Stool with Single Lever Release and Black Nylon Composite Base



961 Options: details & photos on page 35*

Product Info

Casters	add
page 35	varies

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$396	\$396	\$476	\$496	\$509	\$555	\$580	\$596	\$603	\$638	\$691

Specs:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	33.50	85.09	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Overall Height with backrest -high	39.00	99.06	Height Release hardware	single lever		Shipping Weight ⁸	20	9.07
Optional Footring	18.00	45.72						

Adjustable Height Seat^{††12} add

	inch	cm		inch	cm
5 [five] \$	96x-xx	\$0	Seat Height -low	18.50	46.99
7 [seven]	96x-7-xx	\$38	Seat Height -low	20.50	52.07
10 [ten]	96x-10-xx	\$48	Seat Height -low	24.50	62.23

	inch	cm		inch	cm
Seat Height -high	24.00	60.98	Seat Height -high	28.25	71.76
Seat Height -high	28.25	71.76	Seat Height -high	35.00	88.90

T-Style Armrest [†]	add	Procedure Arm ^{†2}	add	Footring [†]	add	Backrest ^{†2}	add
not available		961PA-xx	\$450	961-FR-xx	\$110	Full 961B-xx	\$175
						Crescent 961CB-xx	\$315

958 Industrial Stool with Single Lever Release, Footring, 5 Leg Black Tubular Steel Base, and Glides



Product Info

958
Options:
details & photos on page 35*

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$610	\$610	\$690	\$710	\$723	\$769	\$794	\$810	\$817	\$852	\$905
†		inch	cm	inch		cm	yd		m	
Specs:	Seat Diameter	16.00	40.64	Base Diameter <i>5 leg</i>	24.00	60.96	COM Requirement ¹⁰		1	0.92
	Overall Height <i>with backrest -low</i>	33.00	83.82	Height Adjustment <i>method</i>	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
	Overall Height <i>with backrest -high</i>	39.00	99.06	Height Release <i>hardware</i>	single lever		Shipping Weight ⁸		23	10.43
	Footring	20.00	50.80	Adjustable Height Seat ^{††12}		add	inch		cm	
	5 [five] §	958-xx	\$0	Seat Height -low	17.50	44.45	Seat Height -high		23.00	58.42
	7 [seven]	958-7-xx	\$38	Seat Height -low	19.50	49.53	Seat Height -high		26.50	67.31
	10 [ten]	958-10-xx	\$48	Seat Height -low	23.50	59.69	Seat Height -high		34.00	86.36
T-Style Armrest [†]		add	Procedure Arm ^{†2}	add	Footring [†]	add	Backrest ^{†2}		add	
not available			958PA-xx	\$450	standard	\$0	Full 958B-xx		\$175	
							Crescent 958CB-xx		\$315	

951 Physician Stool with Single Lever Release and Black Nylon Composite Base



Product Info

951
Options:
details & photos on page 35*

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$395	\$395	\$475	\$495	\$508	\$554	\$579	\$595	\$602	\$637	\$690
†		inch	cm	inch		cm	yd		m	
Specs:	Seat Diameter	16.00	40.64	Base Diameter <i>with casters</i>	24.00	60.96	COM Requirement ¹⁰		1	0.92
	Overall Height <i>with backrest -low</i>	33.00	83.82	Height Adjustment <i>method</i>	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
	Overall Height <i>with backrest -high</i>	39.00	99.06	Height Release <i>hardware</i>	single lever		Shipping Weight ⁸		20	9.07
	Optional Footring	18.00	45.72	Adjustable Height Seat ^{††12}		add	inch		cm	
	5 [five] §	958-xx	\$0	Seat Height -low	18.50	46.99	Seat Height -high		24.00	60.98
	7 [seven]	958-7-xx	\$38	Seat Height -low	20.50	52.07	Seat Height -high		28.25	71.76
	10 [ten]	958-10-xx	\$48	Seat Height -low	24.50	62.23	Seat Height -high		35.00	88.90
Casters		add	T-Style Armrest [†]	add	Procedure Arm ^{†2}	add	Footring [†]	add	Backrest ^{†2}	
page 35		varies	not available		951PA-xx	\$450	951-FR-xx	\$110	Full 951B-xx	
									Crescent 951CB-xx	
									\$315	

940 SERIES Physician Exam Sit / Stand Seating

941 Urethane Sit / Stand Stool with Single Lever Release and Black Nylon Composite Base



Product Info

941
Options:
details & photos on page 35*

list										
\$377	available in black urethane only									
†		inch	cm	inch		cm	yd		m	
Specs:	Seat Width	12.00	30.48	Base Diameter <i>with glides</i>	24.00	60.96	COM Requirement ¹⁰		1	0.92
	Seat Depth	12.00	30.48	Optional Footring	20.00	50.80	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
	Overall Height <i>with backrest -low</i>	21.00	53.34	Height Adjustment <i>method</i>	pneumatic		Shipping Weight ⁸		21	9.53
	Overall Height <i>with backrest -high</i>	31.50	80.01	Height Release <i>hardware</i>	single lever					
	Adjustable Height Seat ^{††12}		add	inch		cm	inch		cm	
	5 [five]	941-5-xx	\$0	Seat Height -low	15.50	39.37	Seat Height -high		20.75	52.71
	7 [seven]	941-7-xx	\$0	Seat Height -low	17.25	43.82	Seat Height -high		24.75	62.87
	10 [ten] §	941-xx standard	\$0	Seat Height -low	21.00	53.34	Seat Height -high		31.50	80.01
Casters		add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	
page 35		varies	not available		not available		941-FR-xx	\$110	not available	

933 Physician Stool with Manual Spin Lift Adjustment and Polished Aluminum Base



[Product Info](#)

[details & photos on page 35*](#)

Options:
[page 35](#) varies

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$476	\$476	\$556	\$576	\$589	\$635	\$660	\$676	\$683	\$718	\$771
Specs:	†		inch	cm	inch		cm	yd		m	
Seat Diameter	16.00		40.64	Base Diameter with casters		24.00	60.96	COM Requirement ¹⁰		1	0.92
Overall Height with backrest -low	33.00		83.82	Seat Height -low		19.50	49.53	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
Overall Height with backrest -high	39.00		99.06	Seat Height -high		23.00	58.42	Shipping Weight ⁸		20	9.07
Height Release hardware	manual spin			Height Adjustment method		manual spin lift					
Options:	Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add	
	page 35	varies	not available		not available		not available		Full 933B-xx	\$175	
									Crescent 933CB-xx	\$315	

932 Physician Stool with Manual Spin Lift Adjustment and Brushed Aluminum Base with Toe Caps



[Product Info](#)

[details & photos on page 35*](#)

Options:
[page 35](#) varies

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$467	\$467	\$547	\$567	\$580	\$626	\$651	\$667	\$674	\$709	\$762
Specs:	†		inch	cm	inch		cm	yd		m	
Seat Diameter	16.00		40.64	Base Diameter with casters		24.00	60.96	COM Requirement ¹⁰		1	0.92
Overall Height with backrest -low	32.00		81.28	Seat Height -low		18.50	46.99	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
Overall Height with backrest -high	38.00		96.52	Seat Height -high		22.00	55.88	Shipping Weight ⁸		23	10.43
Height Release hardware	manual spin			Height Adjustment method		manual spin lift					
Options:	Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add	
	page 35	varies	not available		not available		not available		Full 932B-xx	\$175	
									Crescent 932CB-xx	\$315	

931 Physician Stool with Manual Spin Lift Adjustment and Black Nylon Composite Base



[Product Info](#)

[details & photos on page 35*](#)

Options:
[page 35](#) varies

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$382	\$382	\$462	\$482	\$495	\$541	\$566	\$582	\$589	\$624	\$677
Specs:	†		inch	cm	inch		cm	yd		m	
Seat Diameter	16.00		40.64	Base Diameter with casters		24.00	60.96	COM Requirement ¹⁰		1	0.92
Overall Height with backrest -low	33.00		83.82	Seat Height -low		19.50	49.53	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
Overall Height with backrest -high	39.00		99.06	Seat Height -high		23.00	58.42	Shipping Weight ⁸		22	9.98
Height Release hardware	manual spin			Height Adjustment method		manual spin lift					
Options:	Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add	
	page 35	varies	not available		not available		not available		Full 931B-xx	\$175	
									Crescent 931CB-xx	\$315	

928 Urethane Seat Industrial Stool with Footring and Black Tubular 5-Leg Steel Base with Glides



[Product Info](#)

[details & photos on page 35*](#)

Options:
[page 35](#) varies

	list										
	\$532	i available in black urethane only									
Specs:	†		inch	cm	inch		cm	Height Release hardware		single lever	
Seat Diameter	14.00		35.56	Footring		20.00	50.80	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76
Overall Height -low	15.00		38.10	Base Diameter with glides		24.00	60.96	Shipping Weight ⁸		23	10.43
Overall Height -high	19.75		50.17	Height Adjustment method		pneumatic					
Options:	Adjustable Height Seat ^{††12}		add	inch		cm	inch		cm		
	5 [five] [§]	928-xx	\$0	Seat Height -low		15.00	38.10	Seat Height -high		19.75	50.17
	7 [seven]	928-7-xx	\$38	Seat Height -low		16.50	41.91	Seat Height -high		23.25	59.06
	10 [ten]	928-10-xx	\$48	Seat Height -low		20.50	52.07	Seat Height -high		31.00	78.74
Options:	T-Style Armrest [†]		add	Procedure Arm [†]		add	Footring [†]		add	Backrest [†]	
	not available			928PA-xx		\$450	standard		\$0	not available	

878 Lab Stool with D-Handle Release, Footring, and Black Tubular 5-Leg Steel Base with Glides



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$831	\$831	\$911	\$931	\$944	\$990	\$1,015	\$1,031	\$1,038	\$1,073	\$1,126

Specs:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	38.00	96.52	Backrest Height	15.75	45.72	Height Release hardware	D-ring	
Overall Height with backrest -high	48.00	121.92	Base Diameter with glides	24.00	60.96		lb	kg
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸	32	14.51

878 Options:

[details & photos on page 35*](#)

[Product Info](#)

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm		
7 [seven]	878-7-xx	\$0	19.25	48.90	Seat Height-high	26.50	67.31
10 [ten] [§]	878-xx standard	\$0	23.50	59.69	Seat Height-high	34.00	86.36

T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
not available		not available		standard	\$0	Crescent 878CB-xx	\$315

873 Lab Stool with D-Handle Release, Footring, and Polished Aluminum Base



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$831	\$831	\$911	\$931	\$944	\$990	\$1,015	\$1,031	\$1,038	\$1,073	\$1,126

Specs:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	37.00	93.98	Backrest Height	15.75	45.72	Height Release hardware	D-ring	
Overall Height with backrest -high	47.00	119.38	Base Diameter with glides	24.00	60.96		lb	kg
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸	28	12.70

873 Options:

[details & photos on page 35*](#)

[Product Info](#)

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm		
7 [seven]	873-7-xx	\$0	20.50	52.07	Seat Height-high	27.50	69.85
10 [ten] [§]	873-xx standard	\$0	24.50	62.23	Seat Height-high	35.00	88.90

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		standard	\$0	Crescent 873CB-xx	\$315

872 Lab Stool with D-Handle Release, Footring, and Brushed Aluminum Base with Toe Caps



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$864	\$864	\$944	\$964	\$977	\$1,023	\$1,048	\$1,064	\$1,071	\$1,106	\$1,159

Specs:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92
Overall Height with backrest -low	36.00	91.44	Backrest Height	15.75	45.72	Height Release hardware	D-ring	
Overall Height with backrest -high	46.00	116.84	Base Diameter with glides	24.00	60.96		lb	kg
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸	28	12.70

872 Options:

[details & photos on page 35*](#)

[Product Info](#)

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm		
7 [seven]	872-7-xx	\$0	20.50	52.07	Seat Height-high	27.50	69.85
10 [ten] [§]	872-xx standard	\$0	23.50	59.69	Seat Height-high	34.00	86.36

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		standard	\$0	Crescent 872CB-xx	\$315

871 Lab Stool with D-Handle Release, Footring, and Black Nylon Composite Base



Product Info

details & photos on page 35*

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$786	\$786	\$866	\$886	\$899	\$945	\$970	\$986	\$993	\$1,028	\$1,081

Specs:	inch		cm		inch		cm		yd		m	
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92				
Overall Height with backrest -low	37.00	93.98	Backrest Height	15.75	45.72	Height Release hardware	D-ring					
Overall Height with backrest -high	47.00	119.38	Base Diameter with glides	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED	lb		kg			
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic			350	158.76				
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸	27	12.25				

871 Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm	inch	cm
7 [seven]	871-7-xx	\$0	19.50	49.53	Seat Height -low	27.50	69.85
10 [ten] [§]	871-xx standard	\$0	24.50	62.23	Seat Height -low	35.00	88.90

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		standard	\$0	Crescent 871CB-xx	\$315

868 Lab Stool with Single Lever Release, Footring, 5 Leg Black Tubular Steel Base, and Glides



Product Info

details & photos on page 35*

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$827	\$827	\$907	\$927	\$940	\$986	\$1,011	\$1,027	\$1,034	\$1,069	\$1,122

Specs:	inch		cm		inch		cm		yd		m	
Seat Diameter	18.00	45.72	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92				
Overall Height with backrest -low	38.00	96.52	Backrest Height	15.75	45.72	Height Release hardware	single lever					
Overall Height with backrest -high	48.00	121.92	Base Diameter with glides	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED	lb		kg			
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic			350	158.76				
Footring	20.00	50.80	Backrest Adjustment direction	lateral		Shipping Weight ⁸	32	14.51				

868 Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm	inch	cm
7 [seven] [§]	868-7-xx	\$0	19.50	49.53	Seat Height -low	27.50	69.85
10 [ten]	868-xx standard	\$0	23.50	59.69	Seat Height -low	34.00	86.36

T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
not available		not available		standard	\$0	Crescent 868CB-xx	\$315

CB | optional crescent backrest* | standard full backrest



863 Lab Stool with Single Lever Release, Footring, and Polished Aluminum Base



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$827	\$827	\$907	\$927	\$940	\$986	\$1,011	\$1,027	\$1,034	\$1,069	\$1,122

Specs:	†		inch		cm		inch		cm		yd		m		
	inch	cm	inch	cm	inch	cm	inch	cm	yd	m	yd	m	yd	m	
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92							
Overall Height with backrest -low	37.00	93.98	Backrest Height	15.75	45.72	Height Release hardware	single lever								
Overall Height with backrest -high	47.00	119.38	Base Diameter with glides	24.00	60.96			lb	kg						
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76							
Footring	18.00	45.72	Backrest Adjustment direction	lateral		Shipping Weight ⁸	32	14.51							

863 Options:
[details & photos on page 35*](#)

Adjustable Height Seat ^{††12}			add	inch		cm		inch		cm		
7 [seven]	863-7-xx	\$0	Seat Height -low	20.50	52.07	Seat Height -high	27.50	69.85				
10 [ten] [§]	863-xx standard	\$0	Seat Height -low	24.50	62.23	Seat Height -high	35.00	88.90				

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		standard	\$0	Crescent 863CB-xx	\$315

[Product Info](#)

862 Lab Stool with Single Lever Release, Footring, and Brushed Aluminum Base with Toe Caps



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$856	\$856	\$936	\$956	\$969	\$1,015	\$1,040	\$1,056	\$1,063	\$1,098	\$1,151

Specs:	†		inch		cm		inch		cm		yd		m		
	inch	cm	inch	cm	inch	cm	inch	cm	yd	m	yd	m	yd	m	
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92							
Overall Height with backrest -low	36.00	91.44	Backrest Height	15.75	45.72	Height Release hardware	single lever								
Overall Height with backrest -high	46.00	116.84	Base Diameter with glides	24.00	60.96			lb	kg						
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76							
Footring	18.00	45.72	Backrest Adjustment direction	lateral		Shipping Weight ⁸	28	12.70							

862 Options:
[details & photos on page 35*](#)

Adjustable Height Seat ^{††12}			add	inch		cm		inch		cm		
7 [seven]	862-7-xx	\$0	Seat Height -low	19.50	49.53	Seat Height -high	27.50	69.85				
10 [ten] [§]	862-xx standard	\$0	Seat Height -low	23.50	59.69	Seat Height -high	34.00	86.36				

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		standard	\$0	Crescent 862CB-xx	\$315

[Product Info](#)

861 Lab Stool with Single Lever Release, Footring, and Black Nylon Composite Base



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$779	\$779	\$859	\$879	\$892	\$938	\$963	\$979	\$986	\$1,021	\$1,074

Specs:	†		inch		cm		inch		cm		yd		m		
	inch	cm	inch	cm	inch	cm	inch	cm	yd	m	yd	m	yd	m	
Seat Diameter	16.00	40.64	Backrest Width	9.75	24.77	COM Requirement ¹⁰	1	0.92							
Overall Height with backrest -low	37.00	93.98	Backrest Height	15.75	45.72	Height Release hardware	single lever								
Overall Height with backrest -high	47.00	119.38	Base Diameter with glides	24.00	60.96			lb	kg						
Seat Foam thickness	3.00	7.62	Height Adjustment method	pneumatic		Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76							
Footring	18.00	45.72	Backrest Adjustment direction	lateral		Shipping Weight ⁸	27	12.25							

861 Options:
[details & photos on page 35*](#)

Adjustable Height Seat ^{††12}			add	inch		cm		inch		cm		
7 [seven]	861-7-xx	\$0	Seat Height -low	20.50	52.07	Seat Height -high	27.50	69.85				
10 [ten] [§]	861-xx standard	\$0	Seat Height -low	24.50	62.23	Seat Height -high	35.00	88.90				

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	not available		not available		standard	\$0	Crescent 861CB-xx	\$315

[Product Info](#)

848DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Black Tubular 5-Leg Steel Base with Glides



[details & photos on page 35*](#)

[Product Info](#)

list
\$722

i available in black urethane only

Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	37.00	93.98	Base Diameter with casters	24.00	60.96	Height Adjustment seat & back tilt	dual lever			
Overall Height with backrest -high	47.50	120.65	Height Adjustment method	pneumatic				lb	kg	
Footring	20.00	50.80	Backrest Adjustment direction	height			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	45	20.41		

848DC Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm
7 [seven]	848DC-7-xx	\$0		Seat Height-low	17.50 44.45
10 [ten] [§]	848DC-xx standard	\$0		Seat Height-low	21.00 53.34

T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
848DC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

843DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Polished Aluminum Base



[details & photos on page 35*](#)

[Product Info](#)

list
\$733

i available in black urethane only

Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96	Height Adjustment seat & back tilt	dual lever			
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic				lb	kg	
Footring	20.00	50.80	Backrest Adjustment direction	height			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	45	20.41		

843DC Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm
7 [seven]	843DC-7-xx	\$0		Seat Height-low	18.50 46.99
10 [ten] [§]	843DC-xx standard	\$0		Seat Height-low	22.50 57.15

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	843DC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

842DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Brushed Aluminum Base with Toe Caps



[details & photos on page 35*](#)

[Product Info](#)

list
\$766

i available in black urethane only

Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	37.50	97.79	Base Diameter with castes	24.00	60.96	Height Adjustment seat & back tilt	dual lever			
Overall Height with backrest -high	48.00	124.46	Height Adjustment method	pneumatic				lb	kg	
Footring	20.00	50.80	Backrest Adjustment direction	height			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	45	20.41		

842DC Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm
7 [seven]	842DC-7-xx	\$0		Seat Height-low	17.50 44.45
10 [ten] [§]	842DC-xx standard	\$0		Seat Height-low	21.50 54.61

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	842DC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

841DC Self Skin Ergonomic Chair with Seat Tilt, Back Tilt, Footring, and Black Nylon Composite Base



[Product Info](#)

[details & photos on page 35*](#)

list
\$726

available in black urethane only

Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96	Height Adjustment seat & back tilt	dual lever			
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic			lb	kg		
Footring	20.00	50.80	Backrest Adjustment direction	height			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	43	19.50		

841DC Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm	inch	cm		
7 [seven]	841DC-7-xx	\$0		Seat Height-low	18.50	46.99	Seat Height-high	25.50	57.15
10 [ten] [§]	841DC-xx standard	\$0		Seat Height-low	22.50	57.15	Seat Height-high	33.00	83.82

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	841DC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

848SC Self Skin Ergonomic Chair with Footring and Black Tubular 5-Leg Steel Base with Glides



[Product Info](#)

[details & photos on page 35*](#)

list
\$660

available in black urethane only

Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	37.50	95.25	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever			
Overall Height with backrest -high	47.50	120.65	Height Adjustment method	pneumatic			lb	kg		
Footring	20.00	50.80	Backrest Adjustment direction	lateral			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	39	17.69		

848SC Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm	inch	cm		
7 [seven]	848DC-7-xx	\$0		Seat Height-low	17.50	44.45	Seat Height-high	24.50	62.23
10 [ten] [§]	848DC-xx standard	\$0		Seat Height-low	21.00	53.34	Seat Height-high	31.50	80.01

T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
848SC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

843SC Self Skin Ergonomic Chair with Footring and Polished Aluminum Base



[Product Info](#)

[details & photos on page 35*](#)

list
\$670

available in black urethane only

Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever			
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic			lb	kg		
Footring	20.00	50.80	Backrest Adjustment direction	lateral			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	39	17.69		

843SC Options:

Adjustable Height Seat ^{††12}	add	inch	cm	inch	cm	inch	cm		
7 [seven]	843SC-7-xx	\$0		Seat Height-low	18.50	46.99	Seat Height-high	25.50	57.15
10 [ten] [§]	843SC-xx standard	\$0		Seat Height-low	22.50	57.15	Seat Height-high	33.00	83.82

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	843SC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

842SC Self Skin Ergonomic Chair with Footring and Brushed Aluminum Base with Toe Caps



[Product Info](#)

list
\$704

available in black urethane only

Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	37.50	95.25	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever			
Overall Height with backrest -high	48.00	121.92	Height Adjustment method	pneumatic			lb	kg		
Footring	20.00	50.80	Backrest Adjustment direction	lateral			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	39	17.69		

842SC Options:
[details & photos on page 35*](#)

Adjustable Height Seat ^{††12}			add		inch		cm	
7 [seven]	842SC-7-xx	\$0	Seat Height -low	17.50	44.45	Seat Height -high	24.50	62.23
10 [ten] [§]	842SC-xx standard	\$0	Seat Height -low	21.50	54.61	Seat Height -high	32.00	81.28

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	842SC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

841SC Self Skin Ergonomic Chair with Footring, and Black Nylon Composite Base



[Product Info](#)

list
\$598

available in black urethane only


Specs:	inch		cm		inch		cm			
Seat Width	18.00	45.72	Backrest Width	16.00	40.64	Seat & Back substrate	steel			
Seat Depth	17.00	43.18	Backrest Height	12.00	30.48	Seat & Back cover material	urethane			
Overall Height with backrest -low	38.50	97.79	Base Diameter with casters	24.00	60.96	Height Release hardware	single lever			
Overall Height with backrest -high	49.00	124.46	Height Adjustment method	pneumatic			lb	kg		
Footring	20.00	50.80	Backrest Adjustment direction	lateral			Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76	
						Shipping Weight ⁸	38	17.24		

841SC Options:
[details & photos on page 35*](#)

Adjustable Height Seat ^{††12}			add		inch		cm	
7 [seven]	841SC-7-xx	\$0	Seat Height -low	18.50	46.99	Seat Height -high	25.50	57.15
10 [ten] [§]	841SC-xx standard	\$0	Seat Height -low	22.50	57.15	Seat Height -high	33.00	83.82

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	841SC-T-xx	\$187	not available		standard	\$0	Full standard	\$0

T | adjustable urethane t-style armrests option*



838 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Black Tubular 5-Leg Steel Base with Glides



[Product Info](#)

838 Options:
[details & photos on page 35*](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,024	\$1,024	\$1,104	\$1,124	\$1,137	\$1,183	\$1,208	\$1,224	\$1,231	\$1,266	\$1,319

† Specs:	inch		cm		inch		cm		yd		m	
Seat Width	20.00	50.80	Backrest Width	19.00	48.26	COM Requirement ¹⁰		1	0.92			
Seat Depth	18.00	45.72	Backrest Height	16.00	40.64	Height Release hardware		dual lever				
Overall Height with backrest -low	39.00	99.06	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76			
Overall Height with backrest -high	49.50	125.73	Height Adjustment method	pneumatic		Shipping Weight ⁸		45	20.41			
Footring	20.00	50.80	Backrest Adjustment direction	height								

Adjustable Height Seat ^{††12} add			inch		cm		inch		cm	
7 [seven]	838-7-xx	\$0	Seat Height -low	19.50	49.53	Seat Height -high	26.50	67.31		
10 [ten]	838-xx standard	\$0	Seat Height -low	22.50	57.15	Seat Height -high	33.00	83.82		

T-Style Armrest [†] add	Procedure Arm [†] add	Footring [†] add	Backrest [†] add
838-T-xx \$187	not available	standard \$0	Full standard \$0

833 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Polished Aluminum Base



[Product Info](#)

833 Options:
[details & photos on page 35*](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,024	\$1,024	\$1,104	\$1,124	\$1,137	\$1,183	\$1,208	\$1,224	\$1,231	\$1,266	\$1,319

† Specs:	inch		cm		inch		cm		yd		m	
Seat Width	20.00	50.80	Backrest Width	19.00	48.26	COM Requirement ¹⁰		1	0.92			
Seat Depth	18.00	45.72	Backrest Height	16.00	40.64	Height Release hardware		dual lever				
Overall Height with backrest -low	39.00	99.06	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76			
Overall Height with backrest -high	49.50	125.73	Height Adjustment method	pneumatic		Shipping Weight ⁸		45	20.41			
Footring	20.00	50.80	Backrest Adjustment direction	height								

Adjustable Height Seat ^{††12} add			inch		cm		inch		cm	
7 [seven]	833-7-xx	\$0	Seat Height -low	20.00	50.80	Seat Height -high	27.00	68.58		
10 [ten]	833-xx standard	\$0	Seat Height -low	24.00	60.96	Seat Height -high	34.50	87.63		

Casters add	T-Style Armrest [†] add	Procedure Arm [†] add	Footring [†] add	Backrest [†] add
page 35 varies	833-T-xx \$187	not available	standard \$0	Full standard \$0

832 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Brushed Aluminum Base with Toe Caps



[Product Info](#)

832 Options:
[details & photos on page 35*](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,013	\$1,013	\$1,093	\$1,113	\$1,126	\$1,172	\$1,197	\$1,213	\$1,220	\$1,255	\$1,308

† Specs:	inch		cm		inch		cm		yd		m	
Seat Width	20.00	50.80	Backrest Width	19.00	48.26	COM Requirement ¹⁰		1	0.92			
Seat Depth	18.00	45.72	Backrest Height	16.00	40.64	Height Release hardware		dual lever				
Overall Height with backrest -low	38.00	96.52	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76			
Overall Height with backrest -high	48.50	123.19	Height Adjustment method	pneumatic		Shipping Weight ⁸		45	20.41			
Footring	20.00	50.80	Backrest Adjustment direction	height								

Adjustable Height Seat ^{††12} add			inch		cm		inch		cm	
7 [seven]	832-7-xx	\$0	Seat Height -low	19.00	48.26	Seat Height -high	26.00	66.04		
10 [ten]	832-xx standard	\$0	Seat Height -low	23.00	58.42	Seat Height -high	33.50	85.09		

Casters add	T-Style Armrest [†] add	Procedure Arm [†] add	Footring [†] add	Backrest [†] add
page 35 varies	832-T-xx \$187	not available	standard \$0	Full standard \$0

831 Ergonomic Lab Chair with Seat Tilt, Back Tilt, Footring, and Black Nylon Composite Base



Product Info

Options: 831 details & photos on page 35*

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$950	\$950	\$1,030	\$1,050	\$1,063	\$1,109	\$1,134	\$1,150	\$1,157	\$1,192	\$1,245

Specs:	inch		cm		inch		cm		yd		m	
Seat Width	20.00	50.80	Backrest Width	19.00	48.26	COM Requirement ¹⁰		1		0.92		
Seat Depth	18.00	45.72	Backrest Height	16.00	40.64	Height Release hardware		dual lever				
Overall Height with backrest -low	39.00	99.06	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		350		158.76		
Overall Height with backrest -high	49.50	125.73	Height Adjustment method	pneumatic		Shipping Weight ⁸		45		20.41		
Footring	20.00	50.80	Backrest Adjustment direction	height								

Adjustable Height Seat ^{††12}			add		inch		cm		inch		cm	
7 [seven]	831-7-xx	\$0	Seat Height -low	20.00	50.80	Seat Height -high		27.00		68.58		
10 [ten]	831-xx standard	\$0	Seat Height -low	24.00	60.96	Seat Height -high		34.50		87.63		

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	831-T-xx	\$187	not available		standard	\$0	Full standard	\$0

810 SERIES

Laboratory Seating

818 Ergonomic Lab Chair with Footring, and Black Tubular 5-Leg Steel Base with Glides



Product Info

Options: 818 details & photos on page 35*

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$891	\$891	\$971	\$991	\$1,004	\$1,050	\$1,075	\$1,091	\$1,098	\$1,133	\$1,186

Specs:	inch		cm		inch		cm		yd		m	
Seat Width	20.00	50.80	Backrest Width	16.00	40.64	COM Requirement ¹⁰		1		0.92		
Seat Depth	18.00	45.72	Backrest Height	18.00	45.72	Height Release hardware		single lever				
Overall Height with backrest -low	41.50	105.41	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		350		158.76		
Overall Height with backrest -high	52.00	132.08	Height Adjustment method	pneumatic		Shipping Weight ⁸		45		20.41		
Footring	20.00	50.80	Backrest Adjustment direction	lateral								

Adjustable Height Seat ^{††12}			add		inch		cm		inch		cm	
7 [seven]	818-7-xx	\$0	Seat Height -low	19.50	49.53	Seat Height -high		26.50		67.31		
10 [ten]	818-xx standard	\$0	Seat Height -low	22.50	57.15	Seat Height -high		33.00		83.82		

T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
818-T-xx	\$187	not available		standard	\$0	Full standard	\$0

813 Ergonomic Lab Chair with Footring, and Polished Aluminum Base



Product Info

Options: 813 details & photos on page 35*

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$924	\$924	\$1,004	\$1,024	\$1,037	\$1,083	\$1,108	\$1,124	\$1,131	\$1,166	\$1,219

Specs:	inch		cm		inch		cm		yd		m	
Seat Width	20.00	50.80	Backrest Width	16.00	40.64	COM Requirement ¹⁰		1		0.92		
Seat Depth	18.00	45.72	Backrest Height	18.00	45.72	Height Release hardware		single lever				
Overall Height with backrest -low	42.00	106.68	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED		350		158.76		
Overall Height with backrest -high	52.50	133.35	Height Adjustment method	pneumatic		Shipping Weight ⁸		43		19.50		
Footring	20.00	50.80	Backrest Adjustment direction	lateral								

Adjustable Height Seat ^{††12}			add		inch		cm		inch		cm	
7 [seven]	813-7-xx	\$0	Seat Height - low	20.00	50.80	Seat Height -high		27.00		68.58		
10 [ten]	813-xx standard	\$0	Seat Height -low	24.00	60.96	Seat Height -high		34.50		87.63		

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	813-T-xx	\$187	not available		standard	\$0	Full standard	\$0

812 Ergonomic Lab Chair with Footring and Brushed Aluminum Base with Toe Caps



[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$938	\$938	\$1,018	\$1,038	\$1,051	\$1,097	\$1,122	\$1,138	\$1,145	\$1,180	\$1,233

Specs: †

	inch	cm		inch	cm		yd	m
Seat Width	20.00	50.80	Backrest Width	16.00	40.64	COM Requirement ¹⁰	1	0.92
Seat Depth	18.00	45.72	Backrest Height	18.00	45.72	Height Release hardware	single lever	
Overall Height with backrest -low	41.00	104.14	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Overall Height with backrest -high	51.50	130.81	Height Adjustment method	pneumatic		Shipping Weight ⁸	43	19.50
Footring	20.00	50.80	Backrest Adjustment direction	lateral				

Options: **812**

[details & photos on page 35*](#)

Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	812-7-xx	\$0		19.00	48.26	Seat Height -high	26.00	66.04	
10 [ten] [§]	812-xx standard	\$0		23.00	58.42	Seat Height -high	33.50	85.09	

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	812-T-xx	\$187	not available		standard	\$0	Full standard	\$0

811 Ergonomic Lab Chair with Footring, and Black Nylon Composite Base



[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$832	\$832	\$912	\$932	\$945	\$991	\$1,016	\$1,032	\$1,039	\$1,074	\$1,127

Specs: †

	inch	cm		inch	cm		yd	m
Seat Width	20.00	50.80	Backrest Width	16.00	40.64	COM Requirement ¹⁰	1	0.92
Seat Depth	18.00	45.72	Backrest Height	18.00	45.72	Height Release hardware	single lever	
Overall Height with backrest -low	42.00	106.68	Base Diameter with casters	24.00	60.96	Weight Capacity ¹ EVENLY DISTRIBUTED	350	158.76
Overall Height with backrest -high	52.50	133.35	Height Adjustment method	pneumatic		Shipping Weight ⁸	42	19.05
Footring	20.00	50.80	Backrest Adjustment direction	lateral				

Options: **811**

[details & photos on page 35*](#)

Adjustable Height Seat ^{††12}	add	inch	cm	Seat Height -low	inch	cm	Seat Height -high	inch	cm
7 [seven]	811-7-xx	\$0		19.50	49.53	Seat Height -high	26.50	67.31	
10 [ten] [§]	811-xx standard	\$0		24.00	60.96	Seat Height -high	34.50	87.63	

Casters	add	T-Style Armrest [†]	add	Procedure Arm [†]	add	Footring [†]	add	Backrest [†]	add
page 35	varies	811-T-xx	\$187	not available		standard	\$0	Full standard	\$0



xx1 | black nylon composite base

xx2 | brushed aluminum base with toe caps



xx3 | polished aluminum base

xx8 | black tubular five-leg steel base and glides



748 Hydraulic Surgical Stool with Silver Column, Foot Release and Polished Aluminum Base



[Product Info](#)

[details & photos on page 35*](#)

Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,178	\$1,178	\$1,258	\$1,278	\$1,291	\$1,337	\$1,362	\$1,378	\$1,385	\$1,420	\$1,473

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height - low	35.50	90.17	Height Adjustment method	hydraulic		+ FREIGHT * Weight Capacity ¹ EVENLY DISTRIBUTED	lb	kg
Overall Height - high	42.50	107.95	Rotation Lock	on seat			350	158.76
Height Release hardware	foot control		Powder Coat Finish	silver		Shipping Weight ⁸	45	20.41

Options: 748

Casters	add	T-Style Armrest [‡]	add	Procedure Arm [‡]	add	Footring [‡]	add	Backrest [‡]	add
page 35	varies	not available		not available		not available		Full 748B-xx	\$175
								Crescent 748CB-xx	\$315

744 Hydraulic Surgical Stool with Black Column, Foot Release and Black Aluminum Base



[Product Info](#)

[details & photos on page 35*](#)

Specs: †

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,178	\$1,178	\$1,258	\$1,278	\$1,291	\$1,337	\$1,362	\$1,378	\$1,385	\$1,420	\$1,473

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement ¹⁰	1	0.92
Overall Height - low	35.50	90.17	Height Adjustment method	hydraulic		+ FREIGHT * Weight Capacity ¹ EVENLY DISTRIBUTED	lb	kg
Overall Height - high	42.50	107.95	Rotation Lock	on seat			350	158.76
Height Release hardware	foot control		Powder Coat Finish	black		Shipping Weight ⁸	45	20.41

Options: 744

Casters	add	T-Style Armrest [‡]	add	Procedure Arm [‡]	add	Footring [‡]	add	Backrest [‡]	add
page 35	varies	not available		not available		not available		Full 744B-xx	\$175
								Crescent 744CB-xx	\$315

HSS500 Surgeon Stool with Foot-Controlled Hydraulics, 3-Spindle Base, and 4 Dual Tente® Casters



[Product Info](#)

list
\$6,420

[i](#) available in black urethane only

Specs: †

	inch	cm		inch	cm		degrees	
Overall Width ³	25.50	64.77	Base Width ³	24.00	60.96	Armrest Rotation	360	
Overall Depth ³	32.25	81.92	Base Depth ³	28.50	72.39	Armrest Support	180	
Overall Height - low	21.00	53.34	Seat Height - low	21.50	54.61	Backrest Tilt	20	
Overall Height - high	28.00	71.12	Seat Height - high	28.50	72.39	+ FREIGHT * Weight Capacity ¹ EVENLY DISTRIBUTED	lb	kg
Armrest Width ³	4.75	12.07	Backrest Width	14.00	35.56		300	136.08
Armrest Length ³	12.00	30.48	Backrest Length	8.25	20.96	Shipping Weight ⁸	74	33.57
Armrest Height - low	7.00	17.78	Backrest Height - low	10.50	26.67	Tente® Casters four sets - dual	3.00	7.62
Armrest Height - high	15.00	38.10	Backrest Height - high	12.50	31.75			

HSS500 Options:

Saddle Seat	add	Wedge Seat	add	Water Drop Seat	add	Reverse Water Drop	add
HSS500-S	\$0	HSS500-W	\$0	HSS500-WD	\$0	HSS500-RWD	\$0

³ variance: ± 0.5 inch

IN601P Pivot-Arm Clinical Care Recliner with 3 Positions, Trendelenburg, and Upholstered Arms



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$4,493	\$4,493	\$4,972	\$5,092	\$5,172	\$5,448	\$5,599	\$5,692	\$5,733	\$5,943	\$6,263

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	34.50	87.63	Seat Width	22.50	57.15	COM Requirement ¹⁰	7	6.40
Overall Width side table up	54.00	137.16	Seat Depth	19.00	48.26		inch	cm†
Overall Depth	38.00	96.52	Seat Height	22.00	55.88	Fold Down Side Tables	10 x 19	25 x 48
Overall Depth trendelenburg	72.00	182.88	Arm Width	3.00	7.62		lb	kg
Overall Height	51.50	130.81	Width between arms	25.00	63.50	Weight Capacity ¹ EVENLY DISTRIBUTED	400	181.44
Back Width	30.50	77.47	Arm Height from floor	28.50	72.39		Shipping Weight ⁸	197
Back Height above seat	30.00	76.20	Tente® Casters	5.00	12.70			

Product Info

details & photos on page 34*

IN601 Options:

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended†	add	Folding Tables	add
None IN601Pxx	\$0	Mount FR601IVMNT ⁵	\$64	Skirt HRSxx	\$78	Legrest IN601PLxx ^{6S}	\$0	Fold Down Side [§]	\$0
Both IN601PHxx	\$1,740	Pole 2101A	\$64	Cover HRC601xx	\$149	Legs	not available	Swing-Away SAFT20	\$447
Heat Only IN601PH0xx	\$1,494			Pillow FR5659xxUS	\$98				

IN601U Pivot-Arm Clinical Care Recliner with 3 Positions, Trendelenburg, and Urethane Capped Arms



com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$4,666	\$4,666	\$5,145	\$5,265	\$5,345	\$5,621	\$5,772	\$5,865	\$5,906	\$6,116	\$6,436

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	34.50	87.63	Seat Width	22.50	57.15	COM Requirement ¹	7	6.40
Overall Width side table up	54.00	137.16	Seat Depth	19.00	48.26		inch	cm†
Overall Depth	38.00	96.52	Seat Height	22.00	55.88	Fold Down Side Tables	10 x 19	25 x 48
Overall Depth trendelenburg	72.00	182.88	Arm Width	3.00	7.62		lb	kg
Overall Height	51.50	130.81	Width between arms	25.00	63.50	Weight Capacity ¹⁰ EVENLY DISTRIBUTED	400	181.44
Back Width	30.50	77.47	Arm Height from floor	28.50	72.39		Shipping Weight ⁸	197
Back Height above seat	30.00	76.20	Tente® Casters	5.00	12.70			

Product Info

details & photos on page 34*

IN601U Options:

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended†	add	Folding Table	add
None IN601Uxx	\$0	Mount FR601IVMNT ⁵	\$64	Skirt HRSxx	\$78	Legrest IN601ULxx ⁶	\$0	Fold Down Side [§]	\$0
Both IN601UHxx	\$1,740	Pole 2101A	\$64	Cover HRC601xx	\$149	Legs	not available	Swing-Away SAFT20	\$447
Heat Only IN601UH0xx	\$1,494			Pillow FR5659xxUS	\$98				



IN597G Ortho-Biotic II Recliner with Infinite Back Positioning, Trendelenburg, and Wood Capped Arms



[Product Info](#)

[details & photos on page 34*](#)

Heat & Massage [†]		add
None	IN597Gxx	\$0
Both	IN597GHxx	\$1,851
Heat Only	IN597GH0xx	\$1,605

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$4,789	\$4,789	\$5,268	\$5,388	\$5,468	\$5,744	\$5,895	\$5,988	\$6,029	\$6,239	\$6,559
	inch		cm		inch		cm		yd		m
Overall Width	33.00	83.82	Seat Width	23.00	58.42	COM Requirement ¹⁰		6	5.49		
Overall Width <i>side table up</i>	45.00	114.30	Seat Depth	20.00	50.80			inch	cm [†]		
Overall Depth	37.00	93.98	Seat Height	21.00	53.34	Fold Down Side Table ⁴⁵		10 x 19	25 x 48		
Overall Depth <i>trendelenburg</i>	73.00	185.42	Arm Width	3.50	8.89						
Overall Height	46.00	116.84	Width <i>between arms</i>	24.00	60.96	+ FREIGHT*		lb	kg		
Back Width	23.00	58.42	Arm Height <i>from floor</i>	26.75	67.95	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76		
Back Height <i>above seat</i>	30.00	76.20	Tente® Casters	4.00	10.16	Shipping Weight ⁸		210	95.25		

IN597G Options:

IV Mounting	add	Headrest [†]	add	Extended	Folding Table	add
Mount FR597GIVMNT ⁵	\$48	Skirt HRSxx	\$78	Legrest not available	Left Side Single FR597GTBLMNT ⁴⁵	\$96
Pole 2101A	\$64	Cover not available		Legs not available	Swing-Away SAFT20	\$447
		Pillow standard [§]	\$0			

IN597P Ortho-Biotic II Recliner with Infinite Back Positioning, Trendelenburg, and Upholstered Arms



[Product Info](#)

[details & photos on page 34*](#)

Heat & Massage [†]		add
None	IN597Pxx	\$0
Both	IN597PHxx	\$1,851
Heat Only	IN597PH0xx	\$1,605

Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$4,646	\$4,646	\$5,125	\$5,245	\$5,325	\$5,601	\$5,752	\$5,845	\$5,886	\$6,096	\$6,416
	inch		cm		inch		cm		yd		m
Overall Width	33.00	83.82	Seat Width	23.00	58.42	COM Requirement ¹⁰		6	5.49		
Overall Width <i>side table up</i>	45.00	114.30	Seat Depth	20.00	50.80			inch	cm [†]		
Overall Depth	37.00	93.98	Seat Height	21.00	53.34	Fold Down Side Table		10 x 19	25 x 48		
Overall Depth <i>trendelenburg</i>	73.00	185.42	Arm Width	3.50	8.89						
Overall Height	46.00	116.84	Width <i>between arms</i>	24.00	60.96	+ FREIGHT*		lb	kg		
Back Width	23.00	58.42	Arm Height <i>from floor</i>	26.75	67.95	Weight Capacity ¹ EVENLY DISTRIBUTED		350	158.76		
Back Height <i>above seat</i>	30.00	76.20	Tente® Casters	4.00	10.16	Shipping Weight ⁸		210	95.25		

IN597P Options:

IV Mounting	add	Headrest [†]	add	Extended	Folding Table	add
Mount FR597GIVMNT ⁵	\$48	Skirt HRSxx	\$78	Legrest not available	Left Side Single FR597GTBLMNT ⁴⁵	\$96
Pole 2101A	\$64	Cover not available		Legs not available	Swing-Away SAFT20	\$447
		Pillow standard [§]	\$0			



IN587W Extra Wide Clinical Care Recliner with 3 Positions and Trendelenburg



Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$3,592	\$3,592	\$4,071	\$4,191	\$4,271	\$4,547	\$4,698	\$4,791	\$4,832	\$5,042	\$5,362
	inch		cm		inch		cm		yd		m
Overall Width	33.00	83.82	Seat Width		24.00	60.96	COM Requirement¹⁰		6	5.49	
Overall Width side table up	53.50	135.89	Seat Depth		20.50	52.07			inch		cm†
Overall Depth	45.50	115.57	Seat Height		22.00	55.88	Fold Down Side Table		10 x 19	25 x 48	
Overall Depth trendelenburg	78.00	198.12	Arm Width		4.50	11.43					
Overall Height	50.00	127.00	Width between arms		23.25	59.06	+ FREIGHT *		lb		kg
Back Width	30.50	77.47	Arm Height from floor		28.00	71.39	Weight Capacity¹ EVENLY DISTRIBUTED		450	204.12	
Back Height above seat	33.50	85.09	Tente® Casters		5.00	12.70	Shipping Weight⁸		143	64.86	

IN587W Options:

Product Info

details & photos on page 34*

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended	add	Folding Table	add
None IN587Wxx	\$0	Mount 5782RG	\$54	Skirt HRSxx	\$78	Legrest	not available	Fold Down Side ⁵	\$0
Both IN587WHxx	\$1,740	Pole 2101A	\$64	Cover HRC587xx	\$128	4 inch ¹¹ FR5699G	\$40	Swing-Away SAFT10	\$447
Heat Only IN587WH0xx	\$1,493			Pillow FR5659xxUS	\$98				

IN587WD Extra Wide Drop-Arm Clinical Care Recliner with 3 Positions and Trendelenburg



Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$3,900	\$3,900	\$4,379	\$4,499	\$4,579	\$4,855	\$5,006	\$5,099	\$5,140	\$5,350	\$5,670
	inch		cm		inch		cm		yd		m
Overall Width	35.75	90.81	Seat Width		24.00	60.96	COM Requirement¹⁰		6	5.49	
Overall Width side table up	58.25	147.96	Seat Depth		20.50	52.07			inch		cm†
Overall Depth	44.00	111.76	Seat Height		22.00	55.88	Fold Down Side Table		10 x 19	25 x 48	
Overall Depth trendelenburg	78.00	198.12	Arm Width		4.50	11.43	Tente® Casters		5.00	12.70	
Overall Height	50.00	127.00	Width between arms		26.50	67.31	+ FREIGHT *		lb		kg
Back Width	30.50	77.47	Drop Arm Height from floor, arm up		27.50	69.85	Weight Capacity¹ EVENLY DISTRIBUTED		450	204.12	
Back Height above seat	33.50	85.09	Drop Arm Height from floor, arm down		21.50	54.61	Shipping Weight⁸		145	65.77	

IN587WD Options:

Product Info

details & photos on page 34*

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended	add	Folding Table	add
None IN587WDxx	\$0	Mount 5782RG	\$54	Skirt HRSxx	\$78	Legrest	not available	Fold Down Side ⁵	\$0
Both IN587WDHxx	\$1,740	Pole 2101A	\$64	Cover HRC587xx	\$128	Legs	not available	Swing-Away SAFT10	\$447
Heat Only IN587WDH0xx	\$1,493			Pillow FR5659xxUS	\$98				



IN566G Deluxe Clinical Care Recliner with 3 Positions and Trendelenburg



[Product Info](#)

[details & photos on page 34*](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$3,086	\$3,086	\$3,565	\$3,685	\$3,765	\$4,041	\$4,192	\$4,285	\$4,326	\$4,536	\$4,856

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	31.25	79.38	Seat Width	20.50	52.07	COM Requirement ¹⁰	6	5.49
Overall Width side table up	51.50	130.81	Seat Depth	19.00	48.26			
Overall Depth	37.50	95.25	Seat Height	22.25	56.52	Fold Down Side Table	10 x 19	25 x 48
Overall Depth trendelenburg	75.00	190.50	Arm Width	4.50	11.43			
Overall Height	49.50	125.73	Width between arms	21.50	54.61	+ FREIGHT * Weight Capacity ¹ EVENLY DISTRIBUTED	400	181.44
Back Width	27.50	69.85	Arm Height from floor	28.50	72.39			
Back Height above seat	32.00	81.28	Tente® Casters	5.00	12.70	Shipping Weight ⁸	140	64.41

Options:

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended	add	Folding Table	add
None IN566Gxx	\$0	Mount 5782RG	\$54	Skirt HRSxx	\$78	Legrest	not available	Fold Down Side ⁵	\$0
Both IN566GHxx	\$1,666	Pole 2101A	\$64	Cover HRC566xx	\$134	4 inch ¹¹ FR5699G	\$40	Swing-Away SAFT10	\$447
Heat Only IN566GH0xx	\$1,420			Pillow FR5659xxUS	\$98				

IN566DG Deluxe Drop-Arm Clinical Care Recliner with 3 Positions and Trendelenburg



[Product Info](#)

[details & photos on page 34*](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$3,518	\$3,518	\$3,997	\$4,117	\$4,197	\$4,473	\$4,624	\$4,717	\$4,758	\$4,968	\$5,288

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	32.50	82.55	Seat Width	20.50	52.07	COM Requirement ¹⁰	6	5.49
Overall Width side table up	52.50	133.35	Seat Depth	19.00	48.26			
Overall Depth	38.00	96.52	Seat Height	22.50	57.15	Fold Down Side Table	10 x 19	25 x 48
Overall Depth trendelenburg	75.00	190.50	Arm Width	4.50	11.43			
Overall Height	50.50	128.27	Width between arms	22.00	55.88	+ FREIGHT * Weight Capacity ¹ EVENLY DISTRIBUTED	400	181.44
Back Width	27.50	69.95	Drop Arm Height from floor, arm up	28.00	71.12			
Back Height above seat	30.00	76.20	Drop Arm Height from floor, arm down	21.50	54.61	Shipping Weight ⁸	142	64.41

Heat & Massage†	add	IV Mounting	add	Headrest†	add	Extended	add	Folding Table	add
None IN566DGxx	\$0	Mount 5782RG	\$54	Skirt HRSxx	\$78	Legrest	not available	Fold Down Side ⁵	\$0
Both IN566DGHxx	\$1,666	Pole 2101A	\$64	Cover HRC566xx	\$134	Legs	not available	Swing-Away SAFT10	\$447
Heat Only IN566DGH0xx	\$1,420			Pillow FR5659xxUS	\$98				



vinyl: Oatmeal 8602

IN577RG Multi-Care Recliner with 3 Positions and Trendelenburg



Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$2,962	\$2,962	\$3,441	\$3,561	\$3,641	\$3,917	\$4,068	\$4,161	\$4,202	\$4,412	\$4,732
	inch		cm		inch		cm		yd		m
Overall Width	31.50	80.01	Seat Width	21.50	54.61	COM Requirement¹⁰		6	5.49		
Overall Width side table up	52.00	132.08	Seat Depth	19.50	49.53						
Overall Depth	40.50	102.87	Seat Height	20.50	52.07	Fold Down Side Table		10 x 19	25 x 48		
Overall Depth trendelenburg	73.50	186.69	Arm Width	4.50	11.43						
Overall Height	48.00	121.92	Width between arms	21.50	54.61						
Back Width	21.50	54.61	Arm Height from floor	27.50	69.85						
Back Height above seat	34.00	86.36	Tente® Casters	5.00	12.70						
								Weight Capacity¹ EVENLY DISTRIBUTED	350	158.76	
								Shipping Weight⁸	129	58.51	

IN577RG Options:

Product Info

details & photos on page 34*

Heat & Massage [‡]	add	IV Mounting	add	Headrest [‡]	add	Extended	add	Folding Table	add
None IN577RGxx	\$0	Mount FR5782RG	\$54	Skirt HRSxx	\$78	Legrest not available		Fold Down Side [§]	\$0
Both IN577RGHxx	\$1,666	Pole 2101A	\$64	Cover HRC577xx	\$115	4 inch ¹¹ FR5699G	\$40	Swing-Away SAFT10	\$447
Heat Only IN577RGH0xx	\$1,420			Pillow FR5659xxUS	\$98				



600 SERIES

Blood Draw / Phlebotomy Seating

620 Upholstered Bariatric Blood Draw / Phlebotomy Chair with Adjustable Height, Padded, Arm Rest



Specs: †

	com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
	\$1,720	\$1,720	\$1,920	\$1,970	\$2,003	\$2,118	\$2,181	\$2,220	\$2,237	\$2,324	\$2,458
	inch		cm		inch		cm		yd		m
Overall Width	34.00	86.36	Seat Width	24.50	62.23	COM Requirement¹⁰		2.5	2.29		
Overall Depth	29.00	73.66	Seat Depth	18.50	46.99						
Overall Height	37.50	95.25	Seat Height	22.00	55.88	Fold Down Side Table		10 x 19	25 x 48		
Arm Height - low	27.00	62.23	Back Width	25.50	64.77						
Arm Height - high	32.00	81.28	Back Height	17.00	43.18						
						Powder Coat Finish	foghorn				
								Weight Capacity¹ EVENLY DISTRIBUTED	800	362.87	
								Shipping Weight⁸	55	24.95	

620 Options: *

Single Flip Arm [‡]	add	Dual Flip Arms [‡]	add	Wallsaver Wireform Accessory	add
620-xx	\$0	620-DA-xx	\$172	620 WALL SAVER [§] — order separately *	\$66



Product Info



460 High / Low Power Exam Table with Power Assisted Back, Hand Control, and Stirrups



Product Info

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$7,984	\$7,984	\$8,224	\$8,284	\$8,324	\$8,462	\$8,537	\$8,584	\$8,604	\$8,709	\$8,869

Specs: †

	inch	cm					yd	m
Overall Width	27.00	68.58	Pelvic Tilt	standard	Power Assisted			back
Overall Length <i>leg rest extended</i>	72.00	182.88	Pull Out Leg Rest	standard	UL 60601-1			yes
Overall Height ³ <i>- low</i>	18.50	46.99	Treatment Tray	integrated				
Overall Height ³ <i>- high</i>	35.50	90.17	Stirrups	five position continuous lock	COM Requirement ¹⁰	3	2.74	
Paper Roll Holder Widths Available	18.00	45.72	Base	dual pedestal	FREIGHT			
	19.00	48.26	Drive	DC electromechanical	Weight Capacity ¹ EVENLY DISTRIBUTED	600	272.16	
	20.00	50.80	Powder Coat Finish	foghorn	Shipping Weight ⁸	395	179.7	

Options: ★

460 Fully Loaded [‡]	add	460 Foot Switch [‡]	add	460 Assist Arms [‡]	add	460 Physician Storage Drawer [‡]	add
460-FL-xx	\$1,839	460-FS-xx	\$512	460-AA-xx	\$722	460-PS-xx Drawer Only	\$397
460-xx Base Model						460-PS-DW-xx With Warmer	\$608



420 Exam Table with 3 Pass-Through Drawers, 2 Storage Drawers, and Stirrups



Product Info

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,875	\$2,875	\$3,115	\$3,175	\$3,215	\$3,353	\$3,428	\$3,475	\$3,495	\$3,600	\$3,760

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	27.00	68.58	Drawer Width ⁷	18.00	45.72	COM Requirement ¹⁰	3	2.74
Overall Length <i>leg rest extended</i>	72.00	182.88	Drawer Length ⁷	16.50	41.91			
Overall Height	32.00	81.28	Drawer Height ⁷	4.00	10.16			
Foot Step Width	18.00	45.72	Pelvic Tilt	standard	FREIGHT			
Foot Step Length	17.00	43.18	Pull Out Leg Rest	standard	Foot Step Weight Capacity ¹ EVENLY DISTRIBUTED	800	362.87	
Paper Roll Holder Widths Available	18.00	45.72	Treatment Tray	integrated	Table Weight Capacity ¹ EVENLY DISTRIBUTED	500	226.80	
	19.00	48.26	Stirrups	five position continuous lock	Shipping Weight ⁸	335	151.95	
	20.00	50.80	Powder Coat Finish	foghorn				

Options: ★

420 Base Model	add	420 Integrated Electrical Outlet [‡]	add
420-xx \$		420-E0-xx Outlet Only	\$281
		420-DW-xx with Drawer Warmer	\$669



415 Exam Room Treatment Table with Shelf, Adjustable Back, and Contour Top



Product Info

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,332	\$2,332	\$2,572	\$2,632	\$2,672	\$2,810	\$2,885	\$2,932	\$2,952	\$3,057	\$3,217

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	Weight Capacity ¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90						
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn		Shipping Weight ⁸	190	90.72

Options:

No Drawers [‡]	add	Single Pass-Through Drawer [‡]	add	Dual Pass-Through Drawer [‡]	add
415-xx [§]	\$0	415D-xx	\$223	415DD-xx	\$446
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

414 Exam Room Treatment Table with Shelf and Adjustable Back



Product Info

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,311	\$2,311	\$2,551	\$2,611	\$2,651	\$2,789	\$2,864	\$2,911	\$2,931	\$3,036	\$3,196

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	Weight Capacity ¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90						
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn		Shipping Weight ⁸	200	90.72

Options:

No Drawers [‡]	add	Single Pass-Through Drawer [‡]	add	Dual Pass-Through Drawer [‡]	add
414-xx	\$0	414D-xx [§]	\$223	414DD-xx	\$446
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

413 Exam Room Treatment Table with Shelf and Flat Top



Product Info

Specs:

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,956	\$1,956	\$2,196	\$2,256	\$2,296	\$2,434	\$2,509	\$2,556	\$2,576	\$2,681	\$2,841

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹⁰	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45	Weight Capacity ¹ EVENLY DISTRIBUTED	500	226.80
Drawer Length ⁷	19.25	48.90						
Drawer Height ⁷	5.00	12.70	Powder Coat Finish	foghorn		Shipping Weight ⁸	200	90.72

Options:

No Drawers [‡]	add	Single Pass-Through Drawer [‡]	add	Dual Pass-Through Drawer [‡]	add
413-xx	\$0	413D-xx	\$223	413DD-xx [§]	\$446
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

41xD Single Drawer Option:



41xDD Double Drawer Option:



412 Exam Room Treatment Table with Contoured Adjustable Back



Product Info

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,102	\$2,102	\$2,342	\$2,402	\$2,442	\$2,580	\$2,655	\$2,702	\$2,722	\$2,827	\$2,987

Specs:

	†		inch		cm		inch		cm		yd		m	
Overall Width	28.00	71.12					18.00	45.72			COM Requirement ¹⁰		4	3.66
Overall Length	76.00	193.04					19.00	48.26						
Overall Height	32.00	81.28					20.00	50.80						
Drawer Width ⁷	17.75	45.09					1.75	4.45						
Drawer Length ⁷	19.25	48.90												
Drawer Height ⁷	5.00	12.70												
							Powder Coat Finish		foghorn					

412 Options:

No Drawers [†]	add	Single Pass-Through Drawer [†]	add	Dual Pass-Through Drawer [†]	add
412-xx [§]	\$0	412D-xx	\$223	412DD-xx	\$446
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

411 Exam Room Treatment Table with Adjustable Back



Product Info

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$2,079	\$2,079	\$2,319	\$2,379	\$2,419	\$2,557	\$2,632	\$2,679	\$2,699	\$2,804	\$2,964

Specs:

	†		inch		cm		inch		cm		yd		m	
Overall Width	28.00	71.12					18.00	45.72			COM Requirement ¹⁰		4	3.66
Overall Length	76.00	193.04					19.00	48.26						
Overall Height	32.00	81.28					20.00	50.80						
Drawer Width ⁷	17.75	45.09					1.75	4.45						
Drawer Length ⁷	19.25	48.90												
Drawer Height ⁷	5.00	12.70												
							Powder Coat Finish		foghorn					

411 Options:

No Drawers [†]	add	Single Pass-Through Drawer [†]	add	Dual Pass-Through Drawer [†]	add
411-xx	\$0	411D-xx [§]	\$223	411DD-xx	\$446
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb

410 Exam Room Treatment Table with Flat Top



Product Info

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$1,546	\$1,546	\$1,786	\$1,846	\$1,886	\$2,024	\$2,099	\$2,146	\$2,166	\$2,271	\$2,431

Specs:

	†		inch		cm		inch		cm		yd		m	
Overall Width	28.00	71.12					18.00	45.72			COM Requirement ¹		4	3.66
Overall Length	76.00	193.04					19.00	48.26						
Overall Height	32.00	81.28					20.00	50.80						
Drawer Width ⁷	17.75	45.09					1.75	4.45						
Drawer Length ⁷	19.25	48.90												
Drawer Height ⁷	5.00	12.70												
							Powder Coat Finish		foghorn					

410 Options:

No Drawers [†]	add	Single Pass-Through Drawer [†]	add	Dual Pass-Through Drawer [†]	add
410-xx	\$0	410D-xx [§]	\$223	410DD-xx [§]	\$446
Shipping Weight ⁸	+0 lb	Shipping Weight ⁸	+10 lb	Shipping Weight ⁸	+20 lb


41xD Single Drawer Option:



41xDD Double Drawer Option:



RT460 Replacement Top with Foot Pad for Intensa 460 Exam Tables




[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$760	\$760	\$1,000	\$1,060	\$1,100	\$1,238	\$1,313	\$1,360	\$1,380	\$1,485	\$1,645

Specs: †	inch		cm		inch		cm		yd		m	
Overall Width	27.50	69.85	Foot Rest Width	19.00	48.26	COM Requirement ¹⁰ <i>no railroad patterns</i>		3	2.74			
Overall Length	60.00	152.40	Foot Rest Length	18.00	45.72	+ FREIGHT				lb	kg	
Paper Roll Holder	included		Mounting Hardware	included		Shipping Weight ⁸		55	24.95			

Order Number[‡] RT460-xx

RT420 Replacement Top with Foot Pad for Intensa 420 Exam Tables




[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$760	\$760	\$1,000	\$1,060	\$1,100	\$1,238	\$1,313	\$1,360	\$1,380	\$1,485	\$1,645

Specs: †	inch		cm		inch		cm		yd		m	
Overall Width	27.50	69.85	Foot Rest Width	19.00	48.26	COM Requirement ¹⁰ <i>no railroad patterns</i>		3	2.74			
Overall Length	60.00	152.40	Foot Rest Length	18.00	45.72	+ FREIGHT				lb	kg	
Paper Roll Holder	included		Mounting Hardware	included		Shipping Weight ⁸		55	24.95			

Order Number[‡] RT420-xx

RT405 Replacement Top with Foot Pad for 405 Exam Tables




[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$613	\$613	\$853	\$913	\$953	\$1,091	\$1,166	\$1,213	\$1,233	\$1,338	\$1,498

Specs: †	inch		cm		inch		cm		yd		m	
Overall Width	27.50	69.85	Foot Rest Width	20.50	52.07	COM Requirement ¹⁰		3	2.74			
Overall Length	60.00	152.40	Foot Rest Length	17.50	44.45	GROUND				lb	kg	
			Hook & Loop Attachment	included		Shipping Weight ⁸		35	15.88			

Order Number[‡] RT405-xx

RT404 Replacement Top with Foot Pad for 404 Exam Tables




[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$613	\$613	\$853	\$913	\$953	\$1,091	\$1,166	\$1,213	\$1,233	\$1,338	\$1,498

Specs: †	inch		cm		inch		cm		yd		m	
Overall Width	27.50	69.85	Foot Rest Width	20.50	52.07	COM Requirement ¹⁰		3	2.74			
Overall Length	60.00	152.40	Foot Rest Length	17.50	44.45	GROUND				lb	kg	
			Hook & Loop Attachment	included		Shipping Weight ⁸		35	15.88			

Order Number[‡] RT404-xx

RT304 Replacement Top with Foot Pad for 304 Exam Tables



Product Info


com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$856	\$856	\$1,096	\$1,156	\$1,196	\$1,334	\$1,409	\$1,456	\$1,476	\$1,581	\$1,741

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	27.50	69.85	Foot Rest Width	20.00	50.80	COM Requirement ¹⁰	3	2.74
Overall Length	55.00	139.70	Foot Rest Length	18.00	45.72	+ FREIGHT*	lb	kg
						Shipping Weight ⁸	55	24.95

Order Number[‡] RT304-xx

RT204 Replacement Top with Foot Pad for 204 Exam Tables



Product Info

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$881	\$881	\$1,121	\$1,181	\$1,221	\$1,359	\$1,434	\$1,481	\$1,501	\$1,606	\$1,766


Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	27.50	69.85	Foot Rest Width	20.00	50.80	COM Requirement ¹⁰	3	2.74
Overall Length	55.00	139.70	Foot Rest Length	18.00	45.72	+ FREIGHT*	lb	kg
Lift Handles	not included		Top Attachment Knobs	not included		Shipping Weight ⁸	55	24.95

Order Number[‡] RT204-xx

Does not include lift handles and top attachment knobs
If the OEM model has a seamless top, request paper roll holder when ordering

RT104 Replacement Top with Foot Pad for 104 Exam Tables



Product Info


com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$938	\$938	\$1,178	\$1,238	\$1,278	\$1,416	\$1,491	\$1,538	\$1,558	\$1,663	\$1,823

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	27.50	69.85	Foot Rest Width	20.00	50.80	COM Requirement ¹⁰	3	2.74
Overall Length	55.00	139.70	Foot Rest Length	18.00	45.72	+ FREIGHT*	lb	kg
Paper Roll Holder	included		Mounting Hardware	included		Shipping Weight ⁸	55	24.95

Order Number[‡] RT104-xx

RT100 Replacement Top with Foot Pad for 100 Exam Tables



Product Info

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$938	\$938	\$1,178	\$1,238	\$1,278	\$1,416	\$1,491	\$1,538	\$1,558	\$1,663	\$1,823

Specs: †

	inch	cm		inch	cm		yd	m
Overall Width	27.50	69.85	Foot Rest Width	20.00	50.80	COM Requirement ¹⁰	3	2.74
Overall Length	55.00	139.70	Foot Rest Length	18.00	45.72	+ FREIGHT*	lb	kg
Paper Roll Holder	included		Mounting Hardware	included		Shipping Weight ⁸	55	24.95

Order Number[‡] RT100-xx

202 Wallsaver, Stacking, Patient / Guest Chair with Slant Arms



[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$564	\$564	\$644	\$664	\$677	\$723	\$748	\$764	\$771	\$806	\$859

Specs: †	inch		cm		inch		cm		yd		m	
Overall Width	22.50	57.15	Seat Width	18.00	45.72	COM Requirement ¹⁰		1	0.91			
Overall Depth	22.00	55.88	Seat Depth	18.00	45.72							
Overall Height	32.50	82.55	Seat Height	18.00	45.72							
Back Width	18.00	45.72	Arm Width Inside	19.50	49.53							
Back Height	12.50	31.75	Arm Height	26.00	66.04							
Steel Tubing	14 gauge		Powder Coat Finish	black								

Order Number[‡] 202-xx

201 Wallsaver, Stacking, Patient / Guest Chair with Straight Arms



[Product Info](#)

com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$538	\$538	\$618	\$638	\$651	\$697	\$722	\$738	\$745	\$780	\$833

Specs: †	inch		cm		inch		cm		yd		m	
Overall Width	23.50	59.69	Seat Width	18.00	45.72	COM Requirement ¹⁰		1	0.91			
Overall Depth	22.00	55.88	Seat Depth	18.00	45.72							
Overall Height	32.50	82.55	Seat Height	18.00	45.72							
Back Width	18.00	45.72	Arm Width Inside	20.00	50.80							
Back Height	12.50	31.75	Arm Height	26.00	66.04							
Steel Tubing	14 gauge		Powder Coat Finish	black								

Order Number[‡] 201-xx

200 Wallsaver, Stacking, Patient / Guest Chair without Arms



[Product Info](#)


com	tier 1	tier 2	tier 3	tier 4	tier 5	tier 6	tier 7	tier 8	tier 9	tier 10
\$504	\$504	\$584	\$604	\$617	\$663	\$688	\$704	\$711	\$746	\$799

Specs: †	inch		cm		inch		cm		yd		m	
Overall Width	19.50	49.53	Seat Width	18.00	45.72	COM Requirement ¹⁰		1	0.91			
Overall Depth	22.00	55.88	Seat Depth	18.00	45.72							
Overall Height	32.50	82.55	Seat Height	18.00	45.72							
Back Width	21.00	53.34	Steel Tubing	14 gauge								
Back Height	12.50	31.75	Powder Coat Finish	black								

Order Number[‡] 200-xx



179 Overbed Table with Heavy-Duty Steel, Low Profile Base, and Raised Edge Thermofoil Dual Top with Cup Holder



179TFx-y-xx
\$1,676

Specs:

	inch	cm		inch	cm
Overall Width	46.00	116.84	Base Width	32.50	82.85
Overall Depth	19.00	48.26	Base Depth	19.00	48.26
Overall Height - low	28.00	71.12	Base Profile Height	1.94	4.93
Overall Height - high	44.00	111.76	Main Top Width	32.00	81.28
Top Depth	17.00	43.18	Flip Top Width	13.50	34.29

i Solid Surface top available - contact customer service



	lb	kg
Weight Capacity ¹ EVENLY DISTRIBUTED	100	45.36
Shipping Weight ⁸	43	19.50

Options:

Thermofoil Top Shape [‡]	add	Vanity with Mirror [‡]	add	Wire Basket [‡]	add
Rectangle 179TFR-y-xx [§]	\$0	179TFx-y-xx-VM [§]	\$210	179TFx-y-xx-BK [§]	\$210
Kidney 179TFK-y-xx	\$135	<i>vanity with mirror and wire basket 177TFx-y-xx-VM-BK</i>			

[Product Info](#)

[finish options page 39*](#)

177 Overbed Table with Heavy-Duty Steel, Low Profile Base and Laminate Top *or* Raised Edge Thermofoil Top



laminate: 177-y-xx
\$1,460

Specs:

	inch	cm		inch	cm
Overall Width	32.50	82.55	Base Width	32.50	82.85
Overall Depth	19.00	48.26	Base Depth	19.00	48.26
Overall Height - low	28.00	71.12	Base Profile Height	1.94	4.93
Overall Height - high	44.00	111.76	Top Width	32.00	81.28
			Top Depth	17.00	43.18

i Solid Surface top available - contact customer service



	lb	kg
Weight Capacity ¹ EVENLY DISTRIBUTED	100	45.36
Shipping Weight ⁸	36	16.33

Options:

Thermofoil Top Shape [‡]	add	Optional [‡]	add	Vanity with Mirror [‡]	add	Wire Basket [‡]	add
Rectangle 177TFR-y-xx [§]	\$135	Cup Holder	\$102	177TFx-y-xx-VM [§]	\$210	177TFx-y-xx-BK [§]	\$210
Kidney 177TFK-y-xx	\$135	<i>vanity with mirror and wire basket 177TFx-y-xx-VM-BK</i>					

[Product Info](#)

[finish options page 39*](#)

175 Overbed Table with Heavy-Duty Steel, Infinite Stop, Low Profile 'H' Base and Thermofoil Top



base: A919000i[§] top: OS3DLOBTF-xx[§]
\$495 \$260

*order separately / ships separately**

Specs:

	inch	cm		inch	cm
Overall Width	32.00	81.28	Base Width	29.30	74.40
Overall Depth	22.00	55.88	Base Depth	18.20	46.20
Overall Height - low	29.00	73.66	Base Profile Height	1.94	4.93
Overall Height - high	43.00	109.20	Top Width	35.00	88.90
Range Stops	automatic/infinite		Top Depth	18.00	45.72



	lb	kg
Weight Capacity ¹ EVENLY DISTRIBUTED	100	45.36
Shipping Weight ⁸	43	19.50

Options:

Thermofoil Top ^{‡*}	add	Talas Cherry Finish	add
Flat Edge OS3DLOBTF-xx [§]	\$0	OS3DLOBTx-TC [§]	\$32
Raised Edge OS3DLOBTR-xx	\$8	<i>Talas Cherry Finish Option at Additional Charge</i>	

[Product Info](#)

[finish options page 39*](#)

Color Variations perceived texture and coloring may vary on finished product

Laminate Top - Rectangular	Thermofoil Top - Rectangular	Thermofoil Top - Kidney	Raised Edge	Cup Holder	Wire Basket	Vanity Drawer w/ Mirror
					10 x 6 x 9 inch 25.40 x 15.24 x 22.86 cm	

102BAR Bariatric Step Stool, Slip Resistant

102BAR
\$220

Specs: †



[Product Info](#)

	inch	cm
Overall Width	20.00	50.80
Overall Length	26.50	67.31
Overall Height	8.00	20.32
Platform Width	14.75	37.47
Platform Length	19.00	48.26

	inch	cm
Step Height	8.00	20.32
Steel Tubing	1.00	2.54
Steel Tubing Finish	epoxy coating	
Frame Finish	gray / black	

Packaged	1/ea
Ships Assembled	yes
 Weight Capacity¹ EVENLY DISTRIBUTED	1000
Shipping Weight⁸	9

100BAR Bariatric Chrome Step Stool with Slip Resistant Mat

100BAR
\$175

Specs: †



[Product Info](#)

	inch	cm
Overall Width	14.00	35.56
Overall Length	17.00	43.18
Overall Height	9.50	24.13
Platform Width	11.00	27.94
Platform Length	14.00	35.56

	inch	cm
Step Height	9.50	24.13
Steel Tubing	1.00	2.54
Steel Tubing Finish	chrome plated	
Frame Finish	gray / black	

Packaged	1/ea
Ships Assembled	yes
 Weight Capacity¹ EVENLY DISTRIBUTED	500
Shipping Weight⁸	11

100 Chrome Step Stool with Slip Resistant Mat

100
\$99

Specs: †



[Product Info](#)

	inch	cm
Overall Width	14.00	35.56
Overall Length	17.00	43.18
Overall Height	9.00	22.86
Platform Width	11.00	27.94
Platform Length	14.00	35.56

	inch	cm
Step Height	9.00	22.86
Steel Tubing	1.00	2.54
Steel Tubing Finish	chrome plated	
Frame Finish	gray / black	

Packaged	2/cs
Ships Assembled	yes
 Weight Capacity¹ EVENLY DISTRIBUTED	250
Shipping Weight⁸	14



500 Pedestal Exam Light with Chrome-Plated Steel Base and Adjustable Goose Neck



500
\$233

Specs: †

	inch	cm			W	hp(M)	
Overall Height	74.00	187.96	Steel Shaft Finish	chrome-plated	Max LED Equivalent	10	0.014
Height Range -low	34.50	87.63	Pharmacy Shade exterior finish	chrome-plated	Max CFL Equivalent	13	0.018
Height Range -high	53.50	135.89	Pharmacy Shade interior finish	white baked enamel	Max Incandescent	60	0.082
Base Diameter	11.00	27.94	Gooseneck Finish	chrome-plated			
Gooseneck Arm	12.00	30.48	Steel Shaft	adjustable			
			Height Adjustment method	SOK-IT-GUARD™	Ships Assembled		no
Power Cord 3-wire grounded	8.00	2.44	Lamp Head adjustable	360°			
UL Approved	yes		Light Bulb Included	no	Shipping Weight ⁸	34	15.42

Product Info

130-2 Chrome-Plated Steel IV Pole with 4-Leg Steel Base, Rubber Caps, and 2-Hooks



130-2
\$119

Specs: †

	inch	cm					
Base Diameter	22.00	55.88					
Overall Height -low	51.00	129.54	Ram Style Hooks	2	Steel Tubing Finish	chrome-plated	
Overall Height -high	90.00	228.6			Height Adjustment	twist-lock	
Rubber Wheels	2.00	5.08	Packages Shipped	2			
Steel Tubing	1.00	2.54	Ships Assembled	no	Shipping Weight ⁸	13	5.90

Product Info

Footnotes:



- * ships freight: request shipping quote from customer service
- ◇ ships ground: no additional charge
- hyperlink - additional product information and ordering
- page link - jumps to referenced page
- ‡ xx in the product number is the placeholder for vinyl / finish choice
- † standard rounding method used
- § representative image shows this option
- * options available at additional cost & may alter weight specs - order requirements noted
- 1. maximum weight capacity - evenly distributed - refer to user manual
- 2. procedure arm is not compatible / available with either backrest option

- 3. variance: ± 0.5 inch (± 1.27 cm)
- 4. supports a single, left side, table only
- 5. requires factory installation include item with purchase
- 6. pullout footplate not compatible / available with extended legrest option
- 7. interior dimensions
- 8. approximate measurement, options may alter spec
- 9. product is not CFFA-HealthCare 201 Certified
- 10. submit required sample size indicated to:
Intensa COM Department, 1810 S. Elm Street, High Point, NC 27260.
- 11. adjust height spec: +4 inch (+10.16 cm)
- 12. replacement cylinders may be ordered separately at additional cost



UNIVERSAL BED SOCKET IV POLE

Easily slides into IV Pole Attachment. Internal rubber bumpers ensure silent operation and reduce shock. Constructed of durable, anodized aluminum tubing. Adjustable height from 35 to 65 inches (88.90 - 165.10 cm).

Item No. 2101A **add \$64**

IV POLE ATTACHMENT

For use with Bed Socket IV Pole. FR601 & FR597 series require factory installation - add item with order.

Item No. FR601IVMNT⁵ (shown left)

■ IN5601 series

Item No. FR597GIVMNT^{4,5}

■ IN597 series

Item No. 5782RG (shown right)

■ IN566 ■ IN577 ■ IN587 series



WALLSAVER ARM CAPS

Plastic border prevents scarring is a standard feature. No changes to the item number

■ IN566 ■ IN577 ■ IN587 series **add \$0**



FULLY UPHOLSTERED ARMS

For this option add **PA** to the end of the item number

■ IN566 ■ IN577 ■ IN587 series **add \$0**

HEADREST COVER

Protects and extends the life of recliners vinyl. Provides added protection from hair oil, hairspray, coloring agents, etc.

Item No. HRC566xxx ■ IN566 series

Item No. HRC587xxx ■ IN587 series

Item No. HRC577xxx ■ IN577 series

Item No. HRC601xxx ■ IN601 series



HEADREST SKIRT

Easily mounts to recliner back and provides added protection to head-section vinyl from hair oil, hairspray, coloring agents, etc.

Item No. HRSxxx **add \$78**

■ IN566 ■ IN577 ■ IN587

■ IN597 ■ IN601 series



FOLD DOWN SIDE TABLES

Smooth surface edges makes cleaning and disinfection easier. Side tables are level with armrest in the up position.

Item No. FR597GTBLMNT^{4,5}

■ IN597 — single left side table only



SWING-AWAY FOLDING TABLE

360 degree swivel table. Scratch resistant table top is durable and easy to clean. Fold down and out of the way when not needed.

Item No. SAFT20-597 **add \$447**

■ IN597 series

Item No. SAFT20 **add \$447**

■ IN601 series

Item No. SAFT10 **add \$447**

■ IN566 ■ IN577 ■ IN587 series



UNIVERSAL HEADREST PILLOW

The unique shape features an angled design that offers the patient a separate level of support depending on whether the small contour or large contour is placed under the patient's head. Hook and loop attachment allows for variable positioning along back of recliner.



Item No. FR5659xxxUS **add \$98**

LEG EXTENSIONS

Fits all round tube frame designs. Raises recliner by four inches (10.16 cm).

Item No. FR5699G **add \$40**

■ IN566 ■ IN577 ■ IN587 series





Order Numbering Multiple Options & Vinyl Codes

Combining Options

Most single options show the full order number. When combining multiple options, pull each option code from the option order number [or as shown below]. To create the full order number, begin with the main item number (972), then each option code (CB) (-FR), a dash (-), seat height, a dash (-), vinyl/finish codes (see right), a dash (-), then wheel choice.

Order Number for 972 stool with crescent back, footring, medium seat height, with soft wheels in Tier 1 Violet vinyl : 972CB-7-FR-VI-SW.

Order Number for 972 stool with procedure arm, footring, with standard medium seat height & hard wheels in Tier 1 Violet vinyl : 972PA-FR-VI.

Customer's Own Material

Due to the nature of our product, some COM material may not perform or upholster to our standards. Intensa will not be responsible for materials that do not meet our standards. All customer material must be approved prior to manufacture. Submit required sample size indicated to: Intensa COM Department, 1810 S. Elm Street, High Point, NC 27260.

Vinyl/Finish Codes

Most vinyl or finish codes follow the product code (see left) and start with a dash (-). See Approved Upholstery Selections [↗](#)

Approved Customer Vinyls: add '-COM' followed by '-manufacturer name-pattern name-color name'.

- 972CB-7-FR-COM-DesignTex-Chasam-Blue-SW

Tier 1 Vinyls: add the two-letter color code. Example using Alabaster.

- 972CB-7-FR-AL-SW

Tiers 2 - 10 Vinyls: add the tier number '-T#' followed by '-manufacturer name-pattern name-color name'.

- 972CB-7-FR-T3-CFStinson-Tessuto-Lime-SW

OBT Finishes: add the base code letter (-y) followed by the two-letter top finish code (-xx). 175 does not have an optional base color.

- 179-G-MC
 - 175-SO
- See Finish Selections [↗](#)

Physician & Laboratory Seating Options

T-Style Arm Rest [†]	add	Rotating Procedure Arm ^{†2}	add	Upholstered Backrest ^{†2}	add	Crescent Backrest ^{†2}	add	Chrome Footring [†]	add
### -T-xx	\$187	### PA-xx	\$450	### B-xx	\$175	### CB-xx	\$315	### -FR-xx	\$110
adjustable, urethane		adjustable, upholstered, rotates 360°		adjustable, black plastic outer cover		adjustable, upholstered		with twist lock feature	



Casters [†] set of five	add
Hard Wheel standard	\$0
Soft Wheel Urethane Tread ### -SW (hard floors)	\$38
Soft Wheel Self Braking ### -SBCSW (roll weighted)	\$40
Soft Wheel Braking ### -BCSW (roll unweighted)	\$50
Stationary Glides ### -GLD (±2 inch height vs caster)	\$0



Gas Cylinders			inch	cm
5 [five] ^N	###-5	Height	5.50	14.00
7 [seven]	###-7	Height	7.87	20.00
10 [ten]	###-10	Height	10.26	27.00



^N not pictured

Additional Info Notices & Links

Pricing & Availability

current at time of publication & subject to change without notice

Contact Customer Service

for most current information, pricing, freight quotes, and support +1.336.884.4096 | intensa@grahamfield.com

Color Variations

perceived texture and coloring may vary on finished product monitor and/or printer settings will affect appearance of color.

[↗](#) California Prop 65 Warnings

[↗](#) Intensa Limited Warranties

[↗](#) Intensa Terms & Conditions

[↗](#) Green Initiative



175 Top Finishes [XX]



AR ARIA



CO COCOBALA



SC SHAKER CHERRY



SO SOLAR OAK



TC TALAS CHERRY



WC WILD CHERRY



WM WINDSOR MAHOGANY

177-179 Top Finishes [XX]



CL CANDLELIGHT



CP CHOCOLATE PEAR



DV DOLCE VITA



HM HARD ROCK MAPLE



MO MEDIUM OAK

177-179 OBT Base Colors [Y]



A ALMOND



G GRAPHITE GREY



N NEVADA

175 Recommended Care & Cleaning^{R,S}



Clean the surface with a damp cloth or sponge and a mild soap or detergent. Rinse thoroughly with warm water and wipe dry.



Difficult stains such as coffee or tea can be removed using a mild household cleaner/detergent and a soft bristled brush, repeating as necessary.



If a stain persists, use a paste of baking soda and water and apply with a soft bristled brush. Light scrubbing for 10 to 20 strokes should remove most stains. Although baking soda is a low abrasive, excessive scrubbing or exerting too much force could damage the surface.



Stubborn stains that resist any of the above cleaning methods may require the use of undiluted household bleach or nail polish remover. Apply to the stain and let stand no longer than two minutes. Rinse thoroughly with warm water and wipe dry. May repeat if stain appears to be going away and the laminate color has not been affected.



After cleaning rinse surfaces using a water-dampened cloth, and dry thoroughly as some cleaning agents are corrosive by nature, and if not rinsed their residue can cause permanent damage

Recommended Household Cleaners^{R,S}

Windex® • Glass Plus® • Fantastik® • Mr. Clean® • Formula 409® • Isopropyl Alcohol

Warning

Prolonged exposure to bleach will cause discoloration. Always rinse surfaces after cleaning! Failure to rinse after cleaning can cause damage; even if a small amount of cleaning solution remains on the surface.

179 - 177 Recommended Care & Cleaning^{R,S}



ALWAYS rinse overbed table with a water-dampened cloth to remove any cleaning agent residue, and dry thoroughly



Recommend using a mild cleaning agent or detergent when cleaning overbed table. Some cleaning agents are corrosive by nature, and if not rinsed their residue can cause permanent damage



If other hospital-grade cleaning agents are used, follow the instructions provided by the manufacturer



After cleaning rinse surfaces using a water-dampened cloth, and dry thoroughly as some cleaning agents are corrosive by nature, and if not rinsed their residue can cause permanent damage

A Few Notes of Caution

- *Acidic or abrasive cleaners* can damage laminate surfaces; do not use them.
- *Hair, textile and food dyes* can cause permanent stains. If spilled, wipe it up immediately with dish-washing detergent or an all-purpose cleaner.
- *Steel wool and other abrasive pads* will damage the laminate face. Do not use them for cleaning.
- *Drain and Toilet bowl cleaners* contain harsh chemicals that can cause permanent damage. If spilled, wipe up immediately, wash surface with soapy water and rinse several times.

^R recommendations do not serve as a product guarantee. use manufacturer instructions ^S the term "safe" herein, pertains to surface effects and not in reference to human health

Intensa TIER 2




AL ALABASTER



AQ AQUA GREEN



AR ARTICHOKE



BL BLACK



BU BLUE



BT BRITISH TAN




DO DOE



DB DUSTY BLUE



GR GRAY



MO MOSS



NA NAVY



NO NORTHWOODS



OG OCEAN GRAY



PM PALM



PL PEARL




PU PLUM



ST SADDLE TAN



SH SHADOW



WH SUPER WHITE



TL TEA LEAF



TE TEAL



VI VIOLET



WG WARM GRAY



WW WEDGEWOOD

Olympus TIER 6



OLY204 DUNE



OLY202 CORIANDER



OLY213 LAUREL



OLY210 ADMIRAL



OLY01 MELON



OLY09 FERN



OLY208 SHOYA

Natural World TIER 6



NTL221 TIMBER



NTL223 MARIGOLD




NTL222 JUNIPER



NTL228 WALNUT

Raffia TIER 4




SCP2270 ACADEMY




SCP2150 SHOW



SCP2230 CINEMA



SCP2210 SILVER SCREEN



SCP2130 FLICK



SCP2190 PREMIER



SCP2310 PRODUCER

Tessuto TIER 3



SEN224 VENETTO



213 BELLINI



SEN220 MALTA



219 MARBLE



SEN222 LIME



229 MURANO



234 PALAZZO



217 SENECA



SEN225 MINK



233 ADRIATIC



228 MARCO




231 MAJESTY



230 TRUFFLE

Galveston TIER 6



GV-004 TWILIGHT



GV-003 CACTUS



GV-007 CREAM



GV-011 OXBLOOD



GV-016 WINCHESTER



GV-010 CEDAR



GV-026 STUCCO



GV-017 DESERT



GV-014 STREAM

Belmont TIER 4



BL-002 CITRUS



BL-009 CLEMINTINE



BL-007 SISAL



BL-004 OPAL



BL-024 TIDAL



BL-016 STONE



BL-000 GRAVEL




BL-027 PATINA



BL-017 DUNE

Icon TIER 3



IC-003 CELERY



IC-010 COFFEE



IC-013 GRASS



IC-034 MARINE



IC-001 SANGRIA



IC-017 SAND



IC-014 OCEAN



IC-019 TERRA



IC-027 ALMOND

Midway TIER 2



MW-007 CEMENT



MW-009 COPPER



MW-017 SAND



MW-023 SPA



MW-003 MOSS



MW-013 ALGAE



MW-004 INDIGO



MW-000 STONE



MW-016 HAZE



MW-027 WHEAT



MW-002 SUNSHINE



MW-014 SKYLIGHT



MW-006 STORM

Reference this chart for Price Tier | Care Codes detailed on page 43 | Link Icon connects to pattern swatch selection, spec sheet, care guide, and more

Stinson®	tier	care	link	Stinson®	tier	care	link	Stinson®	tier	care	link	Stinson®	tier	care	link	Stinson®	tier	care	link
Alamo	2	W	🔗	Core	2	W	🔗	Level 2.0	3	W	🔗	Polygon	6	W	🔗	Smoothie 2.0	3	W-S	🔗
Amplify	2	W-S	🔗	Criss Cross 2.0	3	W-S	🔗	Loophole 2.0	3	W	🔗	Promenade 2.0	5	W	🔗	Spirit Millennium ⁹	3	W	—
Arbor	3	W	🔗	English Pub ⁹	3	W-S	🔗	Lunaria ⁹	6	W	🔗	Raffia 2.0	4	W	🔗	Stat Con ⁹	6	W	🔗
Ardent	2	W	🔗	Essential	2	W-S	🔗	Matrix ⁹	3	W-S	🔗	Rally ⁹	2	W-S	🔗	Stitch in Time ⁹	1	W-S	🔗
Avant	3	W-S	🔗	Fizz 2.0	3	W	🔗	Meteor	6	W	🔗	Rialto ⁹	3	W-S	🔗	Strand ⁹	3	W-S	🔗
Bandwidth	6	W	🔗	Flanders	3	W-S	🔗	Mica 2.0 ⁹	3	W	🔗	Riverstone	3	W	🔗	Stride	3	W	🔗
Blend	2	W	🔗	Flip ⁹	6	W	🔗	Moccasin 2.0	2	W	🔗	Rogue 2 ⁹	2	W-S	🔗	Sumi	6	W-S	🔗
Blink Ombre	5	W	🔗	Frequency 2.0	5	W	🔗	Montana 2.0	2	W	🔗	Round Trip	6	W	🔗	Symphony Classic ⁰	2	W	—
Box Chain	2	W	🔗	Fuse 2.0	2	W-S	🔗	Natural World 2.0	6	W-S	🔗	Route	2	W	🔗	Tendrill 2.0	2	W	🔗
Budget Black ⁰	1	W	—	Gabby	2	W-S	🔗	Naugasoft ⁹	2	W-S	🔗	Runway	3	W-S	🔗	Tessuto 2.0	3	W	🔗
Burkshire	4	W-S	🔗	GeoMet	2	W	🔗	Nishiki 2.0	3	W-S	🔗	Satchel	3	W-S	🔗	Tilt	2	W	🔗
Cassidy 2.0	2	W-S	🔗	Grand Sierra 2.0 ⁰	2	W	—	Orbit 2.0	3	W	🔗	Scout 2.0	3	W-S	🔗	Tinker	6	W	🔗
Cesto	3	W-S	🔗	Grille 2.0	3	W	🔗	Overlay ⁹	3	W-S	—	Script 2.0	4	W	🔗	Touchstone	3	W-S	🔗
Channel 2.0	3	W-S	🔗	Hopsack	3	W	🔗	Oxford 2.0	3	W	🔗	Sea of Clouds	5	W-S	🔗	Vector	4	W	🔗
Chroma ⁹	5	W-S	🔗	Kilkenny Tweed 2.0 ⁰	3	W-S	—	Painted Moon	5	W-S	🔗	Sentiero	3	W-S	🔗	Walk the Line 2.0	6	W	🔗
Cirque 2.0	5	W	🔗	Kusari 2.0 ⁹	3	W	🔗	Patriot Pluss II 2.0 ⁰	2	W	—	Sentry-HC 2.0	3	W	🔗	Watercolor Grid	6	W	🔗
Citadel	3	W-S	🔗	Lara	2	W	🔗	Pebble Trail	5	W	🔗	Shantung 2.0	2	W	🔗	Weft	2	W	🔗
Clean Slate 2.0	3	W	🔗	Laredo ⁹	2	W-S	🔗	Ping Pong	6	W	🔗	Sirocco	6	W	🔗	Whisper 2.0 ⁰	2	W	—
Colorguard ⁰	2	W	—	Lattice	2	W	🔗	Pod	6	W	🔗	Sleek 2.0 ⁹	3	W-S	🔗	Ziggy 2.0	3	W	🔗
Connect 2.0	2	W	🔗																

PVC Free ^F	tier	care	link	PVC Free ^F	tier	care	link	Patty Madden ^F	tier	care	link	Patty Madden ^F	tier	care	link	Patty Madden ^F	tier	care	link
Arioso ⁹	6	W	🔗	Modena 2.0 ⁹	7	W-S	🔗	Affinity 2.0	6	W	🔗	Odeon	9	W	🔗	Vortexia 2.0	9	W	🔗
Attaché ⁹	3	W	🔗	Petal Dots ⁹	8	W-S	🔗	Banyan ⁰	3	W	—	Olympus 2.0	6	W	🔗	Waterwand	6	W	🔗
Bond ⁹	7	W	🔗	Rain Pool ⁹	9	W-S	🔗	Bari 2.0	6	W	🔗	Origami	9	W	🔗	Wicker 2.0	6	W	🔗
Chaos ⁹	10	W	🔗	Savile ⁹	7	W-S	🔗	Bubblerake 2.0	9	W	🔗	Orvilla ⁹	8	W	🔗	Winnly	5	W	🔗
Composite ⁹	10	W-S	🔗	Sedona ⁹	C	W-S	🔗	Cassiopeia ⁹	8	W-S	🔗	Parasail	7	W-S	🔗	Ybor	3	W	🔗
Corium ⁹	7	W-S	🔗	Silicone-Avail ⁹	9	W-S	🔗	Eli 2.0	3	W	🔗	Perspective 2.0	6	W	🔗				
Didi ⁹	6	W-S	🔗	Silicone-Downtown ⁹	9	W-S	🔗	Hashtag	5	W	🔗	Prisma	3	W	🔗				
Effervescence 2.0 ⁹	7	W-S	🔗	Silicone-Uptown ⁹	9	W-S	🔗	Intersect	9	W	🔗	Raya	6	W	🔗				
Etch ⁹	7	W	🔗	Solitaire ⁹	2	W-S	🔗	Juxy 2.0	6	W	🔗	Revolve 2.0	9	W	🔗				
Fiberdrill ⁹	5	W	🔗	Stage ⁹	5	W-S	🔗	La Scalina ⁹	8	W	🔗	Ripple 2.0	9	W	🔗				
Fiddlesticks ⁹	10	W-S	🔗	Stallion ⁹	6	W-S	🔗	Lava 2.0	9	W	🔗	Rivera ⁹	3	W	🔗	Michael Graves ^F	tier	care	link
Flare ⁹	10	W-S	🔗	Staple ⁹	2	W	🔗	Lens ⁹	8	W-S	🔗	Scoprena	9	W	🔗	Entrada ⁹	2	W	🔗
Galileo ⁹	7	W-S	🔗	Sundance ⁹ w/StaKleen™	5	W-S	🔗	Live Wire	6	W	🔗	Swizzle	7	W-S	🔗	Excursion ⁹	1	W	🔗
Impress ⁹ w/StaKleen™	5	W-S	🔗	Tender ⁹	C	W-S	🔗	Martini 2.0	6	W	🔗	Twist	9	W	🔗	Latitude ⁹	1	W	🔗
Katsura ⁹	7	W-S	🔗	Terra ⁹	8	W	🔗	Median	6	W	🔗	Verge	9	W	🔗	Valise ⁹	1	W-S	🔗
Kid ⁹	5	W-S	🔗	Volterra-HC ⁹	C	W-S	🔗												
Mesa ⁹	C	W-S	🔗	Wyatt ⁹	3	W-S	🔗												

🔗 CF Stinson's Complete Selection of Approved Vinyl Upholstery

Mayer Fabrics	tier	care	link	Mayer Fabrics	tier	care	link	Mayer Fabrics	tier	care	link	Mayer Fabrics	tier	care	link	Spradling®	tier	care	link	
Belmont	4	W-S	🔗	Imprint	2	CFC	🔗	Mingle	3	W-S	🔗	Truss	4	W-S	🔗	Guardian	2	S	🔗	
Caressa	3	CFC	🔗	Interlace	4	CFC	🔗	Mode	4	CFC	🔗	Verona	3	W-S	🔗	Savoy	2	S	🔗	
Catwalk	4	CFC	🔗	Key Largo	2	W-S	🔗	Noble	2	CFC	🔗	Zip	2	CFC	🔗	Cassara	2	S	🔗	
Chaps	4	W-S	🔗	Kinsey	2	W-S	🔗	Outback	2	W-S	🔗				Modena	1	S	🔗		
Cirque	4	CFC	🔗	Ludlow	2	W-S	🔗	Ryder	2	W-S	🔗				Masialas	2	S	🔗		
Etch	2	CFC	🔗	Lyric	2	W	🔗	Silverweave	4	CFC	🔗				Dolce ⁹	2	S	🔗		
Explorer	2	W-S	🔗	Matka	3	W-S	🔗	Sketch	3	W-S	🔗				Patriot Plus	2	S	🔗		
Galveston	6	CFC	🔗	Maverick	2	W-S	🔗	Sonnet	4	W-S	🔗				Renaissance	2	S	🔗		
Icon	3	W-S	🔗	Midway	2	W	🔗	Spartan	2	W-S	🔗									

🔗 link to pattern swatch, spec, and care guide F by CF Stinson D product has been renamed 9 product is not CFFA-HealthCare tier C Call for pricing

Guardian TIER 6



1325 ARMOR



1305 CONVOY



1323 MOAT



1319 SIREN



1324 SENTRY



1311 FENCE



1318 STEEL




1304 LATCH

Savoy TIER 6




9517 EUCALYPTUS



9524 BLACK



409 WARM TAUPE



9516 INDIGO


Cassara TIER 6



9611 AMETHYST



9613 CORK



9607 PATINA



9604 INDIGO

Modena TIER 6



2013 TITAN



2022 ELEPHANT



2002 BLUEBIRD



2008 SEAGRASS



2014 INK



2019 NUTMEG



2007 FRESH



2001 FARRO

Masiala TIER 6



321 BLUESTAR



320 DEMIN



311 UNDERWATER



313 WOODLET



304 CALM



306 BLACK WALNUT



309 BERRY BUSH




317 DAMASCUS




Dolce TIER ?

			
6727 METEOR	6717 MOSS	6724 JET	6725 INDIGO

Patriot Plus TIER ?

			
432 IMPERIAL BLUE	8602 OATMEAL	8604 CONCRETE	8567 DOVE
			
454 ROYAL BLUE	427 BLUERIDGE	8584 TAUPE	8516 CHESTNUT
			
8611 TERRA COTTA	8504 SAND	8533 PORT	8626 EUCALYPTUS

Renaissance TIER ?

			
9214 STEEL BLUE	9201 COBBLESTONE	9205 GYPSUM	9208 WINEBERRY
			
9203 PYRITE	9207 ICE BLUE	9210 CYPRESS	9212 BLUE JADE



Color Variations swatches are digital reproductions. perceived texture and coloring may vary on finished product. monitor and printer settings will effect appearance of color.

Fabric Availability current at time of publication & subject to change without notice

Contact Customer Service for the most current availability





NASSIMI

RECOMMENDED

Care & Cleaning for swatch collections on page 37



Regular Cleaning & Maintenance

Clean the soiled area with mild soap and water, then rinse with fresh water and wipe dry with a clean cloth.



Food Stains / Oils

Wipe affected area with a soft cloth or soft bristle brush using a non-abrasive cleaner, for example Formula 409® or Fantastik®. Rinse with fresh water and wipe dry.



Difficult Stains

Clean the soiled area with 50% a solution of isopropyl alcohol and water. It is important to rinse the cleaned area with fresh water after applying the alcohol solution.



Healthcare / Disinfection

Clean the soiled area with a 10% solution of Household Bleach and water. It is important to rinse the cleaned area with fresh water after applying the bleach solution.

Example

Dirt
Dust
Grime

Ketchup
Chocolate
Coffee
Tea
Juice

Eyeshadow/
Mascara
Lipstick
Crayon
Grease
Mustard

Blood
Urine
Betadine

Disinfectants

In addition to cleaning, healthcare facilities maintain disinfecting regimens to reduce the spread of infections. There is a wide variety of available biocides that include bleaches, peroxides and quaternary ammonium products, among others. They should always be used at the recommended dilution, never in concentrated form. To prolong the life of coated fabrics, the disinfectants should be rinsed off after each use. It should be noted that misuse of disinfectants is the major source of surface deterioration.

For Color Transfer

Certain clothing and accessory dyes (such as those used on denim jeans) may migrate to lighter colors. This phenomenon is increased by humidity and temperature and is irreversible. Please check compatibility when using this product in combination with painted or varnished surfaces

Note: The cleaning products recommended were used in a laboratory setting. This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and referenced cleaning agents. GF Health Products, Inc. has no relationship with the manufacturer of these cleaning products.

* This term and any corresponding data refer to the typical performance in the specific tests indicated and should not be construed to imply the behavior of this or any other material under actual fire conditions



RECOMMENDED

Care & Cleaning for swatch collections on pages 41 & 42

SPRADLING®

Clean within 24 hours to avoid permanent stains in the surface. After this time has elapsed, effectiveness drops to 50%.



We recommend cleaning at least once a week to extend the lifespan of the upholstery, which will prevent dirt and contamination buildup and permanent staining. In the event of a spill, wipe off any excess right away with a dry white cloth or fabric, but careful not to spread to remaining surface.



Clean the whole surface in a circular motion using pH neutral liquid hand soap and water (1 part soap to 9 parts water). Check the Approved Cleaners guide for recommended cleaning products. Wipe with a clean damp cloth to get rid of any excess soap. Repeat as needed and then let dry.



If there is any buildup that cannot be removed with soapy water, a mixture of alcohol and water (1 part isopropyl alcohol to 9 parts water) can be used as a last resort. Remember, using alcohol on a regular basis can affect the properties of the product.

† This section of cleaning guide does not apply to Dolce.



Once a month, you can use a soft-bristle brush for a deeper cleaning. Remember, light colors need more care, more often.

🔗 Cleaning Guides specific for each finish technology. [www.spradling.group]

🔗 Guide to Healthcare Surface Disinfection for the Real World. [www.spradling.group]

Disinfectants

There is a wide variety of available biocides that include bleaches, peroxides and quaternary ammonium products, among others. They should always be used at the recommended dilution, never in concentrated form. To prolong the life of coated fabrics, the disinfectants should be rinsed off after each use.

If you need to disinfect the product because it has come into contact with organic material (blood, urine or fecal matter), check the [disinfection guide](#) for appropriate cleaning products. Clean with a clean white cloth or fabric. If you choose to use a bleach solution, only use 0.5% concentration (1 part bleach to 10 parts water). Rinse and wipe with a clean damp cloth to get rid of any excess and dry.

For Color Transfer

Our products don't transfer onto other materials, but some clothing items (like jeans, microfibers and dark fabrics) could transfer their color to the coated fabrics, depending on how they are dry-cleaned. This happens more with dampness and higher temperatures and is irreversible. Neither GF Health Products, Inc. nor Spradling International, Inc.® is responsible for any damages caused in this manner. Please check compatibility when using this product in combination with painted or varnished surfaces.

Note: The cleaning products recommended were used in a laboratory setting. This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and referenced cleaning agents. GF Health Products, Inc. has no relationship with the manufacturer of these cleaning products.

* This term and any corresponding data refer to the typical performance in the specific tests indicated and should not be construed to imply the behavior of this or any other material under actual fire conditions

Stinson **RECOMMENDED** MAYER Care & Cleaning Codes

for swatch collections on pages 37 & 38 | use chart links on page 40 for full details

W

Clean with water-based foam or water-based cleaners only. Do not use solvents. We recommend using a professional cleaning service.

Bleach solutions of up to 10% may be used for more difficult stains.

W-S

For water-based stains use water-based cleaners. For oil-based stains use solvent-based cleaners. We recommend using a professional cleaning service.

Bleach solutions of up to 20% may be used for more difficult stains.

Crypton® Cleaners are recommended to clean many of Stinson's fabrics. These specifically developed cleaners are water-based, pH neutral, biodegradable, easy to use, and do not contain bleach, solvents or harsh chemicals. 🔗 [Crypton® Cleaners](#)

S

Clean the whole surface in a circular motion using pH neutral liquid hand soap and water (1 part soap to 9 parts water). Wipe with a clean damp cloth to get rid of any excess soap. Repeat as needed and then let dry.

* This term and any corresponding data refer to the typical performance in the specific tests indicated and should not be construed to imply the behavior of this or any other material under actual fire conditions

D

After disinfection, hand dish-washing soap and water or Formula 409® may be used with a soft cloth or damp sponge to remove chemical residue.

After cleaning, rinse area with clean water, then dry with a clean, lint-free cloth. Water rinse recommended after each disinfectant application.

CFC

Clean with mild soap and water or Fantastik® or Formula 409® may be used with a soft cloth. More stubborn stains may be eliminated using a cotton cloth soaked in 70-91% of isopropyl alcohol. Wipe without spreading stain to surrounding area. Repeat with vigorous rubbing in a circular motion. For textured material, a soft brush is recommended. Rinse with fresh water and wipe dry.

Note: The cleaning products recommended were used in a laboratory setting. This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and referenced cleaning agents. GF Health Products, Inc. has no relationship with the manufacturer of these cleaning products.

Put a little
COLOR
 in your
 exam
 room
 with **intensa**
 by graham-field



990 SERIES Physician Exam Stools	4
980 SERIES Physician Exam Stools	5
970 SERIES Physician Exam Stools	6
960 SERIES Physician Exam Stools	7
950 SERIES Physician Exam Stools	8
940 SERIES Physician Sit/Stand Stool	8
930 SERIES Physician Exam Stools	9
920 SERIES Physician Exam Stools	9
870 SERIES Laboratory Seating	10
860 SERIES Laboratory Seating	11
840 SERIES Laboratory Seating	13
830 SERIES Laboratory Seating	16
810 SERIES Laboratory Seating	17
740 SERIES Hydraulic Surgical Seating	19
HSS500 SERIES Hydraulic Surgical Seating	19
IN601 SERIES Pivot-Arm Clinical Care Recliners	20
IN597 SERIES Ortho-Biotic™ II Clinical Care Recliners	21
IN587 SERIES Extra Wide Clinical Care Recliners	22
IN566 SERIES Deluxe Clinical Care Recliners	23
IN577 SERIES Multi-Care Recliners	24
600 SERIES Blood Draw / Phlebotomy Seating	24
460 SERIES Exam Room Treatment Tables	25
420 SERIES Exam Room Treatment Tables	25
410 SERIES Exam Room Treatment Tables	26
RT400 SERIES Replacement Exam Table Tops	28
RT300 SERIES Replacement Exam Table Tops	29
RT200 SERIES Replacement Exam Table Tops	29
RT100 SERIES Replacement Exam Table Tops	29
200 SERIES Wallsaver Patient / Guest Seating	30
170 SERIES Overbed Tables	31
100 SERIES Step Stools	32
500 SERIES Exam Light	33
130 SERIES IV Pole	33
Clinical Care Recliner Options	34
Physician & Laboratory Seating Options	35
Order Numbering Multiple Options & Codes	35
Customer's Material Requirements	35
Additional Information Notices & Links	35
Terms & Conditions Limited Warranties Prop 65 Warnings Green Initiative	
Overbed Table Finish Swatches & Care	36
Tier 1 Vinyl Upholstery Standard Collection Swatches	37
Stinson® Vinyl Upholstery Pattern Collection Swatches	38
Mayer Vinyl Upholstery Texture Collection Swatches	39
Approved Vinyl Upholstery Tier, Care, Link, & Swatch Chart	40
Spradling® Vinyl Upholstery High Performance, Green Collection Swatches	41
Spradling® Vinyl Upholstery Performance Collection Swatches	42
Spradling® Vinyl Upholstery Premiere Collection Swatches	42
Vinyl Upholstery Recommended Care, Cleaning & Disinfection	43
Vinyl Upholstery Composition, Finish, & Performance Comparison Chart	44





A Better Healing Environment Better Results



published **November 2025**

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

GF1800088 RevT25

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

© 2015, GF Health Products, Inc. All Rights Reserved | Intensa, Inc. is a wholly owned subsidiary of GF Health Products, Inc. | Stinson and Zeron are registered trademarks of C. F. Stinson, LLC.
| Intensa, Dura-Gel, and SOCK-IT-GUARD are trademarks of GF Health Products, Inc. | Spradling, Neo, EcoSense, PBG+, Permaguard, Permablok, Permablok Plus, and Sureguard are trademarks of Spradling Holdings, Inc.
| Crypton is a registered trademark of Crypton LLC | Tente is a registered trademark of TENTE International | Formula 409 is a registered trademark of The Clorox Company. | Fantastic is a registered trademark of S.C. Johnson & Son, Inc.

