



Laboratory Seating & Casegoods

published *January 2025*





COLOR

It's in our DNA

Intensa is changing the laboratory seating environment by offering one of the largest color stories and incorporating color into the lab space. We have multiple textile collections with extensive color palettes, matching patterns, and non-PVC vinyl. We also accept pre-approved customer material to match your specific design needs.



Our laboratory seating selection includes many designs, price points, and functions to suit most lab environments. Some designs are offered with different base options, single or dual lever controls, clean room/ESD application, and a fold-down back option for ease of storage under the lab desk.



ESSENTIALS
offer the flexibility to change your
laboratory workflow
when progress requires change





research

searching for the best solution for your lab?


LOOK NO FURTHER



ESSENTIALS innovative, modular design lets you easily and affordably select a laboratory casework solution to match your specific workflow, space, chemistry, and aesthetics.

Plus, Essentials' desktop productivity elements enables you to achieve superior ergonomics and staff productivity in the "transaction-intensive" lab environment.

Choose Essentials for efficiency specific to your laboratory.





pharmaceutical

organize your pharmacy inventory like never before
REFILL AS NEEDED



ESSENTIALS organization for your pharmacy enhances “reach & pick” efficiency for filling accuracy and productivity.

Essentials elegantly organizes active filling, compounding, checking, order entry, and dispatch stations by locating supplies, inventory, technology, and pharmaceuticals at your fingertips. Maximize your technicians’ and pharmacists’ active “reach & pick” zone with Essentials.





from workstations to wet labs,

INTENSA HAS YOUR CASEWORK SOLUTIONS



ESSENTIALS modular casework is ideal for your exam rooms, procedure rooms, nurse stations, and almost any critical workrooms.

Only Essentials gives you the design flexibility to elect finishes, surfaces, components, and desktop productivity elements to provide an inviting, efficient healthcare delivery and interaction space for your patients and staff.





education

we demand a lot from our students

FURNISH THEIR LABS TO MEET THE CHALLENGE



Today's **STEM** curriculum and performance learning labs require innovative interiors that motivate educators and students while providing space-efficient and flexible classrooms and labs.

Labs need to be multi-discipline and multi-purpose learning spaces to enable instructional “*swing spaces*” that can go from geometry to computer to engineering instruction labs daily.



840DC SERIES self skin ergonomic laboratory chair
 seat & back tilt | dual lever hardware | height adjustable back rest | footring

- ★ optional lower range adjustable seat height (-7-)
- ★ optional t-style arm rest (-T-)



841DC
 Black Nylon
 Composite Base



842DC
 Brushed Aluminum
 Base with Toe Caps



843DC
 Polished Aluminum
 Base



848DC
 Black Tubular Five-Leg
 Steel Base with Glides



urethane

self skin ergonomic laboratory chair

840 SERIES

All **840 SERIES** chairs are designed for the harshest environments. They have an all-steel inner frame, are wrapped in self-skin urethane, and are impervious to harsh chemicals. They are also impact-resistant, offer Pneumatic Height Adjustment, 20 inch footring, and safely support 350 pounds¹.





840SC SERIES self skin ergonomic laboratory chair
single lever hardware | lateral adjustable back rest | footring

- ★ optional lower range adjustable seat height (-7-)
- ★ optional t-style arm rest (-T-)



841SC
Black Nylon
Composite Base



842SC
Brushed Aluminum Base
with Toe Caps



843SC-T
Polished Aluminum Base
with optional Arm Rest



848SC
Black Tubular Five-Leg
Steel Base with Glides

830 SERIES upholstered ergonomic laboratory chair
 seat & back tilt | dual lever hardware | height adjustable back rest | footring

- ★ optional low range adjustable seat height (-7-)
- ★ optional t-style arm rest (-T-)



831
 Black Nylon
 Composite Base



832
 Brushed Aluminum
 Base with Toe Caps



833
 Polished Aluminum
 Base



838-7
 Black Tubular Five-Leg
 Steel Base with Glides
 and Low Height
 Range Cylinder

lab chair

upholstered ergonomic laboratory chair

810 & 830 SERIES

All **810 & 830 SERIES** chairs offer pneumatic height adjustment and can safely support 350 pounds¹. The seat and back vinyl can be produced in any of the 24 industry standard colors, approved C.F. Stinson coverings, or customer's own material.





810 SERIES upholstered ergonomic laboratory chair
single lever hardware | lateral adjustable back rest | footring

- ★ optional low range adjustable seat height (-7-)
- ★ optional t-style arm rest (-T-)



811
Black Nylon
Composite Base



812
Brushed Aluminum Base
with Toe Caps



813
Polished Aluminum Base



818
Black Tubular Five-Leg
Steel Base with Glides

Color Variations Swatches are digital reproductions.
Perceived texture and coloring may vary on finished product.
Monitor or printer settings will effect appearance of color.



870 SERIES upholstered ergonomic laboratory chair
d-ring handle release hardware | lateral adjustable back rest | footring

- ★ optional low range adjustable seat height (-7-)
- ★ optional crescent style back rest (CB-) *standard back rest pictured*



871
Black Nylon
Composite Base



872
Brushed Aluminum
Base with Toe Caps



873
Polished Aluminum
Base



878
Black Tubular Five-Leg
Steel Base with Glides



lab stool

upholstered ergonomic laboratory stool

860 & 870 SERIES

All **860 & 870 SERIES** stools offer pneumatic height adjustment, lateral back rest adjustment, and can safely support 350 pounds¹. The seat and back vinyl can be produced in a variety of patterns and colors, approved C.F. Stinson coverings, or customer's own material.





860 SERIES upholstered ergonomic laboratory chair
single lever hardware | lateral adjustable back rest | footring

- ★ optional low range adjustable seat height (-7-)
- ★ optional crescent style back rest (CB-) *standard back rest pictured*



861
Black Nylon
Composite Base



862
Brushed Aluminum Base
with Toe Caps



863
Polished Aluminum Base



868
Black Tubular Five-Leg
Steel Base with Glides

Color Variations Swatches are digital reproductions.
Perceived texture and coloring may vary on finished product.
Monitor or printer settings will effect appearance of color.



Seating Features & Options

HW | hard wheel caster — standard



CB | upholstered

SW | soft wheel urethane tread caster option★



B | standard adjustable urethane or upholstered back rest

BCSW | soft wheel braking caster option★



SBCSW | soft wheel self braking caster option★



T | adjustable urethane t-style arm rests option★

GLD | stationary glide — standard 5-leg base



lab seating

features and options

NO TWO LAB SETTINGS ARE ALIKE.

Along with our wide range of colors and pattern that are available, we offer our customers additional options to customize their workspace for optimal performance and design.



FR | standard 20 inch chrome foot rest with twist lock





8x3 | polished aluminum base



8x2 | brushed aluminum base with toe caps



8x1 | black nylon composite base



8x8 | black tubular five-leg steel base with glides



adjustable crescent back rest option★

with black plastic outer cover



EDC | standard electrostatic discharge components

7 | low height range cylinder option★



Laboratory Seating & Casegoods



published **January 2025**

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

GF2400126RevB24

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

© 2017, GF Health Products, Inc. All Rights Reserved | Intensa, Inc. is a wholly owned subsidiary of GF Health Products, Inc. | Intensa is a registered trademark of GF Health Products, Inc.

