

TREATMENT TABLES

413 SERIES EXAM ROOM TABLE

Exam Room Treatment Table with Shelf and Flat Top
Intensa offers a collection of treatment tables to meet any need in the exam room. The tables are offered with an upholstered back and an option for zero, one or two drawers. Intensa offers an extensive collection of vinyls to add color to the treatment table.

413 Features:

- Powder coat frame finish - Foghorn
- Optional one or two pass through drawers
- Flat top

 CF Stinson Pre-Approved Fabrics

 413 Product Info



413

413 Options:

	item number [†]
No Drawers	413-xx
Single Drawer	413D-xx
Double Drawer	413DD-xx

options available at additional cost



413D



413DD

413 Specifications:

	inch	cm		inch	cm		yd	m
Overall Width	28.00	71.12	Paper Roll Holder Widths Available	18.00	45.72	COM Requirement ¹	4	3.66
Overall Length	76.00	193.04		19.00	48.26			
Overall Height	32.00	81.28		20.00	50.80			
Drawer Width ⁷	17.75	45.09	Foam Thickness	1.75	4.45		lb	kg
Drawer Length ⁷	19.25	48.90				Weight Capacity EVENLY DISTRIBUTED	500	226.80
Drawer Height ⁷	5.00	12.70	Powder Coat Finish		foghorn	Shipping Weight	210	95.25

^{*} Ships: Freight request additional shipping quote from Intensa [†] -xx in the product number is the placeholder for vinyl choice ¹ Submit COM sample for approval
² Adjust Shipping Weight +10 lb (+4.54 kg) per drawer ³ Interior Dimensions



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

