

EXAM & TREATMENT TABLES

RT100 SERIES

REPLACEMENT TOP

Replacement Top with Foot Pad for Ritter® by Midmark® 100 Exam Tables

Intensa offers a replacement top program for your existing exam tables. We provide tops that meet certain manufacture model number specifications. The new top is made of completely new foam and wood substrate and can be upholstered in a variety of vinyl colors and patterns. Revitalize/renew your exam rooms without the expense of new tables!

RT100 Features:

- Replaces Ritter® by Midmark® 100 Exam Table
- Paper roll holder and foot pad included
- Includes mounting hardware

 [CF Stinson Pre-Approved Fabrics](#)

 [Replacement Tops Overview](#)


 [RT100 Product Info](#)

RT100 Options:

	Item no.†
Replacement Top	RT100-xx



RT100 Specifications:


	inch	cm	inch	cm	yd	m	
Overall Width	27.50	69.85	Foot Rest Width	20.00	3	2.74	
Overall Length	55.00	139.70	Foot Rest Length	18.00			
Paper Roll Holder	included		Mounting Hardware	included		Shipping Weight	55 lb / 24.95 kg

*Ships: Freight request additional shipping quote from Intensa †-xx in the product number is the placeholder for vinyl choice ¹ Submit COM sample for approval



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

Manufactured in the USA

