

SEATING

740 SERIES SURGICAL STOOL

740 Features:

- Foot lever adjustment release
- Hydraulic height adjustment
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 **CF Stinson Pre-Approved Fabrics**

 **740 Options Detail and Photos**

Product Info:

 **744**  **748**



744




748

740 Options:

standard seat height	item no. [†]	with backrest [†]	crescent back [†]
Black Aluminum Base Black Powder Coat Column	744-xx	744B-xx	744CB-xx
Polished Aluminum Base Silver Powder Coat Column	748-xx	748B-xx	748CB-xx

options available at additional cost

740 Specifications:

	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter <i>with casters</i>	24.00	60.96	COM Requirement¹	1	0.92
Overall Height <i>-low</i>	35.50	90.17	Height Adjustment <i>method</i>	hydraulic			lb	kg
Overall Height <i>-high</i>	42.50	107.95	Rotation Lock	on seat		Weight Capacity EVENLY DISTRIBUTED	350	158.76
Height Release <i>hardware</i>	foot control		Powder Coat Finish	column		Shipping Weight	45	20.41

Adjustable Height Seat^{††}


	inch	cm		inch	cm
Standard 74x-xx			Seat Height -low	19.50	49.53
			Seat Height -high	26.50	67.31

[†] -xx in the product number is the placeholder for vinyl choice ^{††} Replacement Cylinders may be ordered separately at additional cost  ¹ Submit COM sample for approval



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

© 2016, GF Health Products, Inc. All Rights Reserved | Intensa, Inc. is a wholly owned subsidiary of GF Health Products, Inc. | Intensa is a registered trademark of GF Health Products, Inc.



Manufactured in the USA