

## SURGICAL SEATING

# 740 SERIES HYDRAULIC SURGICAL STOOL

### 740 Features:

- Foot lever adjustment release
- Hydraulic height adjustment
- High density fire resistant foam
- 3/4" (1.91 cm) thick composite substrate with steel reinforcement plate

 Approved Upholstery Selections

Product Info:

 744  748




### 740 Options:

standard seat height	item no. <sup>†</sup>	with backrest <sup>‡</sup>	crescent back <sup>‡</sup>
<b>Black Aluminum Base</b> Black Powder Coat Column	744-xx	744B-xx	744CB-xx
<b>Polished Aluminum Base</b> Silver Powder Coat Column	748-xx	748B-xx	748CB-xx

### Adjustable Height Seat<sup>††12</sup>

		inch	cm		inch	cm
5 [five] <sup>§</sup>	74x-xx	Seat Height -low			Seat Height -high	
		19.50	49.53		26.50	67.31

### 740 Specifications:


<sup>†</sup>	inch	cm		inch	cm		yd	m
Seat Diameter	16.00	40.64	Base Diameter with casters	24.00	60.96	COM Requirement <sup>10</sup>	1	0.92
Overall Height -low	35.50	90.17	Height Adjustment method	hydraulic			lb	kg
Overall Height -high	42.50	107.95	Rotation Lock	on seat		Weight Capacity <sup>1</sup> EVENLY DISTRIBUTED	350	158.76
Height Release hardware	foot control		Powder Coat Finish	column		Shipping Weight <sup>8</sup>	45	20.41

\* Ships Freight; request additional shipping quote <sup>†</sup> -xx in the product number is the placeholder for vinyl choice <sup>§</sup> representative image shows this option  
<sup>†</sup> standard rounding method used <sup>10</sup> submit COM sample for approval <sup>1</sup> maximum weight capacity - evenly distributed - refer to user manual  
<sup>8</sup> approximate measurement & options may alter weight specs <sup>12</sup> replacement cylinders may be ordered separately at additional cost   
 \* options available at additional cost - accessories sold separately



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | intensa@grahamfield.com | intensa.net

© 2016, GF Health Products, Inc. All Rights Reserved | Intensa, Inc. is a wholly owned subsidiary of GF Health Products, Inc. | Intensa is a registered trademark of GF Health Products, Inc.

