

LABORATORY

840DC SERIES LABORATORY SEATING

All 840 series chairs are designed for the harshest environments. They have an all-steel inner frame, are wrapped in self-skin urethane, and are impervious to harsh chemicals.

840DC Features:

- Dual lever adjustment release
- Lever operated pneumatic height adjustment
- Standard back rest with height adjustment
- Lever operated seat tilt and back tilt
- Impact-resistant
- Standard 20 inch (50.80 cm) footring
- Steel substrate
- High density fire resistant foam impervious to harsh chemicals
- Available in black urethane only

[840DC Options Detail and Photos](#)

Product Info:


[841DC](#) [842DC](#) [843DC](#) [848DC](#)

840DC Options:

base with standard seat height	item no ¹	with arm rest ²
Black Nylon Composite	841DC-xx	841DC-T-xx
Brushed Aluminum <i>and Toe Caps</i>	842DC-xx	842DC-T-xx
Polished Aluminum	843DC-xx	843DC-T-xx
Five Leg Tubular Steel⁴ <i>Footring and Glides</i>	848DC-xx	848DC-T-xx

options available at additional cost

840 Specifications:

	inch	cm		inch	cm		
Seat Width	18.00	45.72	Back Width	16.00	40.64	Seat & Back substrate	steel
Seat Depth	17.00	43.18	Back Height	12.00	30.48	Seat & Back cover material	urethane
Overall Height with backrest -low	37.50	97.79	Base Diameter with casters	24.00	60.96		lb kg
Overall Height with backrest -high	48.00	124.46	Height Adjustment method	pneumatic		Weight Capacity EVENLY DISTRIBUTED	350 158.76
Height Release hardware	lever		Back Adjustment direction	height		Shipping Weight	45 20.41

Adjustable Height Seat⁺⁺

	inch	cm		inch	cm	
Standard	84xDC-xx	Seat Height -low	22.50	57.15	Seat Height -high	33.00 83.82
7 [seven]	84xDC-7-xx	Seat Height -low	18.50	46.99	Seat Height -high	25.50 64.77

⁺ -xx in the product number is the placeholder for vinyl choice [†] Replacement Cylinders may be ordered separately at additional cost [‡] Specification Variances
¹ Submit COM sample for approval ⁴ 20 inch (50.80cm) footring; casters add 2" (5.08 cm) in height
[♦] Ships: Ground at no additional charge

item	inch	cm	height	ship weight
	inch	cm	lb	kg
841	—	—	-1	-0.45
842	-1	-2.54	—	—
848	+2	+5.08	+4	+1.81

